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### Illumina Complete Long Read Prep with Enrichment, Human

#### DesignStudio Assay Design Tool

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1. Illumina. Illumina Complete Long Read Prep, Human data sheet. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf). Published 2022. Accessed September 22, 2023.
2. Illumina. Comprehensive whole-genome sequencing with Illumina Complete Long Read Prep, Human technical note. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf). Published 2022. Accessed September 22, 2023.
3. Roessler K. Illumina Complete Long Reads software analysis workflow for human WGS. <https://www.illumina.com/science/genomics-research/articles/complete-long-read-software-analysis.html>. Published 2023. Accessed September 22, 2023.
4. Wagner J, Olson ND, Harris L, et al. Curated variation benchmarks for challenging medically relevant autosomal genes. *Nat Biotechnol.* 2022;40(5):672-680. doi:10.1038/s41587-021-01158-1
5. PharmGKB. VIPs: Very Important Pharmacogenes. [pharmgkb.org/vips](https://www.pharmgkb.org/vips). Accessed September 22, 2023.
6. National Library of Medicine. GTR: Genetic Testing Registry. Precision HealthPGx Panel (25 Genes). [ncbi.nlm.nih.gov/gtr/tests/593428/](https://www.ncbi.nlm.nih.gov/gtr/tests/593428/). Updated November 29, 2022. Accessed September 22, 2023.

7. Pratt VM, Everts RE, Aggarwal P, et al. Characterization of 137 Genomic DNA Reference Materials for 28 Pharmacogenetic Genes: A GeT-RM Collaborative Project. *J Mol Diagn.* 2016;18(1):109-123. doi:10.1016/j.jmoldx.2015.08.005
8. Miller DT, Lee K, Abul-Husn NS, et al. ACMG SF v3.1 list for reporting of secondary findings in clinical exome and genome sequencing: A policy statement of the American College of Medical Genetics and Genomics (ACMG). *Genet Med.* 2022;24(7):1407-1414. doi:10.1016/j.gim.2022.04.006
9. Kulski JK, Suzuki S, Shiina T. Human leukocyte antigen super-locus: nexus of genomic supergenes, SNPs, indels, transcripts, and haplotypes. *Hum Genome Var.* 2022;9(1):49. doi:10.1038/s41439-022-00226-5
10. Bekritsky MA, Colombo C, Eberle MA. Identifying genomic regions with high quality single nucleotide variant calling. Published 2021. Accessed August 30, 2023.
11. Illumina. Illumina Human Comprehensive Panel data sheet. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-enrich-hu-comp-panel-data-sheet-m-gl-02191/long-read-hu-comp-panel-data-sheet-m-gl-02191.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-enrich-hu-comp-panel-data-sheet-m-gl-02191/long-read-hu-comp-panel-data-sheet-m-gl-02191.pdf). Published 2024. Accessed January 26, 2024.

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© 2024 Illumina, Inc. All rights reserved.