


Figure 2: Methylation Assay Reproducibility And Differential Methylation Detection



Comparison of methylation profiles between lung cancer and matching normal tissue. The b value (i.e., the methylation ratio measured for all 1,536 CpG sites) obtained from one replicate experiment is plotted against that obtained from another technical replicate experiment. DNA samples derived from lung adenocarcinoma G12022 (center) and its matching normal tissue (left) demonstrate the reproducibility of replicated assays. The right panel shows the comparison between the normal and corresponding adenocarcinoma samples.

3. Technology

Implementation of the methylation profiling assay on the SNP genotyping platform

The GoldenGate Assay for SNP genotyping was adapted for DNA methylation detection. Non-methylated cytosines (C) are converted to uracils (U) when treated with bisulfite, while methylated cytosines remain unchanged. Because the hybridization properties of uracil are similar to that of thymine (T), the detection of the methylation status of a particular cytosine can be carried out following bisulfite treatment by using a genotyping assay for a C/T polymorphism.

Four probes are designed for each CpG site: two allele-specific oligos (ASO) and two locus-specific oligos (LSO). Each ASO-LSO oligo pair corresponds to either the methylated or unmethylated state of the CpG site (Figure 1). The gap between the ASO and LSO allows difficult sequences or ambiguous bases in CpG islands of interest to be avoided. This flexibility is particularly important for methylation studies because of a decrease in sequence complexity as a result of bisulfite treatment. If other CpG sites are present in close proximity to the target CpG site, we make the assumption that they have the same methylation status as the site of interest. This design hypothesis is based on previous reports and our own bisulfite sequencing results, in which a majority (> 90%) of the adjacent CpG sites were shown to be co-methylated or co-unmethylated^{4,5}. Therefore, the GoldenGate Assay for Methylation provides not only a methylation measurement at the targeted CpG site, but also “inferred” methylation information at the adjacent underlying CpG sites.

This design strategy is used in methylation-specific PCR primer design⁶ and other microarray-based DNA methylation analysis⁷. While many CpG sites exist within each CpG island, we only select those for which robust assays can be designed, i.e., assays with a design score above a predetermined threshold.

The EZ DNA methylation kit (Zymo Research, Orange, CA) is used for bisulfite conversion of all DNA samples, according to the manufacturer’s recommendations. For each sample, 500 ng to 2 µg of genomic DNA can be used for each bisulfite conversion. An aliquot of

the converted DNA (corresponding to 250 ng starting gDNA) is then used to assay up to 1,536 CpG sites simultaneously on an array. The subsequent assay procedures are similar to those described previously for standard SNP genotyping¹. For more details, please see the GoldenGate Assay for Methylation System User Guide (Illumina Part No. #11228975).

The methylation status of an interrogated CpG site is determined by calculating the intensity of the methylated (M) and unmethylated (U) alleles. Background intensity computed from a set of negative controls is subtracted from each analytical data point. The ratio of fluorescent signals is then computed from the two alleles, where methylation value (b) = $\text{Max}(M, 0)/[\text{Max}(M, 0) + \text{U}(M, 0) + 100]$.


b provides a continuous measure of levels of DNA methylation in samples, ranging from 0 in the case of completely unmethylated sites to 1 in completely methylated sites. An absolute value is used in the denominator of the formula, as a compensation for any negative values which may arise from global background subtraction (i.e., over-subtraction). A constant bias of 100 is added to the denominator to regularize b when both U and M values are small.

It is important to note that b only provides a relative but not an absolute methylation level measurement in a sample. For example, for any given CpG site, the two ASO-LSO oligo pairs may not hybridize to their corresponding target sequences at exactly the same efficiency; therefore, when $b = 0.5$, it may not mean that methylation level is 50%. Furthermore, b may not be absolutely zero at some assayed CpG sites even when they are not methylated, presumably because of some degree of cross-hybridization of the assay oligos to other non-specific genomic regions. The lack of absolute quantification of this assay, however, does not undermine its capability for differential methylation detection.


Methylation Assay Reproducibility and Differential Methylation Detection

We obtained highly reproducible DNA methylation profiles between technical replicates (Figure 2), with an average r^2 of 0.98 ± 0.02 when the b values were compared. We observed that the standard devia-

Figure 3: Methylation Detection In GDNA Mixtures



A. Female genomic DNA was diluted with male genomic DNA at ratios of 5:95, 10:90, 20:80 and 50:50. Two sets of mixtures were made and measured: M1 (male NA10923)/F1 (female NA10924) and M2 (male NA07033)/F2 (female NA06999). Methylation levels of six X-chromosome linked genes were calculated as the average of several (1-5) CpG sites for each gene.



B. Methylation profiles of individual CpG sites (only four are shown). Error bars represent the standard deviation of b values calculated from four replicate experiments performed with the first set of mixtures (NA10923/NA10924).

tion of the b value obtained for all the 1,536 CpG sites across four replicates was less than 0.06 in 99% of cases. Therefore, we estimate that our method can discriminate levels of methylation (b) that differ by as little as 0.17 ($1.96 \times \sqrt{2} \times 0.06$). The Illumina BeadStudio software can compare the b values between samples or sample groups and perform a differential methylation analysis. The difference in methylation levels between samples in a pair is quantified by a “Diff Score;” a Diff Score equal to 20, 30, or 40 corresponds to a significance p value equal to 0.01, 0.001, or 0.0001, respectively. To estimate the false positive rate for differential methylation detection, we measure methylation levels at 1,536 CpG sites for the same sample in replicate experiments, and usually less than 1% of CpG sites show up with a Diff Score $\geq |20|$.

Sensitivity of the Methylation Assay

To assess the assay sensitivity, we measured the methylation status of six X-linked housekeeping genes—EFNB1, ELK1, FMR1, G6PD, GPC3, and GLA—together with several hundred genes in male and female genomic DNA samples. Female genomic DNA was diluted into male genomic DNA at ratios of 5:95, 10:90, 20:80 and 50:50, prior to bisulfite conversion. Methylation levels from 5:95 and 0:100 mixtures could be unambiguously distinguished from each other (i.e., the maximum b value in the 0:100 mixtures was less than the minimum b value in the 5:95 mixtures; Figure 3A, 3B). This indicates that the GoldenGate Assay for Methylation can detect as little as 2.5% methylation for well-performing CpG sites. On average, the assay can detect change in b of ≥ 0.2 with 95% confidence, which is consistent with our previous estimation.

