Infinium® II Assay Workflow

Illumina's Infinium II Assay provides unlimited multiplexing for whole-genome genotyping applications with manual or automated workflow

INTRODUCTION

The Infinium II Whole-Genome Genotyping Assay (Figure 1) is designed to interrogate a large number of SNPs at unlimited levels of loci multiplexing. Using a single bead type and dual color channel approach, the Infinium II Assay scales genotyping from 10,000 to hundreds of thousands of SNPs per sample. Illumina's optional Laboratory Information Management System (LIMS) and automation ensure positive sample

tracking while reducing time required and labor costs.

ASSAY PROTOCOL

The DNA sample used for this assay is isothermally amplified in an overnight step (Figure 1; Step 1 and Step 2). This amplification has no appreciable allelic partiality. Additionally, a relatively low DNA sample requirement of 750 ng is sufficient to assay over 500,000 SNP loci. The amplified product is then fragmented by a controlled enzy-

matic process that does not require gel electrophoresis (Step 3). After alcohol precipitation and resuspension of the DNA (Step 4), the BeadChip is prepared for hybridization in the capillary flow-through chamber (Step 5); samples are applied to BeadChips and incubated overnight. The amplified and fragmented DNA samples anneal to locus-specific 50-mers (covalently linked to one of over 500,000 bead-types) during the hybridization step (Step 6). One bead type corresponds

FIGURE 1: INFINIUM II ASSAY PROTOCOL


FIGURE 2: INFINIUM II ASSAY WORKFLOW

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* These are total incubations, but several can be nested within each other

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to each allele per SNP locus. After hybridization, allelic specificity is conferred by enzymatic base extension. Products are subsequently fluorescently stained (Step 7). The intensities of the beads' fluorescence are detected by the Illumina BeadArray Reader (Step 8), and are in turn analyzed using Illumina's software for automated genotype calling (Step 9).

Figure 2 shows the estimated hands-on time required for completing the Infinium II Assay using the Illumina BeadStation 500 system. With Illumina's optional Laboratory Information Management System (LIMS) to ensure positive sample tracking, the Infinium II Assay is a robust protocol with a straight forward workflow that can be automated or processed manually.

ADDITIONAL INFORMATION

For more information about Infinium II or other products and services from Illumina, visit our website or contact technical support at the address below.

Illumina, Inc.
Customer Solutions
 9885 Towne Centre Drive
 San Diego, CA 92121-1975
 1.800.809.4566 (toll free)
 1.858.202.4566 (outside the U.S.)
 techsupport@illumina.com
 www.illumina.com

FOR RESEARCH USE ONLY

© 2006 Illumina, Inc.
 Illumina, BeadArray, BeadXpress, CSPro, Sentrix, Array of Arrays, GoldenGate, Infinium, iSelect, DASL, Sentrix, VeraCode, and Making Sense Out of Life, are trademarks or registered trademarks of Illumina. All other names and marks are the property of their respective owners.
 Pub. No. 370-2006-027 07Dec06