


# The value of non-invasive prenatal testing (NIPT).

A supplement for a  
Genetic Counselor's flipbook



# NIPTs use cell-free DNA.




# NIPTs use cell-free DNA.

- The maternal blood sample contains mixture of both maternal and fetal cfDNA.
- Both maternal and fetal fragments are counted and analyzed in NIPT.
- Aneuploidies are detected by comparing the amount of chromosomal material against a set of reference chromosomes.

# MPS — a deeper approach to sequencing.


Maternal and fetal cfDNA



Alignment and counting

```
CCCTTAGCGCTTTAACGTACGTAAAACCCCTT
AACGTACGTAAAACGGGGTCAAAGTTCCC
GACTTAAAATCGGAATCGATGCCCAAACCTT
AATCGATGCCCAAACGGGGTCAAAGTTCCC
```

cfDNA sequenced via  
Massively Parallel Sequencing (MPS)




# MPS — a deeper approach to sequencing.

- Maternal blood sample contains a mixture of both maternal and fetal cfDNA.
- MPS, or Massively Parallel Sequencing, analyzes cfDNA fragments from a maternal blood sample.
- Sequences are aligned to genome to identify the chromosome origin of the fragment.
- Aneuploidies are detected by comparing the amount of chromosomal material for a specific chromosome in a patient sample against a set of reference chromosomes.

# Detecting fetal chromosome abnormalities.

## Chromosomes



# Detecting fetal chromosome abnormalities.

- The veriFi® prenatal test detects trisomies 21, 18, and 13 in both singleton and twin pregnancies.

# Detecting sex chromosome abnormalities.

## Chromosomes






# Detecting sex chromosome abnormalities.

- In singleton pregnancies, sex chromosome aneuploidies can also be detected. This test is optional based on patient and physician preference.
  - Sex chromosome aneuploidies detected include 45,X, 47,XXX, 47,XXY, 47,XYY.
- In twin pregnancies, the presence or absence of Y is detected. This test is optional based on patient and physician preference.

# How the verifi<sup>®</sup> test is reported.

## Results



# Test results.

## There are three possible results for the autosomes:


- Aneuploidy Detected, Aneuploidy Suspected (borderline value), Aneuploidy Not Detected.

## There are two possible results for the sex chromosome analysis in singletons:

- Aneuploidy Detected and Aneuploidy Not Detected.
  - If aneuploidy is not detected, fetal sex will be reported.
  - For twin pregnancies, the presence or absence of Y will be reported.

Patients who receive an Aneuploidy Detected or Suspected result should be offered a follow-up discussion of results and options of invasive testing (CVS or amniocentesis) for confirmation.

# Typical prenatal testing sequence.



# Typical prenatal testing sequence.


## Pre-test Counseling Considerations

- Aneuploidy testing should be an informed choice.
- The verifi prenatal test screens for aneuploidies of chromosomes 21, 18 and 13. Sex chromosome testing is optional.
- Any test offering should be based on clinical context.
- CVS and amniocentesis remain options for all pregnant women.

## Post-test Counseling Considerations

- 3 possible results for each chromosome tested.
  - **No Aneuploidy Detected** – This result is reassuring but does not ensure an unaffected pregnancy.
  - **Aneuploidy Suspected/Borderline Value** – This result is suggestive of aneuploidy. There is a small possibility that the results might not reflect the chromosomes of the fetus but may reflect changes to the placenta or the patient.
  - **Aneuploidy Detected** – This result is consistent with aneuploidy. There is a small possibility that the results might not reflect the chromosomes of the fetus but may reflect changes to the placenta or the patient.
  - Genetic counseling and CVS/amniocentesis are recommended for Suspected and Detected results.
- Open neural tube defect screening is not a part of NIPT testing and should be offered separately.
- All results should be reviewed in the context of the clinical, medical, and family history.

# NIPT profile: the verifi<sup>®</sup> test by the numbers.




What you can't put a number on –  
the reassurance of knowing

# How is the verifi<sup>®</sup> prenatal test different?


- The verifi<sup>®</sup> test is an easy, in-office blood test – only one tube of blood is required.
- It can be performed as early as 10 weeks in pregnancy.
- It tests for trisomies 21, 18 and 13 and has an option for sex chromosome aneuploidies (Monosomy X, XXX, XXY, XYY).
- Test reports are usually available in 2 – 4 business days after sample receipt.
- It has the lowest failure rate among NIPTs (0.1%).

Because of this, the verifi<sup>®</sup> test is able to provide reliable answers about the most common chromosome abnormalities.

# What is mosaicism?



Cell division



● Normal cells


● Cells with genetic change



# What is mosaicism?

- Mosaicism is a biological phenomenon which results from the presence of cells with differing genetic compositions in the same individual. It usually occurs during the process of mitotic cell division.
- Mosaicism can be an underlying biological cause of discordant NIPT results.

# Types of mosaicism.




# Types of mosaicism.


There are three types of mosaicism:

- **Generalized mosaicism** – Presence of two or more karyotypically different cell lines in both the placenta and the fetus. In cases of generalized mosaicism, there is a possibility of a false positive or a false negative NIPT result depending on the origin of the cfDNA.
- **Confined placental mosaicism** – Presence of two or more karyotypically different cell lines that are confined to the placenta and not present in the fetus. In cases of confined placental mosaicism, there is a possibility of a false positive NIPT result.
- **Fetal mosaicism** – Presence of two or more karyotypically different cell lines that are present in the fetus but not present in the placenta. In cases of fetal mosaicism, there is a possibility of a false negative NIPT result.

# Confined Placental Mosaicism.



## Laboratory techniques used for Chorionic Villus Sampling (CVS)



# Confined Placental Mosaicism (CPM).

- Confined placental mosaicism (CPM) is defined as the presence of two or more karyotypically different cell lines that are confined to the placenta and not present in the fetus.
- The cfDNA that is analyzed in NIPT arises from the cytotrophoblasts; thus an NIPT result is most similar to a CVS direct preparation result.
- **Type 1 CPM** – Mosaicism is confined to trophoblastic cells.
- **Type 2 CPM** – Mosaicism is confined to the cells from chorionic villus stroma.
- **Type 3 CPM** – Mosaicism is present in both: trophoblastic and stroma cells.
- Type 1 and Type 3 CPM could result in false positive NIPT results.