


Figure 2: Data Correlation between Automated and Manual Library Preparation Methods



Data generated by Nextera Rapid Capture Exome automated and manual library preparation methods at UNC were compared by selecting a single control sample (NA12878) from the 12-plex pooled libraries described in the experimental design section. The number of unique variants (indels and SNVs) and the corresponding transition/transversion and heterozygous/homozygous ratios for the unique calls identified by each method are indicated to the left of each Venn diagram. As demonstrated, both methods call a very low number of unique indels and SNVs, indicating high correlation between the data sets. Similar results were observed for other sequenced samples (data not shown).

Table 1: Performance Metrics Compared between Manual and Automated Nextera Rapid Capture Methods

| Processing Method | Sample Plexity [*] | Estimated Insert Size (Median Bases) [†] | Number of Reads Aligned to the Target [‡] | Mean Target Coverage [§] | Duplication ^{**} | Target Bases at 10x ^{††} | Target Bases at 20x ^{‡‡} | Total SNVs Called ^{§§} | Total Indels Called ^{¶¶} |
|-----------------------------|-----------------------------|---|--|-----------------------------------|---------------------------|-----------------------------------|-----------------------------------|---------------------------------|-----------------------------------|
| Automated | 12 | 175 | 54,399,184 | 102.96 | 8.32% | 96.09% | 89.99% | 23,476 | 748 |
| UNC (Manual) | 12 | 163 | 104,949,959 | 118.20 | 17.05% | 95.95% | 90.08% | 23,407 | 743 |
| Illumina (Reference) | 12 | 206 | 71,394,650 | 140.40 | 11.10% | 97.90% | 94.50% | 24,968 | 894 |

^{*} Number of pooled samples sequenced

[†] Median base size of the insert fragments

[‡] Total number of reads that uniquely align to regions defined by the target manifest

[§] Average number of regions covered in the target manifest file

^{**} Percent of duplicate reads due to possible PCR artifacts generated during amplification ^{††} Percent of target bases covered by at least 10 bases

^{‡‡} Percent of target bases covered by at least 20 bases

^{§§} Number of SNVs that contain the PASS filter in the target regions

^{¶¶} Number of indels that contain the PASS filter in the target regions

HiSeq 2500 System performance metrics indicate comparable results between manual and automated (using the Tecan Freedom EVO NGS workstation) Nextera Rapid Capture library preparation methods.

Learn More

To obtain the Nextera Rapid Capture automation method for the Tecan Freedom EVO NGS workstation* discussed in this application note, visit www.tecan.com/ngs.

For questions regarding this application note, send inquiries to NGSPrep@tecan.com.

* The application described in this document is provided as-is, is not an Illumina product, and is not accompanied by any rights or warranties. Customers using or adapting this application should obtain any licenses required and materials from authorized vendors. While customer feedback is welcomed, this application is not supported by Illumina Technical Support, Field Application Scientists, or Field Engineers.

[†] "Illumina-qualified" indicates that analysis by Illumina has shown that libraries prepared with the method perform comparably to those prepared manually.

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

FOR RESEARCH USE ONLY

© 2014 Illumina, Inc. All rights reserved.

Illumina, HiSeq, Nextera, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners. Pub. No. 770-2013-061 Current as of 11 November 2014

To learn more about the Nextera Rapid Capture Exome and Expanded Exome kits, visit www.illumina.com/nrc.

References

1. Coriell biobank (www.coriell.org/research-services/biobanking)
2. HiSeq Analysis Software (support.illumina.com/sequencing/sequencing_software/hiseq-analysis-software.ilmn)

