

# BeadStudio GT Module System Information and Benchmark Performance

System information and benchmark performance data to be used as a guide for determining hardware, software, and memory requirements for analysis of Illumina DNA BeadChips.

## INTRODUCTION

This document includes the following information:

- **Hardware, operating system, and memory requirements** for Illumina DNA BeadChip analysis using the BeadStudio Genotyping Module
- **Project creation guidelines** to optimize genotyping projects for performance or project size
- **Report file sizes** to determine storage needs for final reports and locus summary reports
- **Benchmark performance data** for common operations performed on computers with different specifications

## HARDWARE REQUIREMENTS

Illumina recommends the following hardware for your BeadStudio genotyping projects.

TABLE 1: HARDWARE REQUIREMENTS

MINIMUM SYSTEM	RECOMMENDED SYSTEM
Intel Pentium IV or newer processor	Intel Pentium IV or newer processor
32-bit system	64-bit system
4GB RAM	8GB RAM
100+ GB hard drive	100+ GB hard drive

## OPERATING SYSTEM REQUIREMENTS

Illumina recommends the following operating systems for your BeadStudio genotyping projects.

- 32-bit Windows XP
- 64-bit Windows XP
- 32-bit Windows Vista

Illumina has not fully tested BeadStudio on 64-bit Windows 2003 Server, but it has been demonstrated to work.

## PROJECT CREATION GUIDELINES

At project creation, use the following guidelines to optimize a BeadStudio genotyping project for faster performance or for larger project size: For faster performance, use memory-based storage. For larger project size, use file-based storage.

TABLE 2: PROJECT CREATION GUIDELINES

	STORAGE MODE	CALCULATION MODE
Fastest Performance ↓ Larger Project Size	Memory-based	Pre-calculated
	Memory-based	Real-time calculation
	File-based	Pre-calculated
	File-based	Real-time calculation

## REPORT FILE SIZES


Report file sizes vary based on the Illumina BeadChip used and the number of samples in a project.

TABLE 3: FINAL REPORT AND LOCUS SUMMARY REPORT FILE SIZES

BEADCHIP # SAMPLES	FINAL REPORT FILE SIZE (MB)	LOCUS SUMMARY REPORT FILE SIZE (MB)
HumanHap550 BeadChip 120 Samples	1927 (16 per sample)	97
Human1M BeadChip 120 Samples	4650 (39 per sample)	187
Human610-Quad BeadChip 500 Samples	9196 (18 per sample)	108

**BENCHMARK PERFORMANCE**

The following graphs provide benchmark performance data for operations commonly performed in the BeadStudio Genotyping Module on data generated from Illumina BeadChips.



**ADDITIONAL INFORMATION**

For more information about BeadStudio system information, contact Technical Support at the phone number or email at right. To learn more about Illumina's software solutions, please visit [www.illumina.com](http://www.illumina.com).

**Illumina, Inc.**

**Customer Solutions**

9885 Towne Centre Drive  
 San Diego, CA 92121-1975  
 1.800.809.4566 (toll free)  
 1.858.202.4566 (outside North America)  
[techsupport@illumina.com](mailto:techsupport@illumina.com)  
[www.illumina.com](http://www.illumina.com)

**FOR RESEARCH USE ONLY**

© 2008 Illumina, Inc. All rights reserved.  
 Illumina, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSPPro, iScan, and GenomeStudio are registered trademarks or trademarks of Illumina. All other brands and names contained herein are the property of their respective owners.  
 Pub. No. 970-2008-016 Current as of 06 May 2008

