## illumina

# More Robust Reagents for Use in the Whole-Genome DASL® Assay

Reformatted cDNA synthesis reagents restore high levels of reproducibility across degradationsensitive FFPE samples.

#### Introduction

The cDNA-mediated Annealing, Selection, extension, and Ligation (DASL) Assay allows for expression profiling of partially degraded RNA such as that commonly found in formalin-fixed, paraffin-embedded (FFPE) samples. The Whole-Genome DASL (WG-DASL) HT Assay (>29,000 gene targets) offers researchers the opportunity to analyze hundreds to thousands of RNA transcripts in FFPE samples from as little as 50 ng total RNA, delivering highly reproducible results<sup>1</sup>.

Recent improvements to the cDNA synthesis reagents extend the reagent shelf-life, enabling researchers to continue obtaining the expected high reproducibility across samples. Specifically, this change splits the MCS3 reagent, which is a master mix containing buffer and reverse transcriptase enzyme (RTE), into an MCS4 buffer component and a separate RTE. This new format, MCS4 + RTE, produces equivalent performance to the MCS3 reagent throughout its entire shelf life.



RNA from eight different samples used in the WG-DASL HT assay using either MCS3 or MCS4 + RTE reagents.

FFPE samples (Biochain)

UHR = Universal Human Reference RNA (Stratagene) Brn Ref = Human Brain Reference RNA (Ambion)

### **High-Quality Performance**

To evaluate and compare the performance of MCS4 + RTE with that of the MCS3 reagent, WG-DASL HT experiments using both methods were performed on eight different RNA samples, consisting of six Biochain FFPE samples and two intact MAQC samples [Universal Human Reference (UHR) and Human Brain Reference (BrnRef]). FFPE RNA was extracted using the High Pure RNA Paraffin Kit (Roche), with total RNA input of 200 ng for each FFPE sample and 100 ng for each intact sample. BeadChips were scanned on the iScan System and raw data imported and analyzed using GenomeStudio<sup>®</sup> Data Analysis Software. Results are discussed below.

#### **High Reproducibility**

To assess the reproducibility of the new MCS4 + RTE reagents, all eight RNA samples were processed in duplicate and compared with MCS3 data generated in the same experiment. The raw signal intensities for all 29K probes were then compared between sample replicates for each reagent. The results demonstrate equivalent or improved levels of self-reproducibility for MCS4 + RTE as compared with MCS3, across a range of FFPE sample qualities, with r<sup>2</sup>>0.97 for all samples tested (Figure 1).

#### **High Gene Expression Profile Correlations**

Direct comparisons of the raw signal intensities between experiments performed using MCS4 + RTE and MCS3 yielded correlations of  $r^2$ ~0.99 for intact RNAs. Similar comparisons of FFPE RNAs yielded a range of correlations with  $r^2$ ~0.92–0.95. Representative scatterplots are shown in Figure 2.

#### **Equivalent Sensitivity**

Using intact and FFPE samples, a similar level of sensitivity was obtained on assays done with MCS4 + RTE when compared with those using MCS3. Equivalent numbers of genes were also detected in each sample (Figure 3). Additionally, a high degree of overlap was observed between detected transcripts in samples processed with either MCS4 + RTE or MCS3, with ~99% and ~96% of the detected transcripts shared in intact and FFPE samples, respectively (Figure 4).

#### High Cross-Platform Concordance

To independently validate gene expression profiles generated using the MCS4 + RTE reagents in the WG-DASL HT assay, fold-change comparisons with MAQC data derived from mRNA-Seq on the HiScan<sup>®</sup> and TaqMan<sup>2</sup> platforms were performed. Scatterplot comparisons of the WG-DASL HT assay vs. TaqMan (653 transcripts)



— Figure 3: Equivalent Sensitivity







.

and WG-DASL HT assay vs. mRNA-Seq (13,600 transcripts) data yielded robust correlations of  $r^2 \sim 0.88$  and  $r^2 \sim 0.85$ , respectively, indicating a high degree of measurement concordance between the WG-DASL HT assay (MCS4 + RTE reagent) and the other two gene expression platforms (Figure 5).

#### Conclusions

Splitting the MCS3 reagent in the WG-DASL HT assay into two components, MCS4 buffer and RTE, maintains high assay performance and data quality throughout the product shelf life. Because the gene expression profiles of WG-DASL HT assays performed with MCS3 and MCS4 + RTE are similar, samples processed with both reagents can be analyzed in a single experiment.

#### References

- April C, Klotzle B, Royce T, Wickham-Garcia E, Boyaniwsky T et al. (2009) Whole-genome gene expression profiling of formalin-fixed, paraffin-embedded tissue samples. PLoS ONE 4(12): e8162.
- MAQC Consortium, Shi L, Reid LH, Jones WD, Shippy R, et al. (2006) The Microarray Quality Control (MAQC) project shows inter- and intraplatform reproducibility of gene expression measurements. Nat Biotechnol. 24:1151–61.

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

#### FOR RESEARCH USE ONLY

© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseŠpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. Pub. No. 470-2012-002 Current as of 09 May 2012

