

Qualification of DNA Extracted from FFPE Samples

The TruSight Tumor 26 extension-ligation process relies upon the ability of a DNA polymerase to read across a single-stranded template. Therefore, a quantitative PCR approach was taken to qualify the performance of DNA extracted from FFPE tissue. By comparing the ability of FFPE DNA to be amplified relative to that of a non-FFPE reference gDNA, a ΔCq value can be calculated for each sample. This value can be used to predict its performance with TruSight Tumor 26. A non-FFPE reference gDNA Quality Control Template (QCT) is provided with the Quality Control Primers (QCP). It is important to note that the amplifiable mass (ng) of FFPE DNA capable of generating library product will often be only a fraction of the total amount extracted from FFPE tissues. The exact amount of FFPE DNA input will vary according to the quality of the extracted DNA (Table 1).

To make sure that the ΔCq value accurately reflects the ability of a sample to be amplified, the QCT must be used carefully and consistently to maintain a reliable baseline Cq value. Aliquot 5 μ l of the stock QCT into PCR tube strips for long-term storage to improve pipetting accuracy and template performance. Then, each time the assay is run, 5 μ l of QCT from a single aliquot is added to 495 μ l of DEPC water to make the 1:100 dilution needed as input for qPCR. This can then be run according to the following protocol:

- Vortex the diluted QCT and spin briefly at 280 \times g.
- Dilute 1.5 μ l of each FFPE gDNA with 148.5 μ l DEPC water.
- Aliquot 10 μ l of DEPC water for use as the No Template Control (NTC).
- Prepare the qPCR premix using the following volumes per replicate:

2X qPCR Master Mix*	5.0 μ l
0.5 μ M QCP**	1.0 μ l
DEPC water	2.0 μ l
Total volume per well	8.0 μl

* Volume provided reflects use of Bio-Rad iQ SYBR Green Supermix.
** The QCP primers are prepared as mixed 5 μ M stock, then diluted 1:10 to generate a 0.5 μ M working stock.

To calculate the total volume to premix, 8 μ l is multiplied by (# of samples + QCT + NTC)† 3 replicates each + 10% overfill. The final volume can be transferred to a trough or PCR tube strips.

- Use a multichannel pipette to dispense 8 μ l of premix into each well to be used in the qPCR plate (data was generated with the Bio-Rad CFX384 Real-Time System with clear 384-well plates).
- Pipette 2 μ l NTC into each of 3 wells, 2 μ l QCT into each of 3 wells, and 2 μ l FFPE DNA into each of 3 wells per sample. After each dispense, pipet up and down 3X to ensure transfer.
- Seal the plate and spin briefly at 280 \times g.




Figure 1: Qualification of DNA Samples by PCR—DNA samples described in Table 1 were qualified for use in creating libraries for next-generation sequencing using real-time PCR following parameters stated in the text. The graph shows qPCR amplification curves generated by these samples. Purple: Cq threshold. Blue: DNA samples with abnormal amplification profiles that did not yield libraries suitable for sequencing. Green: DNA samples with profiles passing qPCR qualification that can be used to prepare sequencing libraries.

