




Figure 2: OvineSNP50 BeadChip Probe Spacing



The OvineSNP50 BeadChip provides uniform coverage across the entire ovine genome.

BAC-end sequencing that were subsequently validated using Illumina GoldenGate® Genotyping Assays in over 403 animals from 23 breeds.

### High-Quality Data

The 54,241 SNPs on the OvineSNP50 BeadChip were subjected to rigorous functional testing on multiple breeds to ensure strong performance using the Infinium HD assay. High call rates and accurate genotype calls are required for successful whole-genome association studies. Since complex traits often have relatively small gene effects, potential associations can be missed if the assayed SNP in linkage disequilibrium with the SNP of interest has a low call rate or incorrect genotype call. Illumina ensures that every OvineSNP50 BeadChip offers > 99% average call rate across common breeds.

Table 1 shows the results from internal validation testing of the OvineSNP50 BeadChip content using samples provided by collaborators and the Ovine HapMap Consortium<sup>3</sup>. Illumina scientists and collaborators identified and retained 3,079 loci that appear to have an adjacent polymorphism or deletion among the breeds sampled. Although these loci yielded lower call rates when compared to most loci

Table 1: Performance and Specifications

Parameter	Results*	Product Specification
Average Call Rate	> 99.9%	> 99%
Reproducibility	> 99.9%	> 99.9%
Mendelian Inconsistencies	< 0.01%	< 0.1%

\* Based on genotypes from reference samples

on the panel, they were retained because they may provide biologically relevant information for traits of interest and future improvements in the genome build. These performance and content validation results clearly demonstrate the robust data quality delivered by the product. With such high data quality, the OvineSNP50 BeadChip provides researchers the highest accuracy and reliability for interrogating ovine genotypes in *Ovis aries* breeds.

### Illumina Solutions for Genotyping

The OvineSNP50 BeadChip is compatible with the BeadArray™ Reader, iScan, and HiScan™SQ systems. These array scanners feature high-performance lasers and powerful optical systems that enable rapid scan times and precise assay detection. The HiScanSQ system can also perform Illumina sequencing by synthesis chemistry, the world's most widely adopted next-generation sequencing platform.

The convenient modular design enables researchers to easily build out the system for evolving research needs. An optional Laboratory Information Management System (LIMS) is available to accurately and efficiently track samples. Robotic automation capabilities can be added to improve throughput for labs processing large numbers of samples. With the Infinium Assay workflow, data are processed directly into Illumina's GenomeStudio® software to provide streamlined genotype calling, analysis, and reporting. Researchers can also choose to use Illumina's convenient FastTrack Genotyping service to have samples genotyped and data delivered in a format suitable for GWAS or QTL analysis.

Table 2: BeadChip Content Validation

Species	Breed	Samples	Polymorphic Loci	Mean MAF	Median MAF
<i>Ovis aries</i>	Namaqua Afrikaner	17	33631	0.15	0.18
<i>Ovis aries</i>	Ronderib Afrikaner	17	40602	0.21	0.22
<i>Ovis aries</i>	Afshari	37	45593	0.27	0.26
<i>Ovis aries</i>	Altamura	23	46656	0.28	0.27
<i>Ovis aries</i>	Awassi	2	30046	0.25	0.20
<i>Ovis aries</i>	Bangladeshi BGE	24	41429	0.23	0.23
<i>Ovis aries</i>	Barbados Blackbelly	24	43998	0.25	0.25
<i>Ovis aries</i>	Black Headed Mountain	24	44296	0.25	0.25
<i>Ovis aries</i>	Blackface	5	43365	0.30	0.25
<i>Ovis aries</i>	Blackface - Scottish	57	46713	0.28	0.27
<i>Ovis aries</i>	Boreray	19	37318	0.16	0.19
<i>Ovis aries</i>	Brazilian Creole	22	46600	0.30	0.27
<i>Ovis aries</i>	Bundner Oberlander	23	44476	0.26	0.26
<i>Ovis aries</i>	Castellana	22	46970	0.30	0.28
<i>Ovis aries</i>	Changthangi	29	46030	0.28	0.26
<i>Ovis aries</i>	Chios	23	43360	0.24	0.24
<i>Ovis aries</i>	Churra	118	47048	0.28	0.27
<i>Ovis aries</i>	Comisana	24	46747	0.29	0.27
<i>Ovis aries</i>	Composite	16	46316	0.28	0.27
<i>Ovis aries</i>	Cyprus Fat Tail	30	42322	0.23	0.23
<i>Ovis aries</i>	Deccani	24	43681	0.25	0.25
<i>Ovis aries</i>	Dorper	10	45584	0.25	0.24
<i>Ovis aries</i>	African Dorper	21	43151	0.26	0.25
<i>Ovis aries</i>	African White Dorper	6	37987	0.17	0.21
<i>Ovis aries</i>	Dorset	11	43531	0.27	0.26
<i>Ovis aries</i>	Poll Dorset	7	43664	0.29	0.24
<i>Ovis aries</i>	Dorset Horn	21	39029	0.19	0.21
<i>Ovis aries</i>	East Friesian Brown	39	41295	0.22	0.22
<i>Ovis aries</i>	East Friesian White	9	41800	0.22	0.22
<i>Ovis aries</i>	Engadine Red	24	46418	0.29	0.27
<i>Ovis aries</i>	Ethiopian Menz	34	42799	0.24	0.24
<i>Ovis aries</i>	Finnsheep	106	46566	0.28	0.27
<i>Ovis aries</i>	Galway	49	44443	0.26	0.25
<i>Ovis aries</i>	Garole - Bangladeshi	25	40591	0.22	0.23
<i>Ovis aries</i>	Garole - Indian	26	39012	0.21	0.21
<i>Ovis aries</i>	Garut	22	42529	0.25	0.24
<i>Ovis aries</i>	Gulf Coast Native	95	48143	0.30	0.29
<i>Ovis aries</i>	International Mapping Flock	27	45895	0.28	0.27
<i>Ovis aries</i>	Italian Sarda	5	43318	0.30	0.25

Table 2: BeadChip Content Validation (continued)

Species	Breed	Samples	Polymorphic Loci	Mean MAF	Median MAF
<i>Ovis aries</i>	Karakas	18	45205	0.25	0.25
<i>Ovis aries</i>	Katahdin	9	44661	0.22	0.24
<i>Ovis aries</i>	Lacaune	177	47840	0.28	0.27
<i>Ovis aries</i>	Leccese	22	46852	0.30	0.28
<i>Ovis aries</i>	Merino	7	46016	0.29	0.26
<i>Ovis aries</i>	Australian Merino	46	47726	0.29	0.28
<i>Ovis aries</i>	Chinese Merino	23	45938	0.28	0.27
<i>Ovis aries</i>	Macarthur Merino	10	31054	0.10	0.16
<i>Ovis aries</i>	Merinolandschaf	21	45531	0.29	0.27
<i>Ovis aries</i>	Moghani	34	46571	0.28	0.27
<i>Ovis aries</i>	Morada Nova	21	41324	0.24	0.23
<i>Ovis aries</i>	Navaho	1	19302	0.00	0.19
<i>Ovis aries</i>	Norduz	20	44694	0.25	0.25
<i>Ovis aries</i>	Ojalada	24	47412	0.29	0.28
<i>Ovis aries</i>	Old Norwegian Saelsau	15	45479	0.27	0.26
<i>Ovis aries</i>	Qezel	35	47048	0.29	0.27
<i>Ovis aries</i>	Rambouillet	112	46670	0.28	0.27
<i>Ovis aries</i>	Rasa aragonesa	20	48676	0.30	0.29
<i>Ovis aries</i>	Red Maasai	46	43513	0.24	0.24
<i>Ovis aries</i>	Romanov	8	40228	0.19	0.22
<i>Ovis aries</i>	Romney	6	44278	0.25	0.25
<i>Ovis aries</i>	New Zealand Romney	21	44962	0.29	0.26
<i>Ovis aries</i>	Sakiz	22	40451	0.23	0.23
<i>Ovis aries</i>	Santa Ines	46	45914	0.27	0.26
<i>Ovis aries</i>	Sardinian Ancestral Black	20	45451	0.25	0.25
<i>Ovis aries</i>	Soay	432	36728	0.17	0.19
<i>Ovis aries</i>	Spael-coloured	3	35619	0.17	0.21
<i>Ovis aries</i>	Spael-white	29	44293	0.26	0.25
<i>Ovis aries</i>	St Elizabeth	10	47441	0.30	0.27
<i>Ovis aries</i>	Suffolk	9	45002	0.28	0.24
<i>Ovis aries</i>	Irish Suffolk	55	43627	0.25	0.25
<i>Ovis aries</i>	American Suffolk	5	42012	0.20	0.24
<i>Ovis aries</i>	Sumatra	24	41703	0.23	0.23
<i>Ovis aries</i>	Swiss Black-Brown Mountain	22	44961	0.27	0.26
<i>Ovis aries</i>	Swiss Mirror	20	46122	0.28	0.26
<i>Ovis aries</i>	Swiss White Alpine	21	44986	0.26	0.26
<i>Ovis aries</i>	Texel	17	45155	0.26	0.26
<i>Ovis aries</i>	German Texel	43	45987	0.28	0.27
<i>Ovis aries</i>	New Zealand Texel	22	44348	0.27	0.26
<i>Ovis aries</i>	Scottish Texel	80	44521	0.25	0.25

Table 2: BeadChip Content Validation (continued)

Species	Breed	Samples	Polymorphic Loci	Mean MAF	Median MAF
<i>Ovis aries</i>	Indonesian Thin Tail	5	40026	0.20	0.22
<i>Ovis aries</i>	Sumatran Thin Tail	5	38363	0.20	0.21
<i>Ovis aries</i>	Tibetan	42	44578	0.26	0.25
<i>Ovis aries</i>	Valais Blacknose	23	41547	0.24	0.23
<i>Ovis aries</i>	Valais Red	21	40855	0.24	0.23
<i>Ovis aries</i>	Wiltshire	23	36460	0.17	0.19
<i>Ovis aries</i>	All Domestic sheep	2812	49451	0.30	0.28
Other <i>Ovis</i>	Wild sheep - All	116	41392	0.18	0.20
<i>Ovis musimon</i>	Mouflon - Sardinian	28	43053	0.23	0.23
<i>Ovis musimon</i>	Mouflon - European	24	29444	0.10	0.16
<i>Ovis vignei</i>	Wild Urial	4	25577	0.13	0.15
<i>Ovis ammon</i>	Wild Argali	6	17912	0.00	0.07
<i>Ovis dalli</i>	Thinhorn	2	503	0.00	0.00
<i>Ovis canadensis</i>	Bighorn	50	421	0.00	0.00
<i>Ovis canadensis</i>	Wyoming Bighorn	2	366	0.00	0.00
Outgroup	Outgroup	50	11994	0.00	0.06

Table 3: BeadChip Content Sources

Source	OvineSNP50 Probes
Validated SNPs (Sanger)**	557
Illumina RRS	17,042
mtDNA	8
Other	36,634
Total	54,241

\*\*Validated SNPs from a 1536-plex that the consortium did in advance of the SNP50.

Product Summary

Developed through a collaboration between Illumina scientists and leading ovine thought leaders, the OvineSNP50 BeadChip features more than 54,000 evenly spaced SNPs that provide comprehensive coverage of the ovine genome, enabling a diverse range of genetic research applications. As the first high-density whole-genome genotyping array for sheep, the 12-sample OvineSNP50 BeadChip presents a powerful, efficient, and cost-effective tool for a wide variety of genome-wide genetic analysis applications across the majority of sheep breeds worldwide

