

Ancient DNA Speaks Volumes About New World Populations and Migrations

What can an ancient human genome tell you? Quite a lot when you leverage the speed and sensitivity of Illumina's Genome Analyzer and the power of Infinium® BeadChips.

Introduction

"I lived several thousand years ago in what is now Greenland. I had thick black hair with a tendency towards baldness, brown eyes, brown skin, type A+ blood, shovel-formed front teeth, dry ear wax, a stocky build adapted to living in a cold climate, and I am more closely related to people living in Siberia than I am to anyone now living in North America or Greenland." ~ 4,000 year old Saqqaq individual

This was more information than Eske Willerslev, Ph.D. and Morten Rasmussen¹ expected to uncover when they decided to sequence the DNA of a ~ 4,000 year old tuft of permafrost-preserved human hair. Their goal at the onset was to sequence the first complete ancient human nuclear genome and in the process, gain information about one of the first cultures to settle in the New World.

The hair sample had been at Denmark's National Museum since it was discovered in 1986 at a site on the southwestern slope of Greenland. The location and radiocarbon dating placed the individual in the Saqqaq culture which populated Greenland from 2500 BCE until about 800 BCE².

The Benefits of Speed and Quality

Dr. Willerslev and Mr. Rasmussen knew from sequencing the mitochondrial DNA³ that the specimen was of reasonably good quality. While in ancient DNA terms, the sample looked promising, there would still be limitations and obstacles to overcome in sequencing a complete nuclear genome. "We were going from 16,000 base pairs for a mitochondrial genome to around 3 billion base pairs for a nuclear genome, said Dr. Willerslev. "We needed a system that delivered high output, with high quality. Once we verified that we could obtain a lot of useful sequence from a single run, it was clear the Illumina Genome Analyzer was the best choice. It made a resequencing project like this feasible."

It turned out that speed and output would be important for another reason. A team of German scientists was attempting to sequence another hominid ancient genome—that of a Neandertal⁴. To increase the chances of completing their research first, Dr. Willerslev approached the Beijing Genomics Institute (BGI) to enlist its help. BGI is one of the world's leading scientific organizations, operating seven genome research centers in mainland China for scientific collaboration and sequencing services. With help from BGI, the race was on to determine which team would be the first to sequence and publish an ancient human genome.

Accounting for Contamination and Age

With the sequencer and sequencing sites selected, the team turned to solving a critical issue that could impact the quality of the sequence—how to deal with the problem of contamination, both microbial and human. With keratins offering a protective coat around the internal DNA, any microbial contamination on the surface could be easily removed by bleaching and washing the hair before treatment with proteinase K.

"We were more concerned about human contamination, because we needed to know that we were actually sequencing an ancient human, not 90% modern contamination," said Mr. Rasmussen. "One of the beneficial aspects of Illumina's technology is its use of indexing adaptors and primers for library preparation." After the DNA extraction and library builds were performed in a clean laboratory in Copenhagen, the libraries could be indexed enabling the team to know exactly which sequences left the laboratory in Copenhagen. This made it easier to identify any contamination that might have entered the samples during analysis.



Dr. Eske Willerslev is leader of the Center of Excellence in GeoGenetics at the Natural History Museum of the University of Copenhagen, Morten Rasmussen is a Ph.D. candidate at the Center.

To account for DNA fragmentation, the DNA sample-to-adaptor volumes were adjusted for short insert DNA to approximate the 1:10 ratio suggested in Illumina's protocol, and the library protocol was modified to move the gel cutting step after library amplification to ensure the library would not be lost during purification. The amplification step was also adjusted to avoid misreads caused by deamination of cytosine, which accounts for 95% of ancient DNA damage and makes it difficult to distinguish between evolutionary derived substitutions and those caused by DNA damage. "As recommended in Illumina's protocols,

AGAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
TCAACGTACCGTAACGAAACGTATCAATTGAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
CGACGAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
AACGTACCATTAAAGAGCTACCGTGCAACGAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
AGATTGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
GATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT
CGTATCAATTGAGACTAAATATAACGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAAGAAATGATAACAGTAACACACTTCTGTAAACCTT

