


Improved Performance Using the Latest GenomeStudio Software Release


GenomeStudio software is continually optimized to deliver improved performance for the latest whole-genome arrays. Users in both high- and low-throughput labs should consider upgrading to the most recent GenomeStudio version for optimal performance. To demonstrate the improved performance time using the most recent version of the software, 1,000 samples were imported and processed using two versions of GenomeStudio software. As shown, the more recent version of GenomeStudio software (v2011.1) offers improved project creation times (Figure 7) and improved sample re-clustering times (Figure 8) over the preceding version.

Figure 7: Improved Project Creation Times with GenomeStudio v2011.1



Time required to create a project file from 1,000 samples scanned on four whole-genome arrays, processed using A) *.gtc files and B) *.idat files. The project creation time involves calculating sample and SNPs statistics. The data were processed using an 8 GB workstation with Windows 7 OS. For both file types, GenomeStudio v2011.1 delivered improved performance times across all arrays.

Figure 8: Improved Recluster Time with GenomeStudio v2011.1



Time to recluster 1000 samples on the Omni2.5 array using two recent versions of GenomeStudio. The data were process using 32 GB workstation with Windows 7 OS. GenomeStudio v2011.1 delivered significantly better performance time over the previous version of the software.

Summary

The microarray analysis pipeline can be configured at multiple points to optimize efficiency based on the volume of data being processed. By considering multiple factors such as the data file conversion, loci pre-filtering, and system hardware and software requirements, researchers can minimize the sample processing time for high- and low-throughput environments.