




Identifying Chromosomal Abnormalities Using Infinium SNP BeadChips




Harper VanSteenhouse, PhD
Product Manager, Diagnostics
Molecular Cytogenetics



Large & Growing Opportunity in Array-Based Cytogenetics



Array-based cytology market growing double digits (CAGR = 16–20%)




A technology substitution is underway

ILLUMINA IS POSITIONED FOR GROWTH

Building momentum: Arrays as first-line?

ISCA Consensus Statement: AJHG, May 2010



The World's Most Complete Cytogenetics Offering

Genome-Wide SNP BeadChips



- ▶ Multiple arrays for efficient identification and discovery of aberrant regions
- ▶ Flexible software options for any cytogeneticist (GenomeStudio or KaryoStudio)

Next-Gen Sequencing



- ▶ Sequencing to obtain the ultimate molecular karyotype and detection of any variant in any genome at single base resolution
- ▶ Ideal for discovering new aberrations across the genome with precise breakpoint identification


Plans for FDA Submission: iScan Platform

- ▶ Q4 '09 Illumina submitted Pre-IDE for Cytogenetics, followed by additional discussions with FDA
- ▶ Package includes iScan platform, arrays, reagents and software
- ▶ Intended use
 - *Post-natal*
 - *DD/MR/MCA*
 - *Broad, genome-wide coverage*

iScan/KaryoStudio




HumanCytoSNP-12



- Optimized for cytogenetics


HumanOmni1-Quad



- High resolution detection of micro-deletions
- Efficient four-sample format

Detection of a Wide Range of Aberrations

ILLUMINA DNA ANALYSIS AND SEQUENCING PRODUCTS IDENTIFY STRUCTURAL VARIATION




	POLYPLOIDY	ANEUPLOIDY	BALANCED ABERRATION	UNBALANCED ABERRATION	AMPLIFICATION	LOW-LEVEL MOSAICISM	COPY-NEUTRAL LOH
Infinium HD Assay BeadChips	+	+	-	+	+	+	+
Illumina Sequencing	+	+	+	+	+	+	+
Banding	+	+	+	+	-	+	-
FISH/SKY	+	+	+	+	+	+	-
Array-CGH	-	+	-	+	+	-	-

Adapted from Speicher and Carter, 2005

Getting High-Quality Cytogenetic Data is a SNP

- ▶ **Highest data quality**
 - Proven assay and analysis
 - High accuracy and signal-to-noise
 - Consistently high call rates
- ▶ **Most intelligent content**
 - BeadChips target high-value regions
 - And provide genome-wide coverage
- ▶ **Most informative marker panel**
 - High proportion of SNP markers
 - Optimized to identify the most relevant aberrations with high confidence
- ▶ **Streamlined sample to answer workflow**
 - Fast, simple, and proven protocols
 - Automated processing
 - Results visualization




Infinium HD Portfolio of Cytogenetic Products

Arrays

HumanCytoSNP-12

- ▶ Meticulously selected SNPs for whole-genome coverage enhanced in cytogenetic regions and genes
- ▶ Efficient detection of gains, losses, copy-neutral LOH, aneuploidy, and mosaicism




HumanOmni1-Quad

- ▶ High-resolution identification of aberrations and breakpoints across the whole genome
- ▶ Cutting-edge content in high-value genes and functional regions



Software

KaryoStudio™ v1.0



- ▶ Developed in collaboration with molecular cytogeneticists
- ▶ Automated cross-matching and reporting

GenomeStudio™ 2008



- ▶ Numerous options for copy number analysis with graphical, table, and data displays
- ▶ Links to 3rd party software offerings

Instruments

BeadArray Reader




iScan System




HumanCytoSNP-12 BeadChip: Optimized for Efficient Cytogenetic Analysis

- ▶ Streamlined, most informative set of targeted and whole-genome SNP and non-polymorphic markers
- ▶ Uniform picket fence of entire genome (including ~92% of RefSeq genes)
 - 300,000 markers, mostly SNPs
 - 6.2kb median marker spacing yields ~30kb resolution
- ▶ Higher density in cyto high-value regions (~250 for ~40% of genome)
 - All pericentromeres and subtelomeres
 - Sex chromosomes
 - Common regions of interest (e.g., associated with known syndromes)
 - Regions contain ~9000 genes
- ▶ Higher density in ~400 “disease genes”




HumanOmni1-Quad BeadChip: High-Density Backbone Plus High Value Regions

- ▶ Dense, uniform whole-genome backbone
 - Intelligent SNP selection from all 3 phases of HapMap
 - 1.2kb median marker spacing yields ~6kb resolution
 - > 98% of RefSeq genes
- ▶ Most valuable, cutting-edge content
 - Gene-centric SNPs, indels
 - > 600,000 markers within 10kb of a gene
 - 1000+ disease-associated SNPs
 - Selected SNPs from the first phase of the 1000 Genomes Project
- ▶ Rare and Common CNVs
 - discovery and characterization




The Infinium HD Assay




- ▶ Unconstrained Marker Design
 - Freedom to select the best, most informative SNPs then fill-in with intensity-only probes
- ▶ Well-proven
 - High reproducibility (> 99.9%)
 - High call rates (> 99%)
- ▶ Streamlined, automatable
- ▶ No need to run a reference sample
- ▶ High locus selectivity and allele specificity
 - Two-step enzymatic discrimination

SNPs Provide More Information to Detect Copy Number




SNPs Provide Higher Signal-to-Noise Ratios




Log R Ratio Only, More Markers



Same Markers, Including Genotypes



Abnormal NIGMS Cell Repository on HumanCytoSNP-12



Trisomy 21

22q11.2 Homozygous
Deletion

5q34-qter Gain

Example: Detecting Copy-Neutral Abnormalities

- ▶ Platform comparison performed by independent cytogenetics lab
 - Human1M-Duo and Agilent 244A arrays
 - Looking for cytogenetic abnormalities with focus on copy-neutral events
- ▶ 24 pre-characterized samples plus single unrelated male and single unrelated female reference samples
 - 400ng of input DNA
 - Standard quality filters for each platform
 - Illumina: SNP genotype call rate >99% per sample, copy number measurements $\log R$ std dev < 0.2
 - Agilent: DLR spread < 0.3
 - Manufacturer recommended algorithms

Copy-Neutral Cytogenetic Aberrations

Forms of copy-neutral cytogenetic aberrations

▶ Uniparental disomy

- Case in which individual receives two copies of a chromosomal region from one parent, none from the other

▶ Copy-neutral loss of heterozygosity (or “acquired uniparental disomy”)

- Case in which one allele of a gene in a heterozygote is already inactivated and the second, “good” allele is lost without a net change in copy number. This can occur through a gene conversion event in which the chromosome region containing the inactivated allele is used as a template to repair a gap occurring in the corresponding region of the other chromosome


In either case, the absence of a functional allele leaves the individual vulnerable to phenotypes that may be associated to the affected gene(s)

Copy-Neutral Cytogenetic Aberrations are Being Associated with a Rapidly Growing Number of Clinical Phenotypes

Pycnodysostosis	Am J Hum Genet. 1998 Apr;62(4):848-54
Herlitz junctional epidermolysis bullosa	J Invest Dermatol. 2000 Aug;115(2):307-11
Steroid 5alpha-reductase 2 deficiency	J Clin Endocrinol Metab. 2000 Sep;85(9):3147-50
Angelman syndrome	Hum Mol Genet. 1999 Jan;8(1):129-35.
Primary ciliary dyskinesia	Proc Natl Acad Sci U S A. 2002 Aug 6;99(16):10282-6. Epub 2002 Jul 25.
Late-infantile-onset neuronal ceroid lipofuscinosis	Hum Mutat. 2009 Mar 3
Pediatric acute lymphoblastic leukemia	Int J Oncol. 2009 Jun;34(6):1603-12
Prader-Willi syndrome	Endocrine. 2009 May 7
Myelodysplastic syndromes	Leukemia. 2009 Apr 23.
Silver-russell syndrome	Pediatrics. 2009 May;123(5):e929-31. Epub 2009 Apr 13
Angelman syndrome	Res Dev Disabil. 2009 Sep-Oct;30(5):1095-106.
Waldenstrom's macroglobulinemia.	Cancer Res. 2009 Apr 15;69(8):3579-88. Epub 2009 Apr 7
Autism Spectrum Disorders.	Biol Psychiatry. 2009 Mar 17.
Nephropathic cystinosis	Eur J Hum Genet. 2009 Mar 4
Silver-Russell syndrome (SRS)	J Med Genet. 2009 Mar;46(3):192-7. Epub 2008 Dec 9
Galactosaemia	J Inherit Metab Dis. 2008 Oct 29.
Myoclonus-dystonia	Arch Neurol. 2008 Oct;65(10):1380-5
Charcot-Marie-Tooth / Gaucher disease	Neurology. 2008 Mar 18;70(12):976-8
Glioblastomas	Neuro Oncol. 2008 Dec;10(6):995-1003. Epub 2008 Aug 12
Neuroblastoma	BMC Genomics. 2008 Jul 29;9:353
CLL	Blood. 2007 Feb 1;109(3):1202-10. Epub 2006 Oct 19
Head and neck squamous cell carcinoma	Head Neck. 2008 Oct;30(10):1361-83
Esophageal adenocarcinoma	Cancer Res. 2008 Jun 1;68(11):4163-72
MUTYH-associated polyposis carcinomas	J Pathol. 2008 Sep;216(1):25-31
Follicular lymphoma	Clin Cancer Res. 2007 Aug 15;13(16):4777-85

The Combination of B Allele Frequency and Log R Ratio Data Enables the Detection of CNLOH / UPD


- > 30 cases of CNLOH detected in sample cohort on Human1M-Duo BeadChip





Copy-Neutral Features Are Missed by Array CGH

Chromosome


3



Zoom to position




SNP-based Detection Provides More Information for Better Characterization of Chromosomal Aberrations



A profile of chromosome 3 of a cell line derived from a breast tumor.

Mosaicism Detection Benefits from SNP Markers



- Adapted from Nancarrow DJ, Handoko HY, Stark MS, Whiteman DC, Hayward NK (2007) [SiDCoN: a tool to aid scoring of DNA copy number changes in SNP chip data](#). PLoS ONE 2:e1093.

Mechanisms of mosaicism, chimerism and uniparental disomy identified by single nucleotide polymorphism array analysis

Laura K. Conlin^{1,2}, Brian D. Thiel¹, Carsten G. Bonnemann², Livija Medne⁴, Linda M. Ernst⁴, Elaine H. Zackai², Matthew A. Deardorff², Ian D. Krantz², Hakon Hakonarson^{2,3} and Nancy B. Spinner^{1,2,*}

¹Department of Pathology and Laboratory Medicine, ²Department of Pediatrics and ³Center for Applied Genomics, The Children's Hospital of Philadelphia and University of Pennsylvania School of Medicine, Philadelphia, PA 19104, USA and ⁴Department of Pathology, Feinberg School of Medicine, Northwestern University, Chicago, IL 60611, USA

Main conclusions:

- ▶ Frequency of mosaicism higher than previously reported
- ▶ Successful detection of low-level mosaicism
- ▶ Increased risk of UPD with mosaicism
- ▶ Interesting chimeric case

Mechanisms of Mosaicism, Chimerism, and UPD

Laura Conlin, et al. 2010

- ▶ Frequency of mosaicism higher than previously reported
 - ~1% of their total of >2000 samples (~10% of “aberrant” samples)
 - Much higher than reported previously, likely owing to the new abilities to detect

Table 1. Patients with mosaic aneuploidy or chimerism

Patient number	Type of aneuploidy	Mosaic % by array	Mosaic % by G-banded chromosomes	Mosaic % by metaphase FISH	Mosaic % by interphase FISH	Parental chromosome gained or lost	Mitosis/Meiosis	UPD	Tissue	Cell lines
1	Monosomy X	5	6.67				Mitosis	Yes	Blood	45,X/46,XX
2	Monosomy X	25	16.00			Paternal	Mitosis	Yes	Skin	45,X/46,X,r(Y)
3	Monosomy X	30	25.00		22.00	Paternal	Mitosis	Yes	Blood	45,X/46,XY
4	Monosomy X	30	40.00			Maternal	Mitosis	Yes	Blood	45,X/47,XXX
	Trisomy X	70	60.00			Paternal	Mitosis	No	Blood	45,X/47,XXX
5	Monosomy 7	40	0.00		7.00	Paternal	Mitosis	Yes	Blood	45,XY,-7/46,XY
6	Monosomy X	50	42.11				Mitosis	Yes	Blood	45,X/46,X,r(X)
7	Monosomy X	75	6.67				Mitosis	Yes	Blood	45,X/46,X,r(X)
8	Monosomy X	80	75.00				Mitosis	Yes	Blood	45,X/46,X,r(X)
9	Monosomy X	80	76.67				Mitosis	Yes	Blood	45,X/46,X,r(X)
10	Monosomy X	90					Mitosis	Yes	Blood	45,X/46,XX
11	Trisomy 14	5	0.00			Maternal	MI	No	Skin	47,XX, + 14/46,XX
		50	2.56			Maternal	MI	No	Blood	47,XX, + 14/46,XX
12	Trisomy 8	40	100.00				MI	Yes	Blood	47,XY, + 8/46,XY
13	Trisomy 9	50	0.00	2.50	24.00		MI	Yes	Blood	47,XX, + 9/46,XX
14	Trisomy 18	10	15.15				MII	No	Blood	47,XX, + 18/46,XX
15	Trisomy 14	20	10.00			Paternal	MII	Yes	Blood	47,XX, + 14/46,XX
16	Trisomy 8	5		1.50			Mitosis	No	Blood	47,XY, + 8/46,XY
17	Trisomy 9	20	2.00			Paternal	Mitosis	No	Blood	47,XY, + 9/47,XY
18	Trisomy 9	20	35.00				Mitosis	No	Blood	47,XX, + 9/46,XX
19	Trisomy 8	20			12.62		Mitosis	No	Blood	48,XY, + 8, + 19/46,XY
	Trisomy 19	20			14.00		Mitosis	No	Blood	48,XY, + 8, + 19/46,XY
20	Trisomy 21	50	85.00				Mitosis	No	Skin (hypo)	48,XX, + 7, + 21/46,XX
	Trisomy 7	50	85.00				Mitosis	No	Skin (hypo)	48,XX, + 7, + 21/46,XX
	Trisomy 21	60	50.00				Mitosis	No	Skin (hyper)	48,XX, + 7, + 21/46,XX
	Trisomy 7	60	50.00				Mitosis	No	Skin (hyper)	48,XX, + 7, + 21/46,XX
21	Trisomy 17	50					Mitosis	No	Skin (left)	47,XY, + 17/46,XY
		75					Mitosis	No	Skin (right)	47,XY, + 17/46,XY
30	Chimera	20	45.00			Paternal	Fertilization	Yes	Skin (hyper)	46,XX/46,XY


MI, meiosis I; MII, meiosis II.

Mechanisms of Mosaicism, Chimerism, and UPD

Laura Conlin, et al. 2010

► Successful detection of low-level mosaicism

- ~5% limit for mosaic detection. This is way better than anything reported before using aCGH (20% limit).
- The result is that they find many more mosaics than have been reported in aCGH papers (1% rather than 0.2%).




Mechanisms of Mosaicism, Chimerism, and UPD


Laura Conlin, et al. 2010

- ▶ Increased risk of UPD with mosaicism
 - They found 6 cases of UPD, where there was no consanguinity and >20Mb of homozygosity. (see Table 2 below)
 - UPD was present in 60% of cases meiotic mosaic trisomies (trisomy rescue).

Segmental UPD




D Mosaic UPD




Mechanisms of Mosaicism, Chimerism, and UPD

Laura Conlin, et al. 2010


► Parthenogenetic chimera case



Only one X chromosome haplotype



Chimeric isodisomy



The Industry's Best Cytogenetics SNP Arrays

	HumanCytoSNP -12	HumanOmni1 -Quad	300K Competitor	2.7M Competitor
Total Number of Markers for Copy Number Analysis	301,232	1,140,419	310,759	2,761,979
SNP Markers (% of Total)	294,748 (98%)	1,016,423 (90%)	16,666 (5%)	400,103 (14%)
Marker Spacing – Median	6.2kb	1.2kb	9.6kb	0.7kb
SNP Spacing – Mean	9.6kb	2.5kb	177kb	7.4kb
Feature Redundancy – Mean	12x		1x	
Hands on Time – Manual	< 2 hr per 48 samples (< 1hr Automated)		< 3hrs	
Protocol Duration	~40hr		~30hr	
Samples per Chip	12	4	1	1
Probe Length	50-mer		40-mer	
Hybridisation and Enzymatic Discrimination	Yes		No	
Assay Publications	>275		0	

iScan Reader

- ▶ Fast: 3–16 min per sample
- ▶ High resolution: 0.54 micron
- ▶ Dual line - lasers at 532 & 658nm
- ▶ Multi-application support
- ▶ Enhanced dynamic range and limit of detection
- ▶ Full integration with LIMS and Autoloader 2



KaryoStudio: Streamlined Cytogenetics Analysis



Found Regions Displayed in KaryoStudio

Found regions table

Index	Sample	Chr	Start	Stop	Length	Value	Conf	Comment	CNV Index
237	NA19193	12	61880	34659435	34597555	3	10286.52		2
							6576.188		3
							1539.1		1
							543.955		1
							526.9706		13
							329.3838		0
							316.6322		3
							213.1072		6
508	NA18862	7	7595687		97598	3	212.1275		0
289	NA18855	16	5001884		2328	3	203.6422		2
52	NA18501	15	18457255		6427	3	195.9704		3
454	NA19202	7	75981641		44	3	181.9529		1
235	NA19193	7	75981641		642045	3	146.0023		0
239	NA19193	14	105587447		516201	1	134.1853		4
309	NA18855	6	62026396		1472905	2	130.7134		23
479	NA19202	6	77337162		3323985	2	113.9422		27
542	NA18862	11	7412883		9068134	2	108.8162		34
342	NA19143	3	62398116		64521585	2	108.3492		8
121	NA19203	11	113625984		2162742	2	101.2007		32
186	NA19200	1	176112280		1967695	2	99.16006		9
416	NA18517	16			360258	3	97.75482		1
410	NA18862	11			2428073	2	95.14549		37

**Link out to databases:
UCSC, DGV, ENSEMBL, DECIPHER, etc.**


Samples table with QC score

Index	Name	ID	Plate	Well	% Defects
1	NA19144	NA19144	WG0093202-DNA	E02	1.36
2	NA18501	NA18501	WG0093202-DNA	C03	1.47
3	NA19203	NA19203	WG0093202-DNA	F02	1.05
4	NA18861	NA18861	WG0093202-DNA	D03	1.16
5	NA19200	NA19200	WG0093202-DNA	G02	1.51
6	NA19193	NA19193	WG0093202-DNA	E03	5.48
7	NA18855	NA18855	WG0093202-DNA	H02	1.12
8	NA19143	NA19143	WG0093202-DNA	F03	1.01
9	NA18505	NA18505	WG0093202-DNA	A03	0.93
10	NA18517	NA18517	WG0093202-DNA	G03	0.79
11	NA19202	NA19202	WG0093202-DNA	F03	1.42
12	NA18862	NA18862	WG0093202-DNA	F03	0.99

User-defined
Known Regions track

Found Regions

Karyotype View



Mouse-over region for details

Click identified region to zoom to in data plot

Report Generated From Aberrations

Sample and product information

Genome build and versions of software used


Information on found region

Gene information

Aberration display

Chromosome information

Cross-matches to user-defined Known Regions list



Recent Infinium Cytogenetics Publications

(www.illumina.com/publications)

- ▶ Conlin LK, Thiel BD, Bonnemann CG, et al. (2010) [Mechanisms of mosaicism, chimerism and uniparental disomy identified by single nucleotide polymorphism array analysis](#). Hum Mol Genet 19: 1263-75
- ▶ Descipio C, Morrissette JD, Conlin LK, et al. (2010) [Two siblings with alternate unbalanced recombinants derived from a large cryptic maternal pericentric inversion of chromosome 20](#). Am J Med Genet A 152A: 373-82
- ▶ Johnson DS, Gemelos G, Baner J, et al. (2010) [Preclinical validation of a microarray method for full molecular karyotyping of blastomeres in a 24-h protocol](#). Hum Reprod 25: 1066-75
- ▶ Shaikh TH, Gai X, Perin JC, Glessner JT, Xie H, et al. (2009) [High-resolution mapping and analysis of copy number variations in the human genome: A data resource for clinical and research applications](#). Genome Res Epub ahead of print.
- ▶ Beaujard MP, Chantot S, Dubois M, Keren B, Carpentier W, et al. (2009) [Atypical deletion of 22q11.2: Detection using the FISH TBX1 probe and molecular characterization with high-density SNP arrays](#). Eur J Med Genet Epub ahead of print.
- ▶ Gijsbers AC, Lew JY, Bosch CA, Schuurs-Hoeijmakers JH, van Haeringen A, et al. (2009) [A new diagnostic workflow for patients with mental retardation and/or multiple congenital abnormalities: test arrays first](#). Eur J Hum Genet Epub ahead of print.
- ▶ Jackson EM, Sievert AJ, Gai X, Hakonarson H, Judkins AR, et al. (2009) [Genomic analysis using high-density single nucleotide polymorphism-based oligonucleotide arrays and multiplex ligation-dependent probe amplification provides a comprehensive analysis of INI1/SMARCB1 in malignant rhabdoid tumors](#). Clin Cancer Res 15:1923-30.
- ▶ Kamath BM, Thiel BD, Gai X, Conlin LK, Munoz PS, et al. (2009) [SNP array mapping of chromosome 20p deletions: genotypes, phenotypes, and copy number variation](#). Hum Mutat 30:371-8.
- ▶ Fernández L, Nevado J, Santos F, Heine-Suñer D, Martínez-Glez V, et al. (2009) [A deletion and a duplication in distal 22q11.2 deletion syndrome region. Clinical implications and review](#). BMC Med Genet 10:48.
- ▶ Sievert AJ, Jackson EM, Gai X, Hakonarson H, Judkins AR, et al. (2009) [Duplication of 7q34 in pediatric low-grade gliomas detected by high-density single-nucleotide polymorphism-based microarray results in a novel BRAF mutation](#). J Clin Pathol 19:449-58.
- ▶ ... DA, Zhou L, Chan OT, et al. (2008) [Genomic copy number changes in archived ...](#)
- ▶ ... Bertherat J, Stratakis ... highly ...

Cytogenetic Analysis Solutions from Illumina

- ▶ **Infinium HD Assay:** Unconstrained marker design to select the best, highest value SNPs
- ▶ **HumanCytoSNP-12:** Optimized for efficient cytogenetic analysis
- ▶ **HumanOmni1-Quad:** The most comprehensive and cost-efficient high-density tool to find chromosomal aberrations
- ▶ **KaryoStudio:** Streamlined analysis workflow for rapid identification of cytogenetic aberrations
- ▶ **Sequencing:** The ultimate digital karyotype