

Welcome!

Transcriptome sequencing

Goodbye!

Arrays

- Pervasive transcription
- Antisense transcription
- Non-coding RNA
- Promoter associated RNA
- Micro-RNA/mRNA targets
- Enhancer RNAs
- Gene fusions


Abizar Lakdawalla, PhD
alakdawalla@illumina.com

Infinium, BeadXpress, VeraCode, IntelliHyb,
the property of their respective owners.

illumina



Human Genome *circa* 2010



hotly debated. HOTAIR stands out because it is a long piece of RNA that doesn't encode a protein but still does something biologically important¹. "HOTAIR was a gem in a sea [of long RNAs]," says John Rinn, a genome biologist who discovered the RNA while working at Stanford University in California. "It told us little about what the bulk of these things are doing. For that, we can't even see a common trend."


It is hard to comprehend the upheaval that RNA has been causing in molecular biology over the past few years. Once viewed as a passive

something about the sequences that is entic ular biologists. New s can stream out data m predecessors have ma cellular transcripts pe In 2008, this pro-


**"Many transcripts are made that we don't understand. We still don't know what those transcripts do, if anything."
— Ewan Birney**

of the genome of brea


RNA families



Genome Analyzer Transcriptional Analysis




Tag Sequencing




illumina®

Digital Tag mRNA Sample Preparation



Gene Expression



Veitch et al. *BMC Genomics* 2010, 11:124
<http://www.biomedcentral.com/1471-2164/11/124>



RESEARCH ARTICLE


Open Access


Digital gene expression analysis of two life cycle stages of the human-infective parasite, *Trypanosoma brucei gambiense* reveals differentially expressed clusters of co-regulated genes

Nicola J Veitch^{1†}, Paul CD Johnson^{2†}, Urmi Trivedi³, Sandra Terry¹, David Wildridge¹, Annette MacLeod^{1*}



RNA-Seq: Method





Green



Red



Overlay

X	1479
Y	7089
Green	1828
Red	1468

LETM1 | PCDH7 | CNGA1 | UGT2B15 | PTPN13 | CENPE | PDE5A | ANAPC10 | SPOCK3 | MTNR1A


position/search chr4:107,900,001-114,100,000

jump

clear

size 6,200,000 bp.


configure




mRNA-Seq Data is Information Rich

- ▶ mRNA Expression Profiling
- ▶ Alternative Splicing Analysis
- ▶ Analysis of expressed SNPs and mutations
- ▶ Analysis of Allelic-specific Expression
- ▶ Chimeric Transcript Discovery
- ▶ Gene Discovery and Annotation


normal KRT6A






S100A2



KLK6



tumor





Base Position	Cytosine: Band	Sequence(s)							
51165363	51165368	51167459	51169078	51170393	51171589	51172798	51174089	51175272	51176471
Sequence is too dense for display.									
KRT6A									

Base Position	Cytosine: Band	Sequence(s)							
151796793	151800683	151826273	151844463	151868323	151890943	151910133	151925223	151943313	
Sequence is too dense for display.									
S100A2									

Base Position	Cytosine: Band	Sequence(s)							
56160389	56162298	56164439	56167018	56169228	56171438	56163496	56165698	56168008	
Sequence is too dense for display.									
KLK6									

Tissue-specific genes



Read Counts for All RefSeq Genes

GenomeStudio - RNA Sequencing - maqc35


File Edit View Analysis Tools Window Help

Alleles Exon Table **Gene Table** Junction Table


Index	Chromosome	Start	End	Genes	brain35		UHR35	
					Counts	Counts (Raw)	Counts	Counts (Raw)
4366	12	55018930	55020461	IL23A	0.00	4	0.02	37
4367	12	55021649	55040176	STAT2	0.33	2889	0.30	2684
4368	12	55040622	55042850	APOF	0.00	0	0.02	68
4369	12	55097174	55129433	TIMELESS	0.03	279	0.62	5400
4370	12	55130151	55134696	MIP	0.00	7	0.00	2
4371	12	55148568	55151002	SPRYD4	0.06	229	0.15	617
4372	12	55151035	55168448	GLS2	0.36	1890	0.05	259
4373	12	55201993	55269875	RBMS2	0.02	76	0.03	124
4374	12	55275647	55316430	BAZ2A	0.25	4389	0.50	8963
4375	12	55318226	55326119	ATP5B	6.10	22495	11.03	40650
4376	12	55343392	55368345	PTGES3	0.63	2429	1.49	5700
4377	12	55392478	55405593	NACA	0.36	5053	1.05	14601
4378	12	55411631	55432413	PRIM1	0.04	113	0.28	824
4379	12	55443375	55467841	HSD17B6	0.09	269	0.02	66
4380	12	55603884	55614312	SDR-O	0.01	12	0.00	2
4381	12	55631483	55637685	RDH16	0.01	36	0.01	49
4382	12	55674622	55676736	GPR182	0.00	6	0.00	5
4383	12	55678885	55686497	ZBTB39	0.02	289	0.09	1167
4384	12	55690051	55696592	TAC3	0.08	145	0.04	80
4385	12	55708568	55730160	MYO1A	0.00	14	0.01	78
4386	12	55735693	55758810	TMEM194A	0.02	255	0.23	2452
4387	12	55768944	55775459	NAB2	0.53	2812	0.52	2761
4388	12	55775527	55791428	STAT6	0.23	1766	0.76	5966
4389	12	55808549	55893392	LRP1	1.37	40676	0.76	22668
4390	12	55896845	55906493	NXPH4	0.19	681	0.53	1893
4391	12	55909819	55914952	SHMT2	0.35	1460	2.78	11507
4392	12	55914982	55920742	NDUFA4L2	0.45	1043	0.52	1208
4393	12	55923509	55931236	STAC3	0.03	97	0.05	160
4394	12	55933815	55990513	R3HDM2	0.65	5166	0.45	3564
4395	12	56114810	56130876	INHBC	0.00	10	0.02	68

RNA-seq and Microarray Compared


- ▶ Expression levels are shown, as measured by RNA-Seq and tiling arrays for *Saccharomyces cerevisiae* cells. Agree for genes with medium levels of expression, but correlation is very low for genes with either low or high expression levels.



Performance of RNA sequencing is superior to other gene expression profiling methods



Dynamic Range



Accuracy and Sensitivity

illumina®

Overview of mRNA-Seq Assays

- ▶ Start with 1 ug of Total RNA
- ▶ Purify poly-A mRNA
- ▶ Randomly Fragment RNA

Random Priming


- ▶ Random Prime mRNA → cDNA
- ▶ Make 2nd Strand cDNA
- ▶ Ligate Sequencing Adapters

RNA Ligation

- ▶ CIP Treatment
- ▶ PNK Treatment
- ▶ Ligate 3'- Small RNA Adapters
- ▶ Ligate 5'- Small RNA Adapters
- ▶ RT to make cDNA




- ▶ PCR Amplify 15 Cycles
- ▶ Purify, Quantify
- ▶ Grow Clusters
- ▶ Sequence on Genome Analyzer
 - Single or Paired-end Reads
 - 35 to 100 base read lengths

Overlapping of PMPCB and DNAJC2 genes on two strands



Divergent Transcription from Active Promoters
 Amy C. Seila, *et al.*
Science 322, 1849 (2008);
 DOI: 10.1126/science.1162253


Nascent RNA Sequencing Reveals Widespread Pausing and Divergent Initiation at Human Promoters
 Leighton J. Core, *et al.*
Science 322, 1845 (2008);
 DOI: 10.1126/science.1162228




illumina®

CSDN Analysis

COLON-HOMO



ADIPOSE-HET



Differential Allelic Expression

Advances in the identification and analysis of allele-specific expression

Christopher G Bell and Stephan Beck

Address: Medical Genomics, University College London Cancer Institute, Huntley Street, London WC1E 6BT, UK.

Correspondence: Stephan Beck. Email: s.beck@ucl.ac.uk


OPEN ACCESS Freely available online



Transcriptome-Wide Identification of Novel Imprinted Genes in Neonatal Mouse Brain

Xu Wang¹, Qi Sun², Sean D. McGrath³, Elaine R. Mardis³, Paul D. Soloway⁴, Andrew G. Clark^{1*}

¹Department of Molecular Biology & Genetics, Cornell University, Ithaca, New York, United States of America, ²Computational Biology Service Unit, Life Sciences Core Laboratories Center, Cornell University, Ithaca, New York, United States of America, ³The Genome Center at Washington University, Washington University School of Medicine, St. Louis, Missouri, United States of America, ⁴Division of Nutritional Sciences, College of Agriculture and Life Sciences, Cornell University, Ithaca, New York, United States of America



Read Counts for All RefSeq Exons


GenomeStudio - RNA Sequencing - maqc35

File Edit View Analysis Tools Window Help


Alleles **Exon Table** Gene Table Junction Table

Index	Chromosome	Start	End	Genes	brain35		UHR35	
					Counts	Counts (Raw)	Counts	Counts (Raw)
228	1	1259749	1259787	DVL1	0.85	20	0.85	21
229	1	1260521	1261758	DVL1	2.85	3533	2.25	2782
230	1	1263220	1263426	DVL1	2.04	422	1.89	392
231	1	1263512	1263679	DVL1	3.50	588	3.03	509
232	1	1263765	1263896	DVL1	2.24	296	1.78	235
233	1	1264605	1264682	DVL1	1.90	148	1.58	123
234	1	1264825	1264892	DVL1	1.78	121	1.46	99
235	1	1264979	1265055	DVL1	1.75	135	1.74	134
236	1	1265281	1265420	DVL1	2.04	286	1.81	254
237	1	1265490	1265559	DVL1	1.80	126	1.46	102
238	1	1265653	1265746	DVL1	1.53	144	1.23	116
239	1	1266910	1267048	DVL1	2.22	308	2.13	296
240	1	1267296	1267399	DVL1	2.40	250	2.38	247
241	1	1267606	1267727	DVL1	2.43	296	1.89	231
242	1	1267932	1268001	DVL1	1.01	71	0.87	61
243	1	1274139	1274355	DVL1	1.85	401	1.71	370
244	1	1277934	1278872	MXRA8	0.43	407	0.53	498
245	1	1279091	1279171	MXRA8	0.48	39	0.68	55
246	1	1279273	1279349	MXRA8	0.65	50	0.77	59
247	1	1279436	1279475	MXRA8	0.18	7	0.20	8
248	1	1279597	1279752	MXRA8	0.60	94	0.78	121

Transcriptome isoforms –



Making Sense



- ▶ “... studying transcription for a long time and never seen this kind of transcription before ... but we also see a polymerase that appears to be pointing in the wrong direction.” - Seila/Sharp, MIT

“The gene, in other words, is in an identity crisis.” illumina®


Gene Fusions

LETTERS

Transcriptome sequencing to detect gene fusions in cancer

Christopher A. Maher^{1,3*}, Chandan Kumar-Sinha^{1,3*}, Xuhong Cao^{1,2}, Shanker Kalyana-Sundaram^{1,3}, Bo Han^{1,3}, Xiaojun Jing^{1,3}, Lee Sam^{1,3}, Terrence Barrette^{1,3}, Nallasivam Palanisamy^{1,3} & Arul M. Chinnaiyan^{1,2,3,4,5}



Chimera	Chimera Class	Location	5' Gene	3' Gene	Protein	Weight	Detection Methods
<i>BCR-ABL1</i>	Class I: Translocation	22q11.23	BCR, breakpoint cluster region	ABL1, Abelson tyrosine kinase	BCR, Ig heavy chain constant region	145 kDa	Long read, Short read, qRT-PCR, FISH
<i>MRPS10-HPR</i>	Class I: Translocation	6p21.1	MRPS10, mitochondrial ribosomal protein S10	HPR, histone H4 polycomb target 1	MRPS10, mitochondrial ribosomal protein S10	15 kDa	Long read, Short read, qRT-PCR, FISH
<i>MIPOL1-DGKB</i>	Class II: Inter-chromosomal complex	14q13.3-q21.1	MIPOL1, mirror-image polydactyly 1	DGKB, diacylglycerol kinase, beta	MIPOL1, mirror-image polydactyly 1	90 kDa	Long read, Short read, qRT-PCR, FISH
<i>TMPRSS2-ERG*</i>	Class III: Interstitial Deletion	21q22.3	TMPRSS2, transmembrane protease, serine 2	ERG, v-ets erythroblastosis virus E26 oncogene homolog (avian)	ERG, v-ets erythroblastosis virus E26 oncogene homolog (avian)	68 kDa	Long read, Short read, qRT-PCR, FISH
<i>USP10-ZDHHC7*</i>	Class III: Interstitial Deletion	16q24.1	USP10, ubiquitin specific peptidase 10	ZDHHC7, zinc finger, DHHC-type containing 7	ZDHHC7, zinc finger, DHHC-type containing 7	100 kDa	Long read, Short read, qRT-PCR, aCGH
<i>STRN4-GPSN2*</i>	Class IV: Intra-chromosomal complex	19q13.2	STRN4, striatin, calmodulin binding protein 4	GPSN2, glycoprotein, synaptic 2	GPSN2, glycoprotein, synaptic 2	100 kDa	Short read, qRT-PCR
<i>LMAN2-AP3S1</i>	Class IV: Intra-chromosomal complex	5q35.3	LMAN2 lectin, mannose-binding 2	AP3S1, adaptor protein 3 sigma 1	AP3S1, adaptor protein 3 sigma 1	100 kDa	Long read, Short read, qRT-PCR, FISH
<i>HJURP-EIF4E2*</i>	Class IV: Intra-chromosomal complex	2q37.1	HJURP, Holliday junction recognition protein	EIF4E2, eukaryotic translation initiation factor 4E2	HJURP, Holliday junction recognition protein	100 kDa	Long read, Short read, qRT-PCR, FISH
<i>INPP4A-HJURP*</i>	Class II: Intra-chromosomal complex	2q11.2	INPP4A, inositol polyphosphate-4-phosphatase, type 1	HJURP, Holliday junction recognition protein	INPP4A, inositol polyphosphate-4-phosphatase, type 1	100 kDa	Long read, Short read, qRT-PCR, FISH
<i>RC3H2-RGS3</i>	Class IV: Intra-chromosomal complex	9q34	RC3H2, ring finger and CCCH-type zinc finger domains 2	RGS3, regulator of G-protein signaling 3	RC3H2, ring finger and CCCH-type zinc finger domains 2	100 kDa	Long read, Short read, qRT-PCR, FISH
<i>ZNF649-ZNF577</i>	Class V: Read-through	19q13.33	ZNF649, zinc finger protein 649	ZNF577, zinc finger protein 577	ZNF649, zinc finger protein 649	100 kDa	Long read, Short read, qRT-PCR
<i>MBTPS2-YY2*</i>	Class V: Read-through	Xp22.1-p22.2	MBTPS2, membrane-bound transcription factor peptidase, site 2	YY2, yeast Y2 protein	MBTPS2, membrane-bound transcription factor peptidase, site 2	100 kDa	Long read, Short read, qRT-PCR
<i>C19ORF25-APC2</i>	Class V: Read-through	19p13.3	C19ORF25, chromosome 19 open reading frame 25	APC2, adenomatous polyposis coli 2	C19ORF25, chromosome 19 open reading frame 25	100 kDa	Long read, Short read, qRT-PCR
<i>WDR55-DND1</i>	Class V: Read-through	5q31.3	WDR55, WD repeat domain 55	DND1, dead end homolog 1 (zebrafish)	WDR55, WD repeat domain 55	100 kDa	Long read, Short read, qRT-PCR
<i>SLC45A3-ELK4*</i>	Class V: Read-through	1q32.1	SLC45A3, Solute carrier family 45 member 3	ELK4, ETS domain-containing protein	SLC45A3, Solute carrier family 45 member 3	100 kDa	Short read, qRT-PCR




Integrative analysis of the melanoma transcriptome

Michael F. Berger,¹ Joshua Z. Levin,¹ Krishna Vijayendran,^{1,2} Andrey Sivachenko,¹ Xian Adiconis,¹ Jared Maguire,¹ Laura A. Johnson,^{1,2} James Robinson,¹ Roel G. Verhaak,^{1,2} Carrie Sougnez,¹ Robert C. Onofrio,¹ Liuda Ziaugra,¹ Kristian Cibulskis,¹ Elisabeth Laine,³ Jordi Barretina,¹ Wendy Winckler,¹ David E. Fisher,^{4,5} Gad Getz,¹ Matthew Meyerson,^{1,2,6} David B. Jaffe,¹ Stacey B. Gabriel,¹ Eric S. Lander,^{1,7,8} Reinhard Dummer,³ Andreas Gnirke,¹ Chad Nusbaum,¹ and Levi A. Garraway^{1,2,6,9}

¹The Broad Institute of MIT and Harvard, Cambridge, Massachusetts 02142, USA; ²Department of Medical Oncology, Dana-Farber



Foxa2 Binds in Unknown Regions



... two remote peaks in regions containing no known genes. The first peak (unknown #4) lies close to a mouse EST, while the second peak (unknown #5) overlaps and is close to highly conserved sequence regions.

Chromatin signature reveals over a thousand highly conserved large non-coding RNAs in mammals

Mitchell Guttman^{1,2}, Ido Amit¹, Manuel Garber¹, Courtney French¹, Michael F. Lin¹, David Feldser³, Maite Huarte^{1,6}, Or Zuk¹, Bryce W. Carey^{2,8}, John P. Cassidy^{2,8}, Moran N. Cabili⁷, Rudolf Jaenisch^{2,8}, Tarjei S. Mikkelsen^{1,4}, Tyler Jacks^{2,3}, Nir Hacohen^{1,9}, Bradley E. Bernstein^{1,10,11}, Manolis Kellis^{1,5}, Aviv Regev^{1,2}, John L. Rinn^{1,6,11*} & Eric S. Lander^{1,2,7,8*}





Figure 1 | Intergenic K4-K36 domains produce multi-exonic RNAs.

More to small RNA than known microRNAs



More to small RNAs ...

Hairpin RNA induces secondary small interfering RNA synthesis and silencing in *trans* in fission yeast

Femke Simmer^{1*}†, Alessia Buscaino^{1*}, Isabelle C. Kos-Braun^{1‡}, Alexander Kagansky¹, Abdelhalim Boukaba^{1§}, Takeshi Urano², Alastair R.W. Kerr¹ & Robin C. Allshire¹⁺

Cross-mapping and the identification of editing sites in mature microRNAs in high-throughput sequencing libraries

Michiel J.L. de Hoon,¹ Ryan J. Taft,² Takehiro Hashimoto,¹ Mutsumi Kanamori-Katayama,¹ Hideya Kawaji,¹ Mitsuoki Kawano,¹ Mami Kishima,¹ Timo Lassmann,¹ Geoffrey J. Faulkner,³ John S. Mattick,² Carsten O. Daub,¹ Piero Carninci,¹ Jun Kawai,¹ Harukazu Suzuki,¹ and Yoshihide Hayashizaki^{1,4}

Tiny RNAs associated with transcription start sites in animals



Ryan J Taft¹, Evgeny A Glazov², Nicole Cloonan¹, Cas Simons¹, Stuart Stephen¹, Geoffrey J Faulkner¹, Timo Lassmann³, Alistair R R Forrest^{3,4}, Sean M Grimmond¹, Kate Schroder¹, Katharine Irvine¹, Takahiro Arakawa³, Mari Nakamura³, Atsutaka Kubosaki³, Kengo Hayashida³, Chika Kawazu³, Mitsuyoshi Murata³, Hiromi Nishiyori³, Shiro Fukuda³, Jun Kawai³, Carsten O Daub³, David A Hume^{1,5}, Harukazu Suzuki³, Valerio Orlando^{6,7}, Piero Carninci³, Yoshihide Hayashizaki³ & John S Mattick¹

Post-transcriptional regulation of miR-27 in murine cytomegalovirus infection

AMY H. BUCK,^{1,2} JONATHAN PEROT,³ MICHAEL A. CHISHOLM,^{1,2} DIWAKAR S. KUMAR,^{1,2} LEE TUDDENHAM,³ VALÉRIE COGNAT,⁴ LISA MARCINOWSKI,⁵ LARS DÖLKEN,⁵ and SÉBASTIEN PFEFFER³

illumina®

Identification of micro-RNA targets





Henderson, 2009

illumina®

Argonaute HITS-CLIP decodes microRNA-mRNA interaction maps. Chi SW, Zang JB, Mele A, Darnell RB. *Nature* (2009) 460: 479-86.

Covalently crosslink native argonaute protein-RNA gives you Ago-miRNA and Ago-mRNA binding sites.



Ribo-Seq. RNA engaged in translation



Genome-Wide Analysis In Vivo of Translation with Nucleotide Resolution Using Ribosome Profiling

Nicholas T. Ingolia,^{1*} Sina Ghaemmaghani,^{1,2} John R. S. Newman,¹ Jonathan S. Weissman¹

¹Howard Hughes Medical Institute, Department of Cellular and Molecular Pharmacology, University of California, San Francisco, and California Institute for Quantitative Biosciences, San Francisco, CA 94158, USA. ²Present address: Institute for Neurodegenerative Diseases, University of California, San Francisco, San Francisco, CA 94158, USA.

Ribo-Seq

RNA-Seq







RNA Normalization

illumina®

Genome Analyzer Transcriptional Analysis



C₀t Curves and cDNA Normalization



illumina

Double strand nuclease normalization

mRNA-Seq

Purified Total RNA

Poly-A Selection


RNA Fragmentation

cDNA Synthesis*

Adapter Ligation & PCR

Total RNA-Seq

Purified Total RNA



RNA Fragmentation

cDNA Synthesis*


Adapter Ligation & PCR

DSN Normalization

illumina®

DSN Library Normalization Process Overview

- Library is melted at 98°C
- Library is hybridized at 68°C
- DSN Cleavage
- PCR to enrich library
- Grow Clusters and Sequence




From Evrogen.com

Total RNA-Seq Method Overview


- ▶ **RNA from any organism – bacteria, plants, animals**
- ▶ **Less than 100 nanograms of total RNA**
- ▶ **Full-length cDNA coverage of *all* RNA molecules**
- ▶ **Works with low quality RNA – even FFPE preparations**
- ▶ **Normalization protocol can be used with most RNA-Seq protocols**

Comparison of Output from Different Protocols






na^o


Many ncRNAs are often transcribed from the same strand as mRNA



Many non-coding RNAs are seen with *Total* RNA-Seq




DSN Normalization has little effect on most mRNAs



illumina®


Many non-coding RNAs are observed with *Total* RNA-Seq Method

Total RNA vs. Poly-A mRNA



Plot shows gene level count data
for all RefSeq Genes

Fold Change of Total RNA-Seq vs. qRT-PCR



Key Variables:

- RNA is degraded
- No Poly-A Selection
- + DSN Normalization
- 0.1 μ g input total RNA


illumina®



RNA-Seq in FFPE


illumina®

RNA-Seq Protocol Variations




Total RNA-Seq protocol on FFPE Samples

- ▶ Two Cancer Samples
 - Her2 and Basal breast cancer sub-types
- ▶ Two types of RNA preparations
 - Total RNA from fresh/frozen tissue and highly degraded FFPE RNA
- ▶ RIN of the fresh samples are 9.3 and 8.3




Zoom in
➔



illumina®

Collaboration with Chuck Perou - UNC

Pair-wise Comparison of FFPE



FFPE *Total* RNA-Seq vs. Fresh *Total* RNA-Seq

Whole Transcriptome Studies on HiSeq 2000 System

INDIVIDUAL HUMAN TISSUES


16 tissues, 1 lane per tissue

Standard Poly-A mRNA-Seq Library Preps

TOTAL HUMAN TRANSCRIPTOME

Equal Mixture of 16 Human Tissues

Strand-specific mRNA-Seq



Body Map 2.0 Project

Adrenal, Adipose, Brain, Breast, Colon, Heart, Kidney, Liver, Lung, Lymph Node, Ovary, Prostate, Skeletal Muscle, Testis, Thyroid, White Blood Cells


SUMMARY OF RUN METRICS

	Format	PF reads	Average Error Rates (PhiX spike-ins)	Yield
Individual human tissues	2 X 50	2.55 B	0.22%, 0.44%	>120 Gb in <5 days
Individual human tissues	1 X 75	1.27 B	0.37%, 0.52%	>95 Gb in <4 days
Total human transcriptome	1 X 100	1.19 B	0.67%, 0.81%	>115 Gb in <5 days




PCR-free





GAPDH


12




GAPDH

Gene ID | Gene Details | Selected Probe Info | CpG Island | Available Data Tracks


HCT10735 - SE (1pxl=80B)




HCT10736 - SE (1pxl=80B)



HCT10737 - SE (1pxl=80B)



HCT10738 - SE (1pxl=80B)



2




Base Position	21,073,095	21,079,555	21,086,015	21,092,475	21,098,935	21,105,395	21,111,855	21,118,315	21,124,775	21,130,000
Cytogenetic Band	p24.1									
Sequence (+)	No sequence data file found for this chromosome.									
HCT10735										

Gene ID Gene Details Selected Probe Info CpG Island Available Data Tracks

No PCR vs PCR (15 cycles)

HCT10736 Counts (Raw) vs HCT10738 Counts (Raw)





DATA

ANALYSES

Samples/flow cell and cost/sample (all consum.)

	GAlle	GAlx	HiScanSQ	HiSeq 2000
DGE (25bp, 2 million)	60 € 230	125 € 205	125 € 204	250 € 192
RNA-Seq (75x2, 20 million)	6 € 1585	13 € 898	13 € 835	25 € 522
ChIP-Seq (TF)	125 € 252	250 € 231	250 € 229	500 € 219
Virus/BAC clones (100x2, 30x)	2,770 € 209	5,555 € 207	5,555 € 207	11,111 € 206
Bact./GWAS (100x2, 30x)	277 € 242	556 € 224	556 € 222	1,111 € 214
Hu exome (100x2, 30x)	28 € 568	56 € 387	56 € 368	111 € 287
Hu genome (100x2, 30x)	0.28 € 36,444	0.56 € 18,325	0.56 € 16,430	€ 5,616 @300Gb/run

a

Estrogen receptor alpha controls a gene network in luminal-like breast cancer cells comprising multiple transcription factors and microRNAs. Cicatiello L ...
Am J Pathol (2010) 176: 2113-30.

- ▶ 1270 genes respond to beta-estradiol, includes several direct targets of Eralpha, organized in a gene regulation cascade, stemming from ligand-activated receptor and reaching a large number of downstream targets via AP-2gamma, B-cell activating transcription factor, E2F1 and 2, E74-like factor 3, GTF2IRD1, hairy and enhancer of split homologue-1, MYB, SMAD3, RARalpha, and RXRalpha transcription factors.
- ▶ MicroRNAs integral components of this gene regulation network; miR-107, miR-424, miR-570, miR-618, and miR-760 are regulated by 17beta-estradiol along with other microRNAs that can target a significant number of transcripts belonging to one or more estrogen-responsive gene clusters.

Reduced Representation Bisulfite Seq

Genomic DNA



Genome-scale DNA methylation maps of pluripotent and differentiated cells

Alexander Meissner^{1,2,3*}, Tarjei S. Mikkelsen^{1,4*}, Hongchang Gu², Marius Wenig¹, Jacob Hanna¹, Andrey Sivachenko², Xiaolan Zhang¹, Bradley E. Bernstein^{2,5,6}, Chad Nusbaum², David B. Jaffe¹, Andreas Gnirke¹, Rudolf Jaenisch^{1,7} & Eric S. Lander^{2,3,7,8}

BC Cell Line RRBS Data - MMEL1 Gene



BC Cell Line RRBS Data - CD247 Gene

MCF7

ZR-75-1

BT474


MDA-MB-231

MDA-MB-468

T47D

BT20

MCF10A



BC Cell Line RRBS Data - CDH13 Gene

MCF7

ZR-75-1

BT474

MDA-MB-231

MDA-MB-468

T47D

BT20

MCF10A




Correlation of RRBS and mRNA-Seq Data for PRKD1 Gene

chr14:29,027,620-29,554,515 PRKD1 gene methylation for 5' CpG island region been dramatically changed and anti-correlates well with mRNA levels




CDH13 (H-cadherin): increased methylation in lung tumor

Lung





Breast




Study both individual CpG sites and CpG Islands

illumina




Comparison of Methylation Summarized for ~ 25K CpG Islands



illumina®

Reduced Representation Bisulfite Sequencing




Whole Genome Bisulfite Sequencing



illumina®

CpG Island Methylation Changes & Pathways



What are in these pathways that change?

Top 3 Pathways by Enrichment

cAMP-mediated Signaling
($p = 2.09e-8$)

G-Protein Coupled Receptor Signaling
($p = 4.17e-7$)

Axonal Guidance Signaling
($p = 1.18e-6$)

**Cell Death
Related
Pathways**

Implication: increased methylation of cell death related pathways could promote cell growth and cancer progression

Experimental Data Available

- RNA-Seq
 - 50 bp paired-end standard mRNA-Seq
 - 100 bp single-read directional mRNA-Seq
- DNA Methylation - RRBS
- CNV - low pass genomic DNA sequencing
- Small RNA Analysis