





**Table 3: TruSight Tumor 15 Coverage**

Sample ID	Sample Quality	% of Bases ≥ 500x	Amplicon Mean Coverage
FFPE_Colon1	Medium	99.7%	24,219x
FFPE_Colon2	Low	99.9%	20,763x
FFPE_Colon3	Low	99.2%	35,270x
FFPE_Colon4	High	100.0%	18,357x
FFPE_Colon5	High	100.0%	15,769x
FFPE_Melanoma1	Medium	99.7%	32,707x
FFPE_Melanoma2	Low	99.1%	41,640x
FFPE_Melanoma3	High	100.0%	17,285x
FFPE_Melanoma4	Low	95.7%	10,177x
FFPE_Breast1	High	99.1%	15,501x

DNA was extracted from FFPE tumor samples and then 20 ng of input DNA was evaluated using the TruSight Tumor 15 assay and sequenced on the MiniSeq System. Coverage of ≥ 500x is required for accurate identification of mutations at 5% variant frequency. Sample quality was determined by amplification potential of extracted DNA compared to a non-FFPE control sample in a qPCR assay. High quality is indicated by ΔCq of 0–2. Medium quality is indicated by ΔCq value of 2–4. Low quality is indicated by ΔCq of 4–6.

**Table 4: TruSight Tumor 15 Performance with Characterized Horizon Sample**

Gene	Mutation	Reported Frequency	Detected Frequency	Coverage
<i>BRAF</i>	V600E	10.5%	12.3%	55,457x
<i>KIT</i>	D816V	10.0%	10.3%	5463x
<i>EGFR</i>	ΔE746-A750	2.0%	2.1%	3553x
<i>EGFR</i>	L858R	3.0%	4.1%	1761x
<i>EGFR</i>	T790M	1.0%	1.2%	18,927x
<i>EGFR</i>	G719S	24.5%	25.6%	41,805x
<i>KRAS</i>	G13D	15.0%	15.3%	6745x
<i>KRAS</i>	G12D	6.0%	7.2%	6742x
<i>NRAS</i>	Q61K	12.5%	11.2%	13,154x
<i>PIK3CA</i>	H1047R	17.5%	18.8%	21,522x
<i>PIK3CA</i>	E545K	9.0%	7.8%	13,250x

DNA from the HD-C749 formalin-fixed cell line (Horizon Diagnostics) containing known variants was evaluated using the TruSight Tumor 15 assay and sequenced on the MiniSeq System. Variants were analyzed using VariantStudio. HD-C749 showed 100% concordance over 7 different runs.

**Table 5: TruSight Tumor 15 Performance with FFPE Tumor Samples**

Sample	Reported Mutation	Detected Mutation	Detected Frequency	Coverage
FFPE_Colon1	<i>KRAS</i> G12S	<i>KRAS</i> G12S	22.3%	21,134x
FFPE_Colon2	<i>KRAS</i> G12D	<i>KRAS</i> G12D	11.5%	4322x
FFPE_Colon3	<i>BRAF</i> V600E	<i>BRAF</i> V600E	25.5%	140,040x
FFPE_Colon4	<i>KRAS</i> G12V	<i>KRAS</i> G12V	33.4%	5256x
FFPE_Colon5	<i>KRAS</i> G13D	<i>KRAS</i> G13D	33.0%	4156x
FFPE_Melanoma1	<i>BRAF</i> V600E	<i>BRAF</i> V600E	65.7%	106,924x
FFPE_Melanoma2	<i>KRAS</i> G12R	<i>KRAS</i> G12R	4.1%	54,622x
FFPE_Melanoma3	<i>BRAF</i> V600E	<i>BRAF</i> V600E	93.5%	61,838x
FFPE_Melanoma4	<i>BRAF</i> V600K	<i>BRAF</i> V600K	22.2%	8075x
FFPE_Breast1	<i>AKT1</i> E17K	<i>AKT1</i> E17K	37.3%	56,438x

DNA from FFPE tumor samples was extracted and then evaluated using the TruSight Tumor 15 assay and sequenced on the MiniSeq System. Variants were analyzed using VariantStudio. All 10 FFPE samples had 100% variant concordance.

AAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA/ AATCAACGTACCGTAACGAACGTATCAATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA/ AACGACGAAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA/ TTAAGGTACCATTAAGAGCTACCGTGCAAACGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTGCAAACGACGAAAGAATGATA/ AAAAGATTGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA/ AAGATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTGCAAACGAAAGAATGATAAC/ AACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTGCAAACGAAAGAATGATAAC/

## Reliable, Accurate Performance

TruSight Tumor 15 provides the analytical sensitivity and accuracy needed to identify low-frequency variation with confidence in samples of varying quality. High target coverage (at least 93.5% of bases covered at ≥ 500x) provides analytical sensitivity and accuracy required for low-level variant calling (Table 3). TruSight Tumor 15 run on the MiniSeq System enables variant detection in many different sample types, with detection as low as 1% in high-quality DNA (Table 4), and as low as 5% in low quality FFPE samples (Table 5).

## Summary

TruSight Tumor 15 offers a comprehensive workflow solution for the detection of the most common somatic variants found in solid tumors. Developed according to evidence-based guidelines, with input from key opinion leaders, and late-stage pharmaceutical research, the panel enables labs to use the power of NGS technology to focus on the most relevant genes and analyze low-frequency variants from FFPE DNA with confidence. By assessing 15 genes in a single assay, this panel offers a streamlined, economical workflow solution that can be implemented easily by labs using NGS for the first time.

## Ordering Information

Product	Catalog No.
TruSight Tumor 15 Includes library preparation consumables, oligos, and indexes sufficient for 24 samples	OP-101-1002
TruSight Tumor 15 MiSeq Kit Includes library preparation panel, 3 MiSeq v3 Kits, sufficient for 24 samples	OP-101-1001
TruSight Tumor 15 MiniSeq Kit Includes library preparation panel, 3 MiniSeq High Output Kits (300 cycles), sufficient for 24 samples	20005610

## Learn More

For more information about Illumina technology for oncology applications, visit [www.illumina.com/oncology](http://www.illumina.com/oncology).

## References

- Wagle N, Berger MF, Davis MJ, et al. High-throughput detection of actionable genomic alterations in clinical tumor samples by targeted, massively parallel sequencing. *Cancer Discov.* 2012;2:82-93.
- TruSight Tumor 15 Kit Support ([support.illumina.com/sequencing/sequencing\\_kits/trusight-tumor-15-kit/downloads.html](http://support.illumina.com/sequencing/sequencing_kits/trusight-tumor-15-kit/downloads.html)) Accessed 15 February 2016.
- Van Cutsem E, Cervantes A, Nordlinger B, Arnold D; ESMO Guidelines Working Group. Metastatic colorectal cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up. *Ann Oncol.* 2014;25 Suppl 3:iii1-9.
- Voskoboynik M, Mar V, Mailer S, et al. Clinico-pathological Characteristics Associated With BRAFK and BRAFL 597 Mutations in Melanoma. *Pigment Cell Melanoma Res.* 2015. doi: 10.1111/pcmr.12450.
- Kopetz S, Desai J, Chan E, et al. Phase II Pilot Study of Vemurafenib in Patients With Metastatic BRAF-Mutated Colorectal Cancer. *J Clin Oncol.* 2015;33(34):4032-4038.
- Yan L, Zou L, Zhao W, et al. Clinicopathological significance of c-KIT mutation in gastrointestinal stromal tumors: a systematic review and metaanalysis. *Sci Rep.* 2015. doi: 10.1038/srep13718.
- Salipante SJ, Kawashima T, Rosenthal C, et al. Performance and comparison of Illumina and Ion Torrent next-generation sequencing platforms for 16S rRNA-based bacterial community profiling. *Appl Environ Microbiol.* 2014;80:7583-7591.

