

Reducing Residual Risk in CF Carrier Screening - Using the illumina MiSeqDx™ for Cystic Fibrosis Carrier Screening

M. Chelliserry, R.E. Lenta, J.C. San, W.M. Goldstein, and E.J. Peters,
Illumina Inc., 5200 Illumina Way, San Diego, CA 92122 USA


Introduction

Cystic Fibrosis (CF) is one of the most common genetic diseases in the US, with approximately 1 in 29 Caucasian individuals being a carrier for causative CF mutations. While CF is most common in individuals of Caucasian descent, ACOG/ACMG has recommended that carrier screening be offered to all couples who are considering having a child. In order to increase the detection rate of cystic fibrosis carriers and reduce the residual risk among other demographic groups, genotyping-based tests have been developed with larger mutation panels beyond the 23 ACOG/ACMG-recommended mutations. Despite this effort, non-Caucasian groups still have risk mutations that current genotyping panels will not detect. The CFTR2 project (cftr2.com) has systematically examined reported CF-causing mutations and their functional effects. Here we describe the development of a new next generation sequencing based assay, the Illumina MiSeqDx Cystic Fibrosis Carrier Screening Assay, for the detection of CFTR mutations on the MiSeqDx instrument. The test is designed to detect the 162 CFTR mutations listed in the CFTR2 database and is intended to identify an individual's CF carrier status in genomic DNA extracted from whole blood. The results of the test are intended for interpretation by a certified clinical geneticist or equivalent. This test is not indicated for fetal diagnostic testing, for preimplantation testing or for stand-alone diagnostic purposes.

Assay Technology

The assay technology involves targeted amplification of the CFTR gene followed by automated sequencing by synthesis on the MiSeqDx instrument. TruSeq Custom Amplicon workflow is used for targeted amplification of CFTR gene. The process involves hybridization of CF oligos to unfragmented genomic DNA followed by extension and ligation to form DNA templates containing regions of interest flanked by universal primer sequences. Using indexed primers supplied with the kit, the DNA templates are then PCR amplified, pooled into a single tube and sequenced on the MiSeqDx system.

Library Preparation method for MiSeqDx Cystic Fibrosis Assay



- A Hybridization of custom oligonucleotide probes
- B Extension and ligation
- C Addition of indices and sequencing adapters by PCR
- D Final amplicon ready for sequencing with MiSeq

Assay Workflow

The MiSeqDx Cystic Fibrosis Carrier Screening Assay allows up to 48 samples to be processed in less than 48 hours from extracted DNA through completed data analysis. For the library preparation, up to 96 samples can be processed from extracted DNA to normalized samples (ready to be loaded on the sequencing instrument) within 7 hrs with less than 2.5 hrs of hands on time. Up to 48 samples can be pooled and genotyped/sequenced in a single MiSeqDx run.


Sequencing on the MiSeqDx instrument

- ▶ The Illumina MiSeqDx system is a bench top personal sequencer which utilizes Sequencing by Synthesis (SBS) technology. The MiSeqDx has an integrated fluidics architecture and a built-in CPU, which enables cluster generation, sequencing, data analysis, and mutation report generation to be performed on a single instrument.
- ▶ After library preparation, the pooled, normalized and indexed library is loaded on to the MiSeq reagent cartridge which contains all of the reagents required for cluster generation and SBS. The library is first hybridized, then covalently attached through bridge amplification onto the flow cell surface, and amplified to generate millions of clusters which can then be sequenced using SBS.
- ▶ The SBS process uses four fluorescently labeled nucleotides; during each sequencing cycle, a single dNTP with reversible terminator nucleotide.
- ▶ The imaging subcomponent of the instrument consists of two cameras and two LEDs. Each LED is able to capture fluorescence in two channels (530 nm and 660 nm), which together allow for the system to recognize the four base pairs.
- ▶ The Illumina MiSeqDx Cystic Fibrosis System performs paired end 2 x 150 cycle sequencing to allow sequencing for 150 cycles from both directions, and 2x8 cycle sequencing to determine the sequences of both indexes in each cluster.



Sample Indexing/Multiplexing

The PCR primers include index sequences for sample multiplexing. There are 8 unique i5 primers and 12 unique i7 primers. The MiSeqDx Cystic Fibrosis Carrier Screening assay allows up to 48 samples to be multiplexed on a single sequencing flow cell.



i7 Index PCR Primer	Index Sequence	A i5 primers (white caps)	B i7 primers (orange caps)	C AMP plate
A701	ATCAAGAC	AS1	TGACCTT	
A702	ACAGCTGT	AS2	TGCTAGT	
A703	CAGATCCA	AS3	TGTTTCT	
A704	ACAAAGCG	AS4	TGATGCT	
A705	ACCCAGCA	AS5	TGATGCT	
A706	AACCCCTC	AS6	TGATGCT	
A707	CCCAACTC	AS7	TGATGCT	
A708	CACACAC	AS8	TGATGCT	
A709	GAACCCCA	AS9	CTAGTCA	
A710	TGTGACCA	AS10	CTAGTCA	
A711	AGGCTCAA	AS11	TAGGTTC	
A712	AGGACTGG	AS12	TAGACTA	

To demonstrate that sample indexing has no impact on the sequencing result, 470 unique samples were tested across 2 different indices, 15 unique samples across 3 indices, and 3 unique samples across 47 indices. In each the results were 100% reproducible across all indices tested.

Assay Performance

The Illumina MiSeqDx Cystic Fibrosis Carrier Screening Assay had an average call rate of 99.99% when tested on N=400 unrelated blood samples (393 samples had a 100% call rate, 7 samples had a call rate of 99.35% due to No call for PolyTG) and a call rate >99.9% when tested on >1500 HapMap and human variation samples of multiple ethnicities. The performance of the assay was verified for all variants in the panel; for the rare mutations for which no samples were available, the performance was confirmed by using synthetic, plasmid-based DNA samples.

Accuracy and Reproducibility

A set of 47 Coriell DNA samples, chosen to be representative of the different types of mutations within CFTR gene, were tested using the Illumina MiSeqDx Cystic Fibrosis Carrier Screening assay by 3 operators on each of 3 MiSeqDx instruments. The results across 9 MiSeqDx runs indicated excellent reproducibility (100%) and accuracy (100%) when compared to results obtained with bi-directional Sanger sequencing (for all variations except the 2 large deletions) and PCR assay (for the 2 large deletions).

	#Correct Calls	#No calls	#Mis calls	%Reproducibility
Operator 1	7614	0	0	100.0%
Operator 2	7614	0	0	100.0%
Operator 3	7613	1	0	99.987%
Cumulative	6825	1	0	99.99%

MiSeq Data	Sanger			N/N
	Homozygous Reference	Heterozygous	Homozygous Alternate	
Homozygous Reference	67419	0	0	0
Heterozygous	0	944	0	0
Homozygous Alternate	0	0	153	0
N/N	0	0	0	0

Performance of the Assay in sequencing across Challenging regions was verified through testing 1073 samples; this included runs at low (0.1x recommended DNA input) and high (5x recommended DNA input).


Detection of Challenging Regions

- ▶ Two Large Deletions

Deletion name	Length of Deleted Region
CFTR dele2, 3	21kb
CFTR dele22, 23	1.5kb

PolyTG and PolyT regions in Intron 9

Insertion/deletions in homopolymeric regions



Insertion-cum-deletion in the same region

A 22-bp deletion

Mutation	rs ID	cDNA Name (HGVS)
852del22	rs121908804	c.720_741delAGGGAGAATGATGATGAAGTAC

	N	#Correct Calls	#No calls	#Mis calls	%Reproducibility	% Accuracy
Large Deletion 1	104	104	0	0	100.0%	100.00%
PolyTG/PolyT	1073	1060	13*	0	98.8%	100.00%
Indels in homopolymeric regions (394delTT, 2143delT, 2184delA, 3905insT)	56	56	0	0	100.0%	100.00%
Insertion-cum-deletion (2183delAAA-G)	104	104	0	0	100.0%	100.00%

*12/13 No calls were in the same sample tested at 25ng DNA input (0.1x) and 1/13 was in a sample tested at 1250ng (5x) DNA input.

Results Summary

Parameter	Result
Content	162 variations in CFTR2
Assay Time	7 hr run, 2.5 hr hands-on
Total Run Time (DNA to data)	48hrs
Throughput	48
Accuracy	100%
Reproducibility	> 99.99%
Call Rate (per sample)	> 99.99%
DNA Input	250ng
DNA Source	Blood

Conclusions

The Illumina MiSeqDx Cystic Fibrosis Carrier Screening allows genotyping of 162 variations in the CFTR2 panel in a simple, high throughput workflow with a short hands on time. The high accuracy, reproducibility, and the comprehensive nature of the panel will make it a useful tool to determine a subject's CF carrier status. A broad panel of clinically relevant mutations based upon the CFTR2 database will reduce the residual risk associated with current panels. The design of the instrument and the preliminary performance of the Illumina assay make this system appropriate for clinical use.*

* In development, not available for commercial sale.

© 2012 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.