

cnvPartition v3.2.0 CNV Analysis Plug-in for GenomeStudio® Software

1. General Information

The cnvPartition CNV Analysis Plug-in is a software library that works with Illumina's GenomeStudio data analysis software. This document provides information about downloading, installing, and using the cnvPartition CNV Analysis Plug-in.

2. Downloading and Installing the cnvPartition CNV Analysis Plug-in

1. Download the cnvPartition CNV Analysis Plug-in from one of three places:


- The GenomeStudio Portal
- The software downloads section of [iCom](#)
- The [illumina•connect](#) web page

2. Run the setup program for the plug-in.
3. Follow the instructions in the installation prompts.

The setup program places the *.dll and configuration file into this GenomeStudio directory on your computer:

C:\Program Files (x86)\Illumina\GenomeStudio\CNVAlgorithm\cnvPartition

3. Using the cnvPartition Plug-in with GenomeStudio Software

1. Open a genotyping project using GenomeStudio software.
2. [\[Optional\]](#) Adjust the cnvPartition parameters as needed (see Section 5 of this document).
3. Select **Analysis | CNV Analysis**.
The CNV Analysis dialog appears.
4. Select **Create New CNV Analysis | cnvPartition 3.2.0**.
5. In the CNV Analysis Name area, enter a name for this CNV analysis.
6. [\[Optional\]](#) Adjust the cnvPartition parameters as needed (see Section 5 of this document).
7. Click **Calculate New CNV Analysis**.
8. When the analysis is complete, click **OK**. After running the analysis using GenomeStudio software, you can view the results in many ways:
 - In the **Full Data Table**, in the CNV Value and CNV Confidence columns
Use the column chooser  to display the CNV Value and CNV Confidence columns.
 - In the CNV Region Display tool, by selecting **Analysis | Show CNV Region Display**.
 - In the Illumina Genome Viewer (IGV), by selecting **View | CNV Analysis as Bookmarks**.
 - In the Bookmark Viewer, by going to the IGV and selecting **View | Bookmark Viewer**.

4. Configuring the cnvPartition Plug-in

Starting with v1.2.0 of the cnvPartition Plug-in, the parameters can be changed via an editable configuration file. When you install the cnvPartition Plug-in, a configuration file is installed on your computer in this location:

C:\Program Files (x86)\Illumina\GenomeStudio\CNVAlgorithm\cnvPartition

When you modify and save the default configuration file, your changes are preserved for future sessions.

Note: You can also adjust cnvPartition parameters from within GenomeStudio software.

Perform the following steps to adjust cnvPartition Plug-in parameters from within GenomeStudio:

1. In a GenomeStudio genotyping project, go to **Analysis | CNV Analysis**.
The CNV Analysis dialog appears.
2. Select **cnvPartition 3.2.0** from the drop-down menu.
3. In the **Options** area of the CNV Analysis dialog, modify the parameters as needed.

The following table lists cnvPartition parameters, descriptions, and related information for GenomeStudio software.

Parameter Name	Description	Default Value	Required	Visible in GUI?	Applies to	
					LOH Regions	CNV Regions
AdjustYLRR	Set this flag to true to adjust Log R Ratios downward for Y SNPs.	True	True	False	False	True
AveragePloidy AdjustLRR	Set this flag to true to turn on Average Ploidy Adjustment of the Log R Ratio.	False	True	False	False	True
AveragePloidy AdjustLRR Threshold	The threshold for activating the Average Ploidy Adjustment.	2.5	True	False	False	True
ChiSquare Threshold	Minimum Chi Square value needed to detect a region of extended homozygosity. See the DNA Copy Number Analysis Algorithms Technical Note for more details.	23.5	True	False	True	False
Confidence Threshold	CNV regions with a confidence threshold below this value are excluded from analysis results.	35	True	True	False	True
CopyNeutral LOHOnly	If set to true, detects regions of extended homozygosity only if they are copy neutral.	True	True	False	True	False
DetectExtended Homozygosity	Detect regions of extended homozygosity (regions exhibiting LOH but no change in copy number.)	True	True	True	True	False

		Applies to				
Parameter Name	Description	Default Value	Required	Visible in GUI?	LOH Regions	CNV Regions
	<p>Note: The algorithm for detecting regions of homozygosity is similar to that used in the LOH Detector Autobookmarking Plug-in, and is described in detail in the DNA Copy Number Analysis Algorithms Technical Note. This algorithm is run only on regions with a copy number of 2 as detected by cnvPartition.</p> <p>Additionally, the calculated ChiSquare for the homozygous region is used to populate the CNV Confidence value reported in the Full Data table, and is not subject to filtering by the ConfidenceThreshold parameter. Homozygous regions can be filtered independently of CNV regions using the MinHomozygousRegionSize, ChiSquareThreshold and/or MinHomozygoteCount parameters. CNV regions can be filtered independently of homozygous regions using the ConfidenceThreshold and/or MinProbeCount parameters.</p> <p>To run only CNV detection, toggle DetectExtendedHomozygosity to False. To run LOH detection on Copy Number = 2 regions only, set CopyNeutralLOHOnly to 2.</p> <p>To detect only extended homozygous regions, set ConfidenceThreshold to a very high value, such as 10,000,000. Because even phenotypically normal samples contain many regions of homozygosity, you may want to set this threshold higher, to minimize the number of detected regions.</p>					
Exclude IntensityOnly	Set this flag to true if you want to exclude intensity-only probes from your analysis.	False	True	True	True	True
ExcludeProbes WithAuxValue	Enables you to filter a probe by setting its Aux value in the SNP table to -1. All probes with an Aux value of -1 will be excluded.	-1	True	False	True	True
GapSize Threshold	Regions within probe gaps larger than this value are not considered to be within CNV regions. This helps prevent CNVs from being called across large probe gaps, such as centromeres.	1,000,000	True	False	False	True
GcFile	<p>Enables you to specify the full path to the optional GC content file.</p> <p>Note: The GC content file is a tab-delimited file containing two columns for each probe, with no header. The content is organized as follows:</p> <p style="text-align: center;">Probe name<tab>GC content</p> <p>Where the GC content is a numeric value in the range of 0 to 1. The GC content file is not used by default.</p>	[blank]	False	False	False	True
GcWaveAdjust GCBins	The number of GC value bins used in the polynomial-based GC wave adjustment method.	50	True	False	False	True
GcWaveAdjust LRR	Set this flag to true if you want to adjust Log R Ratios for GC waves.	True	True	True	False	True

Parameter Name	Description	Default Value	Required	Visible in GUI?	Applies to	
					LOH Regions	CNV Regions
GcWave Adjustment Window	The sliding window size, in bases, used by the linear-based GC wave adjustment method.	115000	True	False	False	True
GcWaveAdjust Type	Specifies whether to use a polynomial- or linear-based method for the GC wave adjustment. Valid options are Polynomial or Linear.	Polynomial	True	False	False	True
GenTrain Threshold	Probes with a GenTrain score below this value will not be included in the analysis.	0.00	True	True	True	True
Include Mitochondria	Set this flag to True if you want to include mitochondrial chromosomes in this analysis.	False	True	False	True	True
IncludeSex Chromosomes	Set this flag to True if you want to include X, Y, and XY (pseudoautosomal) chromosomes in this analysis.	True	True	True	True	True
LogDiagnostic Info	Log diagnostic info to a file in the application's temp file directory. The log directory location varies depending on your operating system: Windows XP – C:\Documents and Settings\<username>\Local Settings\Application Data\Illumina\cnvPartition Windows 7 – C:\Users\<username>\AppData\Local\Illumina\cnvPartition	False	True	False	True	True
MinHomozygote Count	Minimum number of homozygotes needed to detect a region of extended homozygosity. See the DNA Copy Number Analysis Algorithms Technical Note for more details.	50	True	True	True	False
MinHomozygous RegionSize	Homozygous regions smaller than this are not detected.	1 Mb	True	True	True	False
MinProbeCount	Regions with probe counts smaller than this are removed from analysis results.	3	True	False	False	True
Smoothing Moving AveragePeriod	The smoothing moving average period.	2	True	False	False	True

Parameter Name	Description	Default Value	Required	Visible in GUI?	Applies to	
					LOH Regions	CNV Regions
SmoothLRR	Set this flag to True if you want to smooth the Log R Ratio.	False	True	False	False	True
UseGcFile	Set this flag to True if you want to use an external GC content file for GC wave correction.	False	True	False	False	True
UseSimple HeterosomyFix	This parameter is still in development and should not be used at this time. Keep this value set at False to leave it off.	False	True	False	False	True

5. Technical Support

Direct questions about installing and using the cnvPartition CNV Analysis Plug-in to Illumina Technical Support at techsupport@illumina.com, 1.800.809.4566 (toll-free), or +1.858.202.4566 (outside North America).

6. Version History

Build	Date
v3.2.0 for GenomeStudio Software	12/2/2011
<ul style="list-style-type: none"> Added compatibility with Illumina HumanOmni5M arrays. Simplified statistical methods to improve speed and memory utilization. This might lead to small changes in CNV boundaries and/or detection of lower confidence CNV regions. 	
Build	Date
v3.1.6 for GenomeStudio Software	3/22/2011
<ul style="list-style-type: none"> Added ability to skip intensity-only probes. Added ability to skip probes with Aux value = -1 in the SNP table. Now call CN=2 on the X chromosome for males. Added LRR adjustment for Y chromosome probes. Now ignore LOH on the Y chromosome. GC wave correction can now be toggled on/off from within the CNV Analysis dialog box. Lowered the minimum homozygous region size to 1 Mb. Added parameter GcWaveAdjustType, set to Polynomial by default. 	
Build	Date
v2.4.4 for GenomeStudio and KaryoStudio Software	11/11/2009
<ul style="list-style-type: none"> Improved consistency between 32-bit and 64-bit PCs—Fixed an occasional rounding issue with 32-bit PCs. 	

- Improved consistency between GenomeStudio software and KaryoStudio software—The upcoming release of KaryoStudio also includes this fix.
- Improved detection of LOH regions—There are now fewer false positives. Also, the min LOH region size was changed to 10 Mb (from 1 Mb) for consistency with KaryoStudio software default settings.
- Added Log R Ratio smoothing—This feature is disabled by default (you can enable it via the config file).

The following experimental features have been added. These features are also disabled by default:

- You can now run LOH detection for the entire genome
 - Previously, you could run LOH detection only on regions for which CN=2
 - In the config file, `CopyNeutralLOHOnly` is set to True by default
- Optional adjustment to Log R Ratio for the Y chromosome
 - Y SNPs are clustered using only males, so the Log R Ratio appears as CN=2
 - To lower all Log R Ratios on Y, turn adjustment on
 - If Y chromosome SNP clusters are already adjusted, additional adjustment could result in inconsistent results
 - In the config file, use `AdjustYLRR`
- Support for highly amplified genomes has been added
 - This is a common situation with cancer samples
 - Log R Ratios are adjusted upward based on calculation of average genomic ploidy
 - In the config file, use `AveragePloidyAdjustLRR` and `AveragePloidyAdjustLRRThreshold`
- Optional GC Wave Adjustment has been added
 - This is based on linear regression of LRR vs. GC content in probes
 - In the config file, use `GcWaveAdjustLRR` and `GcWaveAdjustmentWindow`

Build	Date
v2.3.4 for GenomeStudio and KaryoStudio Software	03/09/2009
<ul style="list-style-type: none"> • Increased accuracy for calling smaller-sized CNVs. • Improved logic for detecting CNVs on the sex chromosomes. Based on gender information from the sample sheet, single copies of X chromosomes are ignored for males. 	
Build	Date
v2.2.1 for GenomeStudio and KaryoStudio Software	01/22/2009
<ul style="list-style-type: none"> • <code>cnvPartition</code> v2.2.1 is compatible with Illumina's GenomeStudio and KaryoStudio software. <p>Note: There are two different versions of <code>cnvPartition</code> v2.2 for KaryoStudio—one for 32-bit systems and one for 64-bit systems. Be sure to install the version of <code>cnvPartition</code> that is compatible with your computer.</p> • The <code>cnvPartition</code> algorithm now ignores Y-chromosome probes for females. • The <code>cnvPartition</code> algorithm now ignores probes with NaN or a missing value for the Log R Ratio. • Changed the minimum value for <code>MinHomozygousRegionSize</code> from 10Mb to 5Mb 	

Build	Date
v1.2.0 for BeadStudio	10/10/2008
<ul style="list-style-type: none"> Added configuration file All chromosomes starting with letter "M" are now recognized as mitochondrial. B Allele Frequency is now ignored for Intensity Only probes. Added configuration option to detect extended regions of homozygosity such as copy-neutral LOH (default = True) Added configuration option for the minimum size of a detected region of homozygosity Added configuration option to ignore CNV regions with less than a minimum number of probes defined by the user (default=3) 	
Build	Date
v1.0.2 for BeadStudio	03/11/2008
<ul style="list-style-type: none"> Added Confidence Threshold parameter—any CNV region with a confidence lower than this value will be removed from the CNV output. The recommended default is 35. Added Probe Gap Size Threshold parameter—a CNV region overlapping with a probe gap larger than this value will not be created. Helps prevent CNV regions from being found across centromeres and other large probe gaps. The recommended default is 1,000,000 base pairs. The internal expected Log R Ratio mean and standard deviation values for each copy number value have been changed to better reflect actual Log R Ratio values seen in normal HapMap samples. This results in higher overall accuracy and precision. 	
Build	Date
v1.0.1 for BeadStudio	01/29/2008
<ul style="list-style-type: none"> Improved detection of smaller CNVs Both Log R Ratio and B Allele Frequency are used as inputs. Mitochondrial SNPs can be included in an analysis (optionally). The CNV Confidence score is now calculated. 	
Build	Date
v0.9.3 for BeadStudio	09/04/2007
<ul style="list-style-type: none"> Improved classification of CN = 0 regions (homozygous deletions) Lowered the Min SNPs Per Bin parameter default to ten (from 30) to improve detection of CN = 0 regions (homozygous deletions) 	
Build	Date
v0.9.2 for BeadStudio	07/23/2007
<ul style="list-style-type: none"> Initial release 	

For Research Use Only

© 2008–2011 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSPPro, GenomeStudio, Genetic Energy, HiSeq, HiScan, Eco, TruSeq, and MiSeq are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

The cnvPartition CNV Analysis Plug-in incorporates portions of code from the Cephes Math Library.

Cephes Math Library Release 2.8: June, 2000 Copyright 1984, 1987, 1995, 2000 by Stephen L. Moshier

Contributors: * Sergey Bochkano (ALGLIB project). Translation from C to pseudocode.

See subroutines comments for additional copyrights.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

-Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.

-Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer listed in this license in the documentation and/or other materials provided with the distribution.

-Neither the name of the copyright holders nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.