

Figure 5: Targeted RNA-Seq Analysis in BaseSpace Cloud—View clustered heat maps with the TruSeq Targeted RNA v1.0 App in BaseSpace.

Simplified Bioinformatics

Data analysis with the MiniSeq System requires no informatics expertise or command-line experience. The MiniSeq System features Local Run Manager software, an on-instrument system for creating a run, monitoring status, and analyzing sequencing data. With Local Run Manager, on-instrument data analysis can be automatically performed upon completion of the sequencing run. The data analysis modules generate simple reports for a wide range of sequencing applications. The modular design allows users to install and update individual analysis modules as needed.

In addition, sequencing data generated with the MiniSeq System can be instantly transferred, stored, and analyzed in the BaseSpace computing environment (cloud-based or onsite). BaseSpace RNA-Seq software Apps provide expert-preferred data analysis tools packaged in an intuitive, click-and-go user interface designed for informatics novices. These Apps deliver optimized analysis pipelines that support a range of common RNA-Seq data analysis needs such as alignment, variant calling, and more. The TruSeq Targeted RNA v1.0 App² aligns targeted sequence reads using a banded Smith-Waterman alignment, quantifies the relative expression of genes and isoforms between several samples, and compares abundance across input samples (Figure 5).

For downstream analysis, BaseSpace Apps generate output files that can be directly input into a broad range of data analysis tools. The BaseSpace Environment includes a growing community of developers who use and provide software tools for visualization, analysis, and sharing. This NGS ecosystem provides one of the largest collections of commercial and open-source analysis tools currently available.

Expression Analysis of FFPE Samples

RNA isolated from FFPE samples can vary in quality among different specimens, or within different samples from the same specimen. These samples can prove challenging for conventional expression analysis methods such as RT-PCR or microarray analysis. TruSeq




Figure 6: Targeted Sequencing with FFPE Samples—Gene expression profiles for FFPE breast cancer samples prepared with the TruSeq Targeted RNA expression kit. $R^2 = 0.987$.

Targeted RNA Expression Kits offer an accurate and powerful method for analyzing low-quality or FFPE samples (Figure 6). With the ability to profile 12–2000 targets per sample, and to multiplex up to 384 samples per run, the TruSeq Targeted RNA Expression Kits enables efficient and cost-effective analysis of common FFPE tissues such as tumor samples or tumor-normal paired samples.

Targeted RNA-Seq vs Expression Microarrays

Compared to microarrays, RT-PCR, and other traditional technologies, NGS targeted RNA-Seq offers several advantages. Targeted RNA-Seq allows deep sequencing, which provides higher sensitivity and a virtually unlimited dynamic range. This sensitivity enables superior performance for the detection and analysis of low-abundance, low-quality, or FFPE samples. Furthermore, sequencing provides in-depth, base-by-base coverage over the area of interest, which offers a more comprehensive approach for discovery applications or for the measurement of rare transcripts compared to traditional methods.

Summary

The MiniSeq Targeted RNA Sequencing Solution offers a highly sensitive and accurate method for analyzing specific transcripts of interest. By harnessing the broad dynamic range of NGS sequencing, researchers can obtain more sensitive and accurate measurements for their targeted gene expression studies. Whether looking for the speed of a fixed panel or the flexibility of a custom panel, the MiniSeq Targeted RNA Sequencing Solution delivers high-quality NGS data in a smaller, more accessible platform.

Learn More

For more on DesignStudio and Targeted RNA Expression Experiments, go to www.illumina.com/content/dam/illumina-marketing/documents/products/technotes/truseq_targ_rna_design.pdf.

To learn more about expression analysis of FFPE samples, visit: www.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote-expression-analysis-ffpe-samples.pdf.

For more on the TruSeq Targeted RNA Expression fixed panels, go to: www.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet_truseq_targeted_rna_expression.pdf.

