

Instrument Flexibility

TruSeq Synthetic Long-Read DNA libraries can be sequenced on Illumina HiSeq® 2500 or HiSeq 2000 Systems. The TruSeq method enables researchers to access valuable long-read information by leveraging a single sequencing instrument for multiple applications, without requiring additional specialized equipment. This method provides greater insight into the genome at a fraction of the cost of conventional approaches.

Simple Analysis and Assembly

Push-button analysis in the BaseSpace environment simplifies assembly of long reads. Data can be transferred from an Illumina sequencing instrument to the BaseSpace cloud instantly. Designed for use with the TruSeq Synthetic Long-Read DNA Library Prep and Barcode Kits, the TruSeq Long-Read Assembly App⁵ constructs long sequences from shorter sequencing reads. The intuitive user interface (Figure 5) simplifies data analysis so that researchers can analyze data simply by selecting the project and file destination. The app then processes the short reads, assembles the initial contigs using overlap-based methods, and finally creates contig scaffolds to generate synthetically long reads (Figure 6). The assembled reads are exported in standard FASTQ format, which can be imported directly into downstream assembly tools for further analysis. The TruSeq Long-Read Assembly App enables one-click informatics for novice users, without requiring extensive expertise or infrastructure.


Summary

By combining shorter reads, TruSeq Synthetic Long-Read technology enables researchers to leverage a single sequencer and obtain more accurate long reads than conventional approaches. Designed with biologists in mind, the TruSeq Long-Read Assembly App simplifies bioinformatics so that researchers can spend less time analyzing data and more time focusing on their research. TruSeq Synthetic Long-Read technology provides a comprehensive solution for genome assembly and genome finishing.

References


1. TruSeq | Illumina (www.illumina.com/truseq) Accessed 06 June 2014.
2. Public Data - BaseSpace (basespace.illumina.com/datacentral) Accessed 06 June 2014.
3. Gurevich A, Saveliev V, Vyahhi N, Tesler G (2013) QUAST: quality assessment tool for genome assemblies. *Bioinformatics* 15: 1072–5.
4. Kurtz S, Phillippy A, Delcher AL, Smoot M, Shumway M, et al. (2004) Versatile and open software for comparing large genomes. *Genome Biol* 5: R12.
5. BaseSpace (basespace.illumina.com) Accessed 06 June 2014.

Figure 5: Simplified Analysis



Featuring an intuitive user interface, the TruSeq Long-Read Assembly App simplifies data analysis for any biological researcher, regardless of bioinformatics expertise. Results are displayed in straightforward tabular and graphical formats.

Figure 6: TruSeq Long-Read Assembly App Workflow



The TruSeq Long-Read Assembly App uses sequencing data to construct long reads, delivering results in a separate FASTQ file.

Ordering Information

Product	Catalog No.
TruSeq Synthetic Long-Read DNA Library Prep Kit (4 samples)	FC-126-1001
TruSeq Synthetic Long-Read DNA Barcode Kit (1 sample)	FC-126-1002
TruSeq Synthetic Long-Read DNA Barcode Kit (4 samples)	FC-126-1003
TruSeq Synthetic Long-Read DNA Accessory Kit	FC-126-1004

