

TruSight® Myeloid Sequencing Panel

Using expert-defined content and proven next-generation sequencing technology to identify somatic mutations efficiently and cost-effectively in hematological malignancies.

Highlights

- **Expert-Defined Content**
Designed by a consortium of recognized experts to target 54 genes mutated frequently in myeloid malignancies
- **Streamlined, Comprehensive Method**
Single workflow includes library preparation, sequencing, data analysis, and data annotation
- **Cost-Effective, Time-Efficient Solution**
Assess multiple genes simultaneously for approximately the same cost as a single-gene assay
- **High-Accuracy, High-Sensitivity Analysis**
Limit of detection down to 5% mutant allele frequency with 500x minimum coverage of each region

Introduction

Blood cancers affect more than 1 million people in the United States alone¹. Current methods for assessing these myeloid malignancies can be effective, but are time-consuming and expensive when looking at multiple variants, and may not determine the underlying genetic cause of the disease.

The TruSight Myeloid Sequencing Panel uses next-generation sequencing (NGS) technology to provide a comprehensive assessment of 54 genes (tumor suppressor genes and oncogenic hotspots) in one assay (Table 1). The panel targets mutations with known involvement

in acute myeloid leukemia (AML), myelodysplastic syndrome (MDS), myeloproliferative neoplasms (MPN), chronic myelogenous leukemia (CML), chronic myelomonocytic leukemia (CMML), and juvenile myelomonocytic leukemia (JMML). The result is a single assay for accurate, economical, and rapid profiling of liquid tumors.

Targeted Content Design Strategy

A consortium of recognized experts in blood cancer disorders designed the content for TruSight Myeloid. The sequencing panel specifically targets 54 genes known in the peer-reviewed literature to be frequently mutated in hematological malignancies, focusing on leukemia and myeloproliferative disorders. Targeted genes include those with cited clinical relevance, providing a comprehensive picture of the disease and its progression. Targeted genes include those involved in the MDS-AML continuum.

Superior Coverage

TruSight Myeloid features a highly optimized oligo pool specific for investigating genomic changes associated with hematological malignancies. The panel focuses on ~141 kb of genomic content, consisting of 568 amplicons of ~250 bp in length designed against the human NCBI37/hg19 reference genome. The oligo pool targets 15 full genes (exons only) plus exonic hotspots of an additional 39 genes, providing nearly 100% coverage of all targeted regions (Table 2). This optimized oligo pool provides uniform coverage of the target regions, enabling > 500x coverage for > 95% of amplicons at > 5,000x

Table 1: Gene Regions Assessed by the TruSight Myeloid Sequencing Panel

Gene	Target Region (exon)	Gene	Target Region (exon)	Gene	Target Region (exon)	Gene	Target Region (exon)
ABL1	4-6	DNMT3A	full	KDM6A	full	RAD21	full
ASXL1	12	ETV6/TEL	full	KIT	2, 8-11, 13 + 17	RUNX1	full
ATRX	8-10 and 17-31	EZH2	full	KRAS	2 + 3	SETBP1	4 (partial)
BCOR	full	FBXW7	9 + 10 + 11	MLL	5-8	SF3B1	13-16
BCORL1	full	FLT3	14 + 15 + 20	MPL	10	SMC1A	2, 11, 16 + 17
BRAF	15	GATA1	2	MYD88	3-5	SMC3	10, 13, 19, 23, 25 + 28
CALR	9	GATA2	2-6	NOTCH1	26-28 + 34	SRSF2	1
CBL	8 + 9	GNAS	8 + 9	NPM1	12	STAG2	full
CBLB	9, 10	HRAS	2 + 3	NRAS	2 + 3	TET2	3-11
CBLC	9, 10	IDH1	4	PDGFRA	12, 14, 18	TP53	2-11
CDKN2A	full	IDH2	4	PHF6	full	U2AF1	2 + 6
CEBPA	full	IKZF1	full	PTEN	5 + 7	WT1	7 + 9
CSF3R	14-17	JAK2	12 + 14	PTPN11	3 + 13	ZRSR2	full
CUX1	full	JAK3	13				

mean coverage. This translates into 8 samples per run, providing the required sensitivity and accuracy to call rare variants with confidence.

Table 2: Coverage Details

Cumulative target region size	~ 141 kb
Number of target genes	54
Amplicon size	~ 250 bp
Number of amplicons	568
Recommended mean coverage	5000x
Target minimum coverage	500x
Percent exons covered at 500x	95
Sample Throughput - MiniSeq™ High Throughput Kit	8 samples/run
Sample Throughput - MiSeq® v3 chemistry	8 samples/run

Simple, Integrated Workflow

The TruSight Myeloid Sequencing Panel offers a fully integrated DNA-to-data solution, including a streamlined workflow and automated data analysis with specific variant calling (Figure 1). Starting with 50 ng DNA isolated from blood, bone marrow, or fine needle aspirates (FNA), libraries are generated with the highly multiplexed oligonucleotide probes. Sample-specific indexes are added to each library, enabling pooling of libraries before sequencing for higher throughput. Pooled libraries are loaded into an Illumina desktop sequencing system for automated sequencing and data analysis. These desktop sequencers include the MiniSeq™, MiSeq®, and MiSeqDx (run in RUO mode). Analyzed data can be imported into VariantStudio software for accurate variant annotation, classification, and reporting. The entire process is completed in just 3 days.




Figure 1: Simple, Integrated Workflow—TruSight Myeloid offers integrated library preparation, sequencing, and automated data analysis, creating a streamlined workflow that goes from DNA to data in 3 days.

Optimized Assay Chemistry

The TruSight Myeloid assay begins with hybridizing a highly multiplexed pool of oligo pairs, one upstream, and one downstream of each region of interest (Figure 2). Each oligo contains unique, target-specific sequence, and a universal adapter sequence that is used in a subsequent amplification reaction. A proprietary extension–ligation reaction extends across the region of interest, followed by ligation to unite the 2 probes and yield a library of new templates with common ends. The resulting extension-ligation templates are PCR-amplified, incorporating 2 unique, library-specific indexes. Final reaction products are converted to a single-stranded, adapter-ligated normalized library using a bead-based protocol. The sequence-ready library can be loaded into the MiSeq reagent cartridge ready for sequencing on the MiSeq system without additional processing.




Figure 2: Optimized Assay Chemistry—The TruSight Myeloid assay enables simple, streamlined hybridization and amplification of targeted regions.

The TruSight Myeloid Sequencing Panel can also be used with the high-throughput, easy-to-use NextSeq 500 Desktop Sequencer. Learn more at www.illumina.com/nextseq.



Illumina • 1.800.809.4566 toll-free (US) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2016 Illumina, Inc. All rights reserved. Illumina, TruSight, MiniSeq, MiSeq, BaseSpace, TruSeq, and the pumpkin orange color are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. Pub. No. 1170-2013-004 Current as of 19 January 2016



AGAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC
TCAAACGTACCGTAACGAACGTATCATTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC
CGACGAAAAGAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTAC
AACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC
AGAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC
GATTACTTGATCCACTGATTCAAACGTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC
CGTATCAATTGAGACTAAATATAACGTACCATTAAAGATTACTTGATCCAAGTATTCAAACGTACCGTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCGTTAACCTTAAAGATTACTTGATCCACTGATTCAAC