


Figure 2: Infinium HD Assay Chemistry



values in the industry, which maximize the likelihood of finding true associations for a given phenotype. By strategically selecting the most powerful tag SNPs, Illumina scientists can ensure maximum power to identify associations, while reducing the redundant information on each BeadChip.

Assay Chemistry

The Infinium Assay has a simple, streamlined workflow (Figure 2). This workflow, as well as the data acquisition and analysis workflows, are common across all products. The Infinium Assay protocol features plate-based sample preparation and whole-genome amplification


without PCR or ligation steps, significantly reducing labor and sample handling errors. Following sample DNA hybridization to the Bead-Chip, two-step allele detection provides high call rates and accuracy. Selectivity and specificity are accomplished in two steps. Target hybridization to bead-bound 50-mer oligos ensures high selectivity while enzymatic single-base extension provides powerful specificity. The single-base extension also incorporates a labeled nucleotide for assay readout. The staining reagent is optimized to provide a higher signal, and more balanced intensities between red and green channels. These features contribute to industry-leading accuracy, high call rates, and copy number data with lower background signal.

Genotype Calling

The Infinium Assay produces two-color readouts (one color per allele) for each SNP in a genotyping study. Intensity values for each color channel, A and B, convey information about the allelic ratio at a single genomic locus (Figure 3). Typical studies incorporate values for large numbers of samples (hundreds to tens of thousands) to ensure significant statistical representation. When these values are appropriately normalized and plotted, distinct patterns (or clusters) emerge, in which samples that have identical genotypes at a given locus exhibit similar signal profiles (A and B values) and aggregate into clusters. For diploid organisms, bi-allelic loci are expected to exhibit three clusters (AA, AB, and BB).

Genotype calls are based on information derived from a standard cluster file, which provides statistical data from a representative sample set. This enables genotypes to be called by referencing assay signal intensities against known data for a given locus. Because the call accuracy is tied to the quality of the cluster data, having an efficient and robust clustering algorithm is essential for accurate genotyping. The proven Illumina Gentrain2 algorithm accurately and efficiently identifies cluster patterns of genotyping samples and reports summary statistics. These statistics are used for downstream genotype calling and CNV analysis. The Gentrain2 algorithm is further described at www.illumina.com/documents/products/technotes/technote_gentrain2.pdf.

Figure 3: Typical Diploid Genotyping Plot



Three clusters of points can be seen for this example locus. The red, purple, and blue regions represent the AA, AB, and BB clusters, respectively.

Ordering Information

Product	Catalog No.	Product	Catalog No.
HumanOmniExpress-24 v1.0 DNA Analysis BeadChip Kits		HumanOmniExpress-24+ v1.0 DNA Analysis BeadChip Kits	
48 samples	WG-312-3001	48 sample	WG-312-3005
288 sample	WG-312-3002	288 sample	WG-312-3006
1152 sample	WG-312-3003	1152 sample	WG-312-3007
HumanOmni2.5-8 v1.1 DNA Analysis BeadChip Kits		HumanOmni2.5-8+ v1.1 DNA Analysis BeadChip Kits	
16 sample	WG-312-2511	16 sample	WG-312-2521
48 sample	WG-312-2512	48 sample	WG-312-2522
96 sample	WG-312-2513	96 sample	WG-312-2523
384 sample	WG-312-2514	384 sample	WG-312-2524
HumanOmni5-Quad v1.0 DNA Analysis BeadChip Kits		HumanOmni5-Quad+ v1.0 DNA Analysis BeadChip Kits	
16 sample	WG-311-5001	16 sample	WG-311-5005
48 sample	WG-311-5002	48 sample	WG-311-5006
96 sample	WG-311-5003	96 sample	WG-311-5007
384 sample	WG-311-5004	384 sample	WG-311-5008

Data Analysis Software

Illumina GenomeStudio® Data Analysis Software offers integrated genotyping and copy number tools as well as the Genome Viewer for graphic visualization. GenomeStudio has an open plug-in interface to integrate third-party applications for more downstream data analysis options. Beeline Software provides a direct path to project creation and sample management for large array experiments. The time required for data analysis is reduced by flexible allele calling and data filtering prior to entry into GenomeStudio. Learn more about GenomeStudio and Beeline by visiting support.illumina.com/array/array_software.ilmn. The illumina•Connect program leverages this open architecture and has made numerous plug-ins available to support genotyping and copy number analysis. The illumina•Connect plug-ins are further described and available for download at www.illumina.com/software/illumina_connect.ilmn.

BeadArray Scanners and Automation Systems

Omni micorarrays are compatible with the Illumina iScan and HiScan systems. These cutting-edge array scanners feature high-performance lasers and powerful optical systems that enable rapid scan times and precise assay detection. A convenient modular design enables researchers to easily build out the systems for evolving research needs.

Illumina, Inc. • 5200 Illumina Way, San Diego, CA 92122 USA • 1.800.809.4566 toll-free • 1.858.202.4566 tel • illumina.com

FOR RESEARCH USE ONLY

© 2010–2013 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GALx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, Sentrix, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub. No. 370-2010-027 Current as of 08 October 2013

An optional Laboratory Information Management System (LIMS) is available to accurately and efficiently track samples. For more information on Illumina LIMS, visit www.illumina.com/software/lims.ilmn. Robotic automation capabilities can be added to improve throughput for labs processing large numbers of samples. More information on robotic automation of Illumina microarrays can be found at www.illumina.com/products/iscan_system_automation_option_packages.ilmn.

Illumina FastTrack Microarray Services

Illumina FastTrack Microarray Services are available to analyze samples in a timely fashion, at a reasonable cost, using any Infinium DNA Analysis BeadChip. This option allows researchers to acquire high-quality data for limited studies or before purchasing their own equipment. For more information on Illumina FastTrack Services, visit www.illumina.com/services.ilmn.

Learn More

To learn more about the Omni Array Family, visit www.illumina.com/applications/genotyping/omni_family.ilmn.

