

Simplified Probe Design

Design TruSeq Custom Amplicon probes targeting user-defined regions using DesignStudio, a free, easy-to-use online software tool that provides optimized coverage and estimates total project pricing (Figure 2). After logging on to a personalized account and naming the project, researchers can select coordinates for targeting genomic regions of interest. Probe design is automatically performed by using an algorithm that considers a range of factors, including GC content, specificity, probe interaction, and coverage. Candidate amplicons are visualized and assessed using estimated success scores. Probes can be filtered with user-defined tags, and then added to, or removed from, the design. After visualization and QC, the CAT design can then be ordered along with suggested accessory kits.

Typical design success, depending upon custom-selected content, is 90% or better for desired bases (Table 1). In some cases, higher or lower success rates can be achieved, depending upon the sequence context (eg, regions of homology or %GC content). To use DesignStudio, customers will need an account on the Illumina ecommerce system. Request an account at icom.illumina.com/Account/Register.

Optimized Assay Chemistry

The TruSeq Custom Amplicon v1.5 assay chemistry begins with hybridizing 2 custom-designed probes upstream and downstream of the region of interest (Figure 3). Each probe contains unique, target-specific sequence as well as a universal adapter sequence that is used in a subsequent amplification reaction. The addition of an improved hybridization buffer, OHS2, ensures that the probes are optimally bound across regions of varying GC-richness and improves assay performance in potentially difficult- to-address regions (Figure 4). A proprietary extension-ligation reaction extends across the region of interest, followed by ligation to unite the 2 probes and yield a library of new template molecules with common ends. During this step, an enhanced extension-ligation mix, ELM4, increases library yield and uniformity (Table 2).

While the increase in uniformity may appear modest, many amplicons that are poorly represented when using ELM3 achieve significantly higher read counts with ELM4, enabling higher confidence for variant calling in these regions. Extension-ligation templates are PCR-amplified and 2 unique, sample-specific indexes are incorporated. Final reaction products are converted to a single-stranded, adapter-ligated, normalized library using a bead-based protocol. The sequence-ready library can be stored in the optimized LNS2 storage buffer for up to 3 weeks, or loaded directly onto the MiSeq System without additional processing.[†]




Figure 2: DesignStudio for Custom Probe Design— The Manage Targets Screen shows a Target Region view in DesignStudio software. Easily visualize targeted genomic regions and attempted amplicons to assess design coverage and score. Summary metrics for the entire project are on the left pane, along with project information and user-defined labels for convenient data sorting during the design phase.

Table 1: DesignStudio Parameters for Custom Amplicon

Parameters	Project A	Project B	Project C	Project D
Total Sequence Input (bp)	34,767	69,423	205,441	4,572
Total Number of Amplicons	293	264	1,434	114
Designability Score (% Base Coverage)	98.3	97.0	93.2	91.3

[†] Also compatible with the MiniSeq, NextSeq Series, and HiSeq Series Systems.




Figure 7: Illumina Amplicon Viewer—Streamlined, intuitive analysis toolbars in Amplicon Viewer. Top pane: visualize coverage levels of each amplicon or across entire genome regions (shown) and view quality scores, as well as variant calls in the middle and lower panes, respectively. Left scroll bars allow users to navigate through samples and targets, and track variants.

After the sequencing run on the MiSeq System, data are automatically aligned using the MiSeq Reporter and can be visualized using the Illumina Amplicon Viewer. As shown in Figure 7, large genomic regions containing multiple amplicons can be viewed, or zoomed in to view individual amplicons. Quality scores are tracked, and variant calls are easily and intuitively visualized. Target and sample statistics are managed in simple drop-down menus. From design to analysis, the TruSeq Custom Amplicon user experience is customized and streamlined, and keeps project data highly accessible.

Summary

TruSeq Custom Amplicon v1.5 is a powerful targeted resequencing solution that provides highly multiplexed amplicon sequencing performance with an optimized chemistry workflow and integrated analysis. Paired with the MiSeq System,[§] The TruSeq Custom Amplicon v1.5 assay produces remarkably high-quality data for up to 96 samples in a single run. The robust combination of the TruSeq Custom Amplicon assay and the MiSeq System has driven studies comprising tens of thousands of samples, including hundreds of thousands of amplicons. Using TruSeq Custom Amplicon v1.5 with the MiSeq System enables the highest quality, most reliable, and simplest approach to deep interrogation of genomic regions of interest.

TruSeq Custom Amplicon Specifications

Specification	Value
Minimum input DNA	50 ng
Attempted amplicons per reaction ^a	16–1,536
Content range ^a	2–650 kb
Specificity ^a	> 70%
Coverage uniformity (> 0.2x mean) ^a	> 80%
Species available	Human, mouse, rat, and bovine
Samples types	Fresh, frozen, and FFPE
Amplicon sizes	150, 175, 250, and 425 bp

a. Target values will vary due to custom designs.

§ Also compatible with the MiniSeq, NextSeq Series, and HiSeq Series Systems.

AATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTAA
ACGTACCGTAACGAACTATCAATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGC AACGAAAGAAAGAAATGATAACAGTAACACACTTCTGTTAACCTTAAG
ACGAAAGAAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCG
GTAACCTTAAGAGCTACCGTGC AACGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAAGTAACATTAAAGAGCTACCGTGC AACGAAAGAAATGATAAC
AATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGT
TACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGC AACGAAAGAAATGATAACAGTA
TATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGC AACGAAAGAAATGATAACAGTA

