


Figure 3: Coverage Uniformity by Workflow



Comparison of the sequencing depth versus total number of bases covered demonstrates that coverage uniformity is consistent across all three workflows.

To examine variant calling performance more closely, specificity, sensitivity, and array concordance to the Infinium® HumanOmni2.5M BeadChip® were measured for both insert size conditions (Table 2). Variant calling specificity measurements assume inheritance patterns based on the samples forming a familial trio (two parents and a child). Except for *de novo* mutations in the child, any variant identified in the child should be also called in a minimum of one parent. Conflict rate is based on comparison to expected Mendelian inheritance patterns and is an established proxy for the false positive rate. Sensitivity was measured as the ability to detect a set of well-characterized variants for NA12878, as reported previously⁷. Array concordance was measured as the agreement between SNV calls in the sequencing data versus a curated set of high-confidence calls made using a high-density microarray. These metrics indicate that variant calling is comparable for both insert size libraries.

Workflow Effects on Coverage Uniformity and Variant Calling

To examine the impact of the new WGS workflows on coverage uniformity, sequencing data from all three workflows (WGS PCR-Free Analysis Workflow 2.0, the new WGS PCR-Free Analysis Workflow 3.0, and the new WGS Low-Input Analysis Workflow v3.0) were compared. By plotting the sequencing depth versus the total number of bases covered for each data set, the data show that uniformity is comparable across all workflows (Figure 3).

Workflow-related differences in variant calling were determined by comparing various metrics for SNVs and indels, including call rate across all reference positions, the total number of variant calls, the ratio of heterozygous to homozygous variants, the fraction of variants not found in dbSNP, and the transition to transversion ratio (Table 3).

Table 3: SNV and Indel Call Quality Statistics by Workflow

Variant Class	Quality Metric	WGS PCR-Free Analysis Workflow v2.0	New WGS PCR-Free Analysis Workflow v3.0	New WGS Low-Input Analysis Workflow v3.0
	Average Coverage Depth	35.3x	32.5x	30.9x
SNV	Call Rate (%)	95.47	94.78	94.24
	Total SNVs	3,510,092	3,496,305	3,459,373
	Ti/Tv	2.08	2.08	2.09
	Novelty Rate	3.28	3.29	3.19
	Het/Hom	1.6	1.62	1.63
Indel	Total Indels	368,466	523,681	487,965
	Novelty Rate	6.24	6.72	6.32
	Het/Hom	1.95	2.03	2.45

Metrics calculated as an average across the CEPH trio NA12891, NA12892, and NA12878. Total quality results depend on the average coverage depth.

Call Rate: Percent of non-N reference genome in which a reference or non-reference call was made for both alleles

Total SNVs: Total number of SNVs that have 'Pass' value in the FILTER key of the VCF file

Ti/Tv: Transition to Transversion ratio of SNV calls

Novelty Rate: Percent of SNVs or indels not found in dbSNP132

Het/Hom Rate: Heterozygous to Homozygous ratio of SNV or indel calls

Total Indels: Total number of indels that have 'Pass' value in the FILTER key of the VCF file

Table 4: Specificity, Sensitivity, and Array Concordance by Workflow

Detection Metric	Variant Class	WGS PCR-Free Analysis-Workflow v2.0	New WGS PCR-Free Analysis Workflow v3.0	New WGS Low-Input Analysis Workflow v3.0
Array Concordance	SNV	99.42%	99.43%	99.43%
	Indel	41.72%	52.92%	50.58%
Sensitivity	SNV	92.87%	93.55%	93.58%
	Indel	99.87%	99.88%	99.88%
Specificity	SNV	99.87%	99.88%	99.88%
	Indel	98.55%	98.19%	97.62%

Array Concordance: Genome concordance with calls from HumanOmni2.5M array calculated as an average across the CEPH trio

Sensitivity: Recovery rate of NA12878 variants previously reported⁷. (95,005 SNVs and 11,403 indels)

Specificity: Mendelian non-conflict rate for the variants called in the CEPH trio.

