

Copy Number Variation (CNV) Analysis with BeadStudio

CNV Analysis is fully supported with BeadStudio software in combination with all Illumina Infinium[®] DNA Analysis BeadChips.

Introduction

Illumina provides many plug-ins designed specifically for copy number variant (CNV) analysis within BeadStudio. The following list describes several plug-ins and where they can be obtained. All of these applications require BeadStudio and internet access. For the latest list of plug-ins, please visit the illumina•Connect website at <http://www.illumina.com/illuminaconnect>.

Table 1: CNV Analysis Options

| Plug-in | Function |
|-----------------------|--|
| cnvPartition | Annotate CNV regions and copy number estimate |
| QuantiSNP | Provide copy number estimates with confidence scores |
| dChip | Provide copy number estimates and LOH scores |
| Homozygosity Detector | Detect extended tracts of homozygosity |
| LOH Score | Provide LOH scores |
| ChromoZone | Annotate CNV regions |

cnvPartition

Summary: Developed by Illumina, this algorithm uses a recursive partitioning approach and provides copy number estimates, as well as chromosomal bookmarks.

Availability: This plug-in is currently available for BeadStudio. It can be downloaded from the BeadStudio portal or from the illumina•Connect plug-ins webpage.

Use: This plug-in is run from the Genotyping Module.

Output: Unlike ChromoZone, cnvPartition is optimized to use intensity information, and thus provides both copy number estimates and bookmarks.

Performance: 25 seconds per sample¹ (Hap550).

QuantiSNP

Summary: Developed at Oxford University, QuantiSNP has been designed and optimized for Illumina SNP genotyping and CNV data. Based on a Hidden-Markov Model², it provides copy number estimates and bookmarks with confidence scores. QuantiSNP is optimal for finding small-sized CNV regions.

Availability: A plug-in is required to call this algorithm in BeadStudio. This plug-in is available for download from the illumina•Connect website or the BeadStudio portal. The QuantiSNP Windows executable is also required and can be obtained by contacting Dr. Ioannis Ragoussis at (ioannisr@well.ox.ac.uk) or from the Oxford University website at <http://www.well.ox.ac.uk/QuantiSNP/>.

Use: The BeadStudio plug-in outputs the correct data input for QuantiSNP. The plug-in is run by selecting QuantiSNP CNV Analysis from the drop-down list of CNV Analysis algorithms in Analysis | CNV Analysis. QuantiSNP results can then be input into BeadStudio in the form of a CNV Analysis. Bookmarks can also be created within the Illumina Genome Viewer (IGV) based on the CNV Analysis.

Output: QuantiSNP provides copy number estimates with confidence scores (p-values).

Performance: Nine minutes per sample using 10 optimization steps¹ (Hap550), 40–60 minutes for Human1M-Duo.

dChip

Summary: dChip was developed at Harvard (Cheng Li laboratory) and has been used on various types of microarray data³. It provides copy number estimates and LOH scores.

Availability: A compatible custom report plug-in is currently available for BeadStudio. This can be obtained from the illumina•Connect website, or the BeadStudio portal. The setup program will install the plug-in in the correct directory. In addition, dChip must be downloaded from <http://www.hsph.harvard.edu/biostats/complab/dchip/download.htm>.

Use: The custom report plug-in will provide the correct data file for input to dChip. The plug-in is run by selecting dChipGTLInput Reports from the drop-down list of Custom Reports in Analysis | Reports | Report Wizard.

Output: dChip provides both copy number estimates and LOH scores.

Note: Illumina only provides a custom report plug-in that outputs BeadStudio data in a dChip compatible format. Illumina does not currently provide any further support of dChip.

Homozygosity Detector

Summary: Developed by Illumina, this algorithm searches for extended regions of homozygosity (LOH) and does not use intensity information.

Availability: This plug-in is currently available for BeadStudio.

