





- Add 18 µl of master mix to all wells containing dilutions, NTCs, or samples. Pipette up and down to mix, carefully avoiding any bubbles.
- Seal the plate, centrifuge briefly, and place into the Real-Time PCR instrument.
- Choose the appropriate program for your experiment from the proposed protocols. These protocols are just starting point guidance and need to be adjusted based on your reagent provider recommendations.

RNA Template		cDNA Template	
Fast Real-Time PCR	Standard Real-Time PCR	Fast Real-Time PCR	Standard Real-Time PCR
50°C for 15 minutes	50°C for 15 minutes	95°C for 2 minutes	95°C for 2 minutes
95°C for 20 seconds	95°C for 2 minutes	40 cycles of: 95°C for 5 seconds	40 cycles of: 95°C for 15 seconds
40 cycles of: 95°C for 5 seconds	40 cycles of: 95°C for 15 seconds	60°C for 15 seconds	60°C for 1 minute
60°C for 15 seconds	60°C for 1 minute		

- On the Plate Setup window specify your (S)tandards as a 5-point, 10-fold serial dilution on wells A1 through E2 as depicted in Table 1. Also specify the (N)o Template Controls (NTC) on wells F1 and F2. Finally if running (U)known samples, specify them in the remaining wells (A3 through F8) of the plate.

- Click Start Run. While the instrument is running you can visualize the amplification in real-time in the Monitor Run window. The run will take approximately 40 minutes (Fast protocol) or up to 90 minutes (Standard protocol).

### Step 5: Data Analysis

Once the run is finished the software automatically opens the Analyze Data window and performs a basic analysis using auto-baseline and threshold settings. Select the Amplification Plot tab. This view of amplification will show spacing of 3.32 cycles between amplification curves for a well-optimized assay. This is because it mathematically takes 3.32 cycles for a 100%-efficient PCR amplification to increase the number of template molecules 10-fold ( $2^{3.32} = 10$ ) in any given sample.

To extract quantification data from Real-Time PCR amplification curves, the results need to be plotted as a linear regression of the Cq values versus the log of the RNA/cDNA quantities. Such a representation is typically called a standard curve. The Eco software automatically generates a standard curve. To view the curve, select the Results tab in the Analyze Data window (Figure 3).

Standard curves provide a wealth of information about a PCR assay. The slope of the line is a measure of the assay's efficiency. Slopes between -3.1 and -3.6 are considered acceptable (90% and 110% efficient, respectively), while a slope of -3.32 is indicative of 100% efficiency.

The R<sup>2</sup> is a measure of the performance of the assay, and is the coefficient of correlation between the data generated and the results expected under ideal conditions. The R<sup>2</sup> should be greater than 0.99.

Figure 2: dRn Amplification Plots of an Example Assay (plus baseline, plus threshold)





