



"For example, data from Affymetrix test sites, which use multiple short oligonucleotide probes per target with perfect match and mismatch sequences, and Illumina test sites, which use plasma-etched silicon wafers containing beads with long oligonucleotide probes, were remarkably similar in detection and detection consistency, gene list overlap and ratio compression analyses. In other words, the expression patterns generated were reflective of biology regardless of the differences in technology."

Excerpt from "MicroArray Quality Control (MAQC) project shows interplatform reproducibility of gene expression measurements."<sup>1</sup>

Taqman real-time PCR assays in terms of correlations of log ratios for expression signals between different samples. Log expression ratios for Samples A and B were calculated for each platform and the rank-based correlation (R) of these ratios was calculated for each cross-platform comparison. As with the comparisons of differential expression in Figure 41, this comparison showed the results from the Affymetrix and Illumina platforms to be most similar to each other, with an R value of 0.936. Other commercial platforms gave R values ranging from 0.792 and 0.933 when compared to Illumina and Affymetrix (Figure 5b1 and Table S111). Comparisons to Taqman reverse transcription PCR showed the results from the Illumina and Affymetrix platforms to be the most similar to this orthogonal technology, with R values of .902 and .905, respectively. The other whole-genome platforms gave R values ranging between 0.839 and 0.894 when compared to Taqman. These results show the Illumina platform to produce results highly consistent with other commonly used commercial array platforms as well as gene expression technologies that are not based on microarrays.


## Conclusions

The outcome of the MAQC project was a successful comparison of multiple whole-genome gene expression profiles across various commercial platforms at an unprecedented scale. This scale, combined with sequence-based probe mapping and a large number of alternative technology gene expression measurements (such as TaqMan assays), sets this study apart from previous cross-platform comparison studies. Despite the use of very different technologies, a relatively high level of interplatform concordance was observed, especially in the results produced by Illumina BeadChips and Affymetrix GeneChips. Given such high levels of concordance, choosing a gene expression platform should primarily be driven by such factors as technical performance, cost, usability, input requirements] (Figure A), and content quality. The Illumina BeadChips, which were shown to consistently be among the very best performers across the various technical measurements, are also substantially less expensive (Figure B) than all other commercial whole-genome arrays used in this study. With array and reagent costs less than half that of the other commercial arrays, the use of Illumina BeadChips allows for experimental designs to be expanded, yielding more powerful and far-reaching results with the same research budget.

## Reference


1. MAQC Consortium. (2006). The MicroArray Quality Control (MAQC) project shows interplatform reproducibility of gene expression measurements. Nature Biotechnology 24(9), 1151-1161.

Figure A: Input RNA Required



The bar graph above shows the RNA required for input into the platform specified on the X-axis.

Figure B: Materials Cost



The approximate cost for materials to run 20 samples (i.e., per site cost) is shown above with the platform specified on the X-axis.

FOR RESEARCH USE ONLY

AGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC  
 TCAACGTACCGTAACGAACGATCATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTGCAACGAACGAAAAAGAATGATAACAGTAACACACTTCTGTTAACCTT  
 CGACGAAAAGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTAAACGTAACGTAACAGTAAAGATTGATA  
 AACGTACCATTAAGAGCTACCGTTGCAACAGTAACACAGTTCGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTGCAACGAACGAAAAGAATGATA  
 AGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC  
 GATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTGCAACGAACGAAAAGAATGATAACAG  
 CGTATCAATTGAGACTAAATATAACGTACCATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTTGCAACGAACGAAAAGAATGATAACAG