



reads consistently achieve an average mapped depth of approximately 40x. The high coverage uniformity indicates that these two runs enable more accurate calling of variants that are distant from the mean depth. Highly uniform coverage also signifies that a sequencing run achieved a high percentage of bases at  $\geq 30\times$  coverage, while a wider distribution would require additional sequencing.

### Reduced Coverage Gaps

The percentage of bases at low read depth quantifies the number of gaps in a data set. A gap exists if a base or several consecutive bases are either not sequenced at all or sequenced with a read depth below the specified cutoff. In such instances, there will be few or no reads to support variant calls (SNPs or indels) in the gap region, resulting in an increased rate of missed calls. Deep coverage of these traditionally difficult-to-sequence regions enables researchers to have more confidence in the variant calls made. Increased coverage can be especially significant when analyzing somatic mutations in the context of certain phenotypic traits or diseases.


Both TruSeq DNA PCR-Free and TruSeq Nano DNA libraries show an approximate 70% decrease in the number of gaps compared to libraries created using the original TruSeq DNA method. This drastic reduction in coverage gaps stems from the reduced bias provided by TruSeq Nano DNA kits and the elimination of PCR-induced bias by the PCR-Free kits. TruSeq DNA PCR-Free and Nano DNA libraries also demonstrate almost 50% reduction in the total size of gaps. Examples of decreased gaps in coverage are shown in Figure 3, using the human *RNPEPL1* (arginyl aminopeptidase-like 1) and the *ZBTB34* (zinc finger and BTB domain-containing 34) promoters.

### Diversity

Diversity refers to the number of unique fragments generated in a sequencing run. High diversity and few duplicate reads provide greater coverage of the genome, increasing the confidence of variant calls as each base is covered by a greater number of unique fragments.

All three TruSeq DNA kits demonstrate similar depths of coverage and produce low percentages of duplicates (Table 1). Diversity levels, reported in billions of fragments, are dependent on the complexity of a given genome. The data in Table 1 demonstrate that the reduced input requirement offered by TruSeq Nano DNA kits and the elimination of PCR by TruSeq DNA PCR-Free kits continue to provide high-quality data. Coverage depth and diversity are consistent across all three kits, and the percentages of duplicates remain low (below 2.5%). The diversity levels provided by all TruSeq DNA kits (a minimum of 2 billion) are ideal for sequencing a human genome at  $\geq 30\times$  coverage without increasing the number of duplicates.

Figure 3: Reduced Gaps in Coverage



The increased genomic coverage provided by TruSeq Nano DNA and TruSeq DNA PCR-Free libraries results in fewer coverage gaps. The examples shown here demonstrate enhanced coverage of the GC-rich coding regions of the *RNPEPL1* promoter (A) and the *ZBTB34* promoter (B). These data were generated using the Integrative Genomics Viewer developed by the Broad Institute<sup>4</sup>.

Table 1: TruSeq Kits Provide High Diversity Levels

Sample	Coverage Depth	Duplicates	Diversity (billions)
TruSeq DNA Kit (300 bp)	38.64	1.88%	13.22
TruSeq DNA PCR-Free Kit (350 bp)	39.76	1.49%	17.19
TruSeq Nano DNA Kit (350 bp)	39.80	2.43%	10.52

All TruSeq DNA kits generate enough unique molecules to sequence a human genome at  $\geq 30\times$  coverage. These data were generated from a  $2 \times 10^1$  bp sequencing run using four lanes of a HiSeq 2500 instrument in rapid-run mode.



