




Figure 2: Effect of Read Depth on Assay Drop-Out



For both high (>20,000 reads/assay) and low (<1,000 reads/assay) read depth, raw data points were normalized by total counts per sample, converted to RP10K units and log2 transformed. At high read depth, only one assay falls below the threshold of 10 reads, while at low read depth, three assays fall below this threshold.

## Normalization

Normalization is an essential step in gene expression analysis. In microarray data analysis, normalization adjusts for systematic biases to allow comparisons of expression levels between and within samples. For RNA-Seq, normalization helps correct for gene length, count distribution between samples, and differences in sequencing depth. The most common method for normalizing RNA-Seq data is by total counts, where counts are divided by the transcript length and multiplied by millions of mapped reads<sup>3</sup>.

Expression ratio accuracy relies on the use of an internal gene by which the TruSeq Targeted RNA data is normalized by sample. This allows for cross-sample and run comparison by adjusting the data to the counts read in the sample and the internal gene. Good candidates for normalization include genes from RNA-Seq data that are relatively unchanged across samples. As an example, the adjustment of expression ratios using such a gene (PER2) is shown in Figure 3. Invariant expression of genes used for normalization should be verified across all samples being studied.

## Conclusions

Although it is a slight departure from other RNA expression methods that routinely include highly expressed housekeeping genes, the read budget concept provides guidelines to maximize study design and take full advantage of MiSeq output. To most efficiently use your read budget for TruSeq Targeted RNA Expression, it is not recommended to use very highly expressed transcripts. A good practice is to include genes that are low to moderate expressors. Assuming an average of 1,000 counts per assay allows for up to 15,000 reactions per run, and sufficient dynamic range to cover most applications.


illumina, Inc. • 9885 Towne Centre Drive, San Diego, CA 92121 USA • 1.800.809.4566 toll-free • 1.858.202.4566 tel • techsupport@illumina.com • illumina.com

### FOR RESEARCH USE ONLY

© 2013-2014 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CSPRO, DASL, Eco, Genetic Energy, GAlIx, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. Pub. No. 470-2013-001 Current as of 20 November 2014

Figure 3: Normalizing with an Invariant Gene



TREx expression ratios normalized to an example of an invariant gene, PER2 (purple) and not normalized (blue) were compared to RNA-Seq expression ratios. Normalization with this gene results in better correlation.

For experiments requiring a very large dynamic range, decrease the number of samples per run to ensure low expressors are detected. See the TruSeq Targeted RNA Expression User Guide<sup>4</sup> for further sample pooling guidelines.

Internal gene-based normalization is recommended for accurate relative quantification, as long as this gene expression is invariant. This step can be performed manually across all samples. See the MiSeq Reporter User Guide<sup>5</sup> for further instructions on how to normalization by total counts is performed.

## References

1. [http://www.illumina.com/systems/miseq/performance\\_specifications.ilmn](http://www.illumina.com/systems/miseq/performance_specifications.ilmn)
2. <http://www.ncbi.nlm.nih.gov/geo/query/acc.cgi?acc=GSE30611>
3. Mortizavi A, Williams BA, McCue K, Schaeffer L, Wold B (2008) Mapping and quantifying mammalian transcriptomes by RNA-Seq. Nature Methods 5;621–28.
4. [http://support.illumina.com/downloads/truseq\\_targeted\\_rna\\_expression\\_guide\\_15034665.ilmn](http://support.illumina.com/downloads/truseq_targeted_rna_expression_guide_15034665.ilmn)
5. [http://support.illumina.com/documents/documentation/Software\\_Documentation/MiSeqReporter/MiSeqReporter\\_UserGuide\\_15028784\\_J.pdf](http://support.illumina.com/documents/documentation/Software_Documentation/MiSeqReporter/MiSeqReporter_UserGuide_15028784_J.pdf)