

Using a PhiX Control for HiSeq[®] Sequencing Runs

A low-concentration spike-in of Illumina PhiX Control v3 provides quality and calibration controls.

Highlights

- For most libraries, Illumina recommends using PhiX Control v3 (Catalog # FC-110-3001) in a low-concentration spike-in (1%) for improved sequencing quality control
- For unbalanced samples with high AT or GC content, use a dedicated PhiX control lane to improve cross-talk and phasing calculations
- For samples with low diversity, use a dedicated control lane and a high-concentration spike-in (40% or higher) of PhiX to create a more diverse set of clusters

What are PhiX Control Libraries?

Control libraries generated from the PhiX virus serve as an effective control in sequencing runs. Characteristics of the PhiX genome provide several benefits:

- **Small**—PhiX is a small genome, which enables quick alignment and estimation of error rates.
- **Diverse**—The PhiX genome contains approximately 45% GC and 55% AT.
- **Well-Defined**—PhiX has a well-defined genome sequence.

Illumina cluster detection algorithms are optimized around a balanced representation of A, T, G, and C nucleotides. Illumina PhiX Control v3 is a balanced and diverse library that can help mitigate sequencing challenges in unbalanced and low diversity libraries.

The mean insert size of the PhiX v3 library is approximately 375 bp, corresponding to approximately a 500 bp library size if visualized on a Bioanalyzer.

Benefits of Using a PhiX Control

PhiX libraries provide a quality control for cluster generation, sequencing, and alignment, and a calibration control for cross-talk matrix generation, phasing, and prephasing.

Cluster Generation

Illumina PhiX Control v3 is shipped as a ready-to-use 10 nM library that can be used as a positive control in the clustering process. If a problem occurs in sample preparation, PhiX will still generate clusters. This helps you discern whether a lack of clusters is due to sample preparation failure or a failure in the cluster generation process.

Cross-Talk Matrix Generation

During an Illumina sequencing run, the cross-talk due to spectral overlap between the four fluorescently-labeled nucleotides is calculated during template generation (cluster identification) in cycles 1–4.

For proper cross-talk calculation, HiSeq Control Software (HCS) v1.3.8, or higher, requires approximately equal numbers and at least 50,000 clusters/mm² in each of the four bases. Therefore, it is imperative to have a balanced representation of bases at the beginning of each read, excluding the Index Read.

Phasing and Prephasing

During sequencing by synthesis, each DNA strand in a cluster extends by one base per cycle. A small proportion of strands may become out of phase with the current cycle, either falling a base behind (phasing) or jumping a base ahead (prephasing). The phasing and prephasing rates define the fraction of molecules that become phased or prephased per cycle.

Calculation of these rates requires a balanced and random base composition in cycles 2–12. Any library that does not comprise a balanced base composition (e.g., initial cycle indexing, restriction enzyme libraries) should use a control lane.

High GC samples typically show higher phasing rates. However, if the sample has good diversity (for example, whole-genome sequencing libraries with 40–60% GC), it does not require a control lane.

Alignment

Because PhiX has a small, well-defined genome sequence, it is an excellent alignment control. If a lane containing PhiX (control lane or low-concentration spike-in) is designated for your run, Real Time Analysis (RTA) software aligns complete sequences to the PhiX reference beginning after the 25th cycle is accumulated and calculates error rates, providing an indication of sequencing success during the run.

Using a PhiX Control for Unbalanced Samples

Most mammalian genomic or whole-transcriptome RNA samples have a balanced genomic composition (approximately equal proportions of A, T, G, and C). These samples do not require a dedicated control lane to generate accurate matrix and phasing estimations. Unbalanced samples contain genomes with high AT or GC content (less than 40% or greater than 60%). Some examples are Arabidopsis, Plasmodium, some bacteria, and bisulfite conversion studies.

Use a dedicated control lane with unbalanced samples to improve cross-talk and phasing calculations.

