TruSight[™] Hereditary Cancer Panel

NGS panel for cancer research featuring 113 genes associated with genetic risk

- Evaluate germline mutations associated with cancer risk using comprehensive, expert-selected panel content
- Enrich and prepare sequencing-ready libraries in 6.5 hours, with only 2 hours hands-on time
- Achieve excellent coverage uniformity for highly accurate detection of SNVs, indels, and CNVs
- Adjust throughput and sequence 2–256 samples per run on Illumina benchtop sequencing systems



Introduction

Genetic variants play an important role in determining cancer predisposition. The TruSight Hereditary Cancer Panel enables researchers to perform a comprehensive evaluation of the genes in which these variants are located. Developed in collaboration with experts in cancer genomics, the TruSight Hereditary Cancer Panel is a targeted sequencing panel designed to assess germline mutations across 113 genes and 125 single nucleotide polymorphisms (SNPs) for variant identification and polygenic risk scoring.

The assay uses predesigned, ready-to-use oligo probes that cover all exonic regions and 20 bp of flanking intronic regions for each targeted gene. Libraries are prepared using hybrid-capture chemistry integrated with Illumina DNA Prep with Enrichment.* Illumina DNA Prep with Enrichment uses an innovative bead-based chemistry with a simplified, single hybridization step for fast and efficient library preparation. Illumina DNA Prep with Enrichment is compatible with all Illumina benchtop sequencing systems, offering flexibility in experimental design across a wide range of sample throughput (Table 1).1 Combining the speed of Illumina DNA Prep with Enrichment with the MiSeg[™] System, the entire workflow (Figure 1), from sample to data, can be completed in 48 hours.

Table 1: TruSight Hereditary Cancer Panel specifications

| Parameter | Details |
|------------------------------|--|
| Supported sequencing systems | iSeq™ 100 System, MiniSeq™ System, MiSeq System, MiSeqDx System (in research mode), NextSeq™ 550 System, NextSeq 550Dx System (in research mode) |
| Panel size | 403 kb, 113 genes (covering all exons), 125 SNPs (48 ID SNPs and 77 SNPs for polygenic risk score) |
| No. of probes | 10,341 oligo probes |
| Sample type | Genomic DNA, blood ^a , or saliva ^a |
| DNA input | 50-1000 ng DNA |
| Total assay time | 48 hours from DNA to data |
| Library prep time | 6.5 hrs total time, 2 hrs hands-on time |
| Sample throughput | 384 indexes available for variable throughput from 2–256 samples per run at average coverage of 300× (minimum coverage 100×) |
| Samples per tube | 8 enrichments (up to 12 samples per enrichment) |

(Illumina, Catalog no. 20018706).




Figure 1: Fast, flexible NGS workflow— The TruSight Hereditary Cancer Panel uses Illumina DNA Prep with Enrichment library prep chemistry, which integrates library preparation and enrichment steps. A fast, streamlined, and optimized workflow delivers fully enriched libraries in just 6.5 hours. TruSight Hereditary Cancer is compatible with the iSeq 100, MiniSeq, MiSeq Series, and NextSeq Systems.

Illumina DNA Prep was formerly known as Nextera™ DNA Flex Library Prep kit. The two kits use the same tagmentation chemistry and have identical product performance specifications and kit configurations.

Flexibility of throughput with Illumina sequencing systems

The TruSight Hereditary Cancer Panel is compatible with multiple Illumina sequencing systems, providing flexibility and control over experimental design. Users can select instruments or reagent kits according to laboratory needs. Sample throughput can range from 2–256 samples per run (Table 2).

Comprehensive content design

The TruSight Hereditary Cancer Panel includes an extensive list of genes commonly associated with hereditary predisposition to breast, colon, ovarian, and gastric cancers (Figure 2). The content was developed with input and feedback from key opinion leaders on genetic risk assessment from Germany, France, and the United Kingdom. The panel includes 10,341 probes that target 113 genes related to cancer predisposition (Table 3), and evaluated on population studies of cases vs. controls. Also included are 48 SNPs for identity and gender determination purposes, and 77 SNPs for BOADICEA polygenic risk score.^{2,3} Analysis enables the detection of single-nucleotide variants (SNVs), insertions/deletions (indels), and copy-number variants (CNVs) in a single assay (Table 4, Table 5).

Table 3: TruSight Hereditary Cancer Panel gene content

| ACD | DIS3L2 | GREM1 | PIK3CA | SDHD |
|--------|--------------|--------|---------|---------|
| AIP | EPCAM | HOXB13 | PMS2 | SLX4 |
| AKT1 | ERCC1 | KIF1B | POLD1 | SMAD4 |
| APC | ERCC2 | KIT | POLE | SMARCA4 |
| ATM | ERCC3 | LZTR1 | POT1 | SMARCB1 |
| BAP1 | ERCC4 | MAX | PRKAR1A | SMARCE1 |
| BARD1 | ERCC5 | MEN1 | PTCH1 | SPINK1 |
| BLM | FAM175A | MET | PTEN | SPRED1 |
| BMPR1A | FANCA | MITF | RAD50 | STK11 |
| BRCA1 | FANCB | MLH1 | RAD51 | SUFU |
| BRCA2 | FANCC | MRE11A | RAD51B | TERF2IP |
| BRIP1 | FANCD2 | MSH2 | RAD51C | TERT |
| CASR | FANCE | MSH3 | RAD51D | TMEM127 |
| CDC73 | FANCF | MSH6 | RB1 | TP53 |
| CDH1 | FANCG | MUTYH | RECQL4 | TSC1 |
| CDK4 | FANCI | NBN | RET | TSC2 |
| CDKN1B | FANCL | NF1 | RHBDF2 | VHL |
| CDKN2A | FANCM | NF2 | RINT1 | WT1 |
| CEBPA | FH | NSD1 | RUNX1 | XPA |
| CHEK2 | FLCN | NTHL1 | SDHA | XPC |
| CTRC | GALNT12 | PALB2 | SDHAF2 | XRCC2 |
| DDB2 | GATA2 | PDGFRA | SDHB | |
| DICER1 | GPC3 | РНОХ2В | SDHC | |
| | | | | |

a. For the complete list of SNPs included in the panel, visit www.illumina.com/ TruSightHereditaryCancer.

Table 2: Sample batching and output variation between instruments and reagent kits

| Sequencing system ^a | Reagent Kit | Single reads | Output | Run time | Sample plexity ^b |
|--------------------------------|-------------|--------------|--------|----------|-----------------------------|
| iSeq 100 System | 100 i1 | 4M | 1.2 Gb | 19 hr | 2 |
| | v2 Micro | 4M | 1.2 Gb | 19 hr | 2 |
| MiSeq and MiSeqDx Systems | v2 Standard | 15M | 4.5 Gb | 24 hr | 9 |
| | v3 Standard | 25M | 7.5 Gb | 28 hr | 16 |
| MiniOn a Octoberra | Mid Output | 8M | 2.4 Gb | 17 hr | 5 |
| MiniSeq System | High Output | 25M | 7.5 Gb | 24 hr | 16 |
| NextSeg 550 and | Mid Output | 130M | 39 Gb | 26 hr | 80 |
| NextSeq 550Dx Systems | High Output | 400M | 120 Gb | 39 hr | 256 |

a. Theoretical outputs and times for the iSeq 100 and MiniSeq Systems are based on instrument specifications. Internal verification for the TruSight Hereditary Cancer Panel was performed on the MiSeq and NextSeq 550 Systems only.

b. Sample throughput is based on 300× average coverage per sample.

Table 4: Variant detection in Horizon Discovery samples using TruSight Hereditary Cancer Panela,b

| | | | | | | Observed N | MAF at varied | DNA inputs ^c |
|--------|-------|---------|--------------|---------------------|--------------|------------|---------------|-------------------------|
| Sample | Gene | Variant | Variant type | Consequence | Expected MAF | 50 ng | 25 ng | 10 ng |
| | BRCA1 | P871L | SNV | missense mutation | 100% | 100% | 100% | 99.8% |
| | BRCA1 | S1613G | SNV | missense mutation | 50% | 49.8% | 47.7% | 45.8% |
| | BRCA1 | K1183R | SNV | missense mutation | 50% | 45.0% | 43.9% | 44.9% |
| | BRCA1 | K820E | SNV | missense mutation | 50% | 48.1% | 43.6% | 45.6% |
| | BRCA1 | D435Y | SNV | missense mutation | 50% | 42.8% | 46.3% | 44.6% |
| HD793 | BRCA2 | V2466A | SNV | missense mutation | 100% | 99.9% | 100% | 100% |
| | BRCA2 | N289H | SNV | missense mutation | 50% | 39.2% | 40.5% | 40.5% |
| | BRCA2 | N991D | SNV | missense mutation | 50% | 48.6% | 48.1% | 48.0% |
| | BRCA2 | N1784fs | Deletion | frameshift mutation | 50% | 42.2% | 35.7% | 38.9% |
| | BRIP1 | S919P | SNV | missense mutation | 100% | 99.7% | 99.9% | 100% |
| | NBN | E185Q | SNV | missense mutation | 50% | 41.1% | 35.1% | 38.5% |
| | BARD1 | R378S | SNV | missense mutation | 50% | 50.5% | 49.9% | 48.0% |
| | BRCA2 | V2466A | SNV | missense mutation | 100% | 99.9% | 99.9% | 99.8% |
| HD794 | BRCA2 | 12675fs | Insertion | frameshift mutation | 50% | 41.0% | 40.9% | 40.3% |
| | BRIP1 | S919P | SNV | missense mutation | 100% | 99.9% | 100% | 100% |
| | NBN | E185Q | SNV | missense mutation | 100% | 100% | 100% | 100% |

a. Sequencing was performed on the MiSeq System.

Table 5: Variant detection in collaboration samples using TruSight Hereditary Cancer Panela,b

| Sample | Gene | Reference allele | Variant allele | Variant type | Consequence | Rep 1 MAF ^c | Rep 2 MAF° |
|--------|---------------|------------------|----------------|--------------|--------------------|------------------------|------------|
| 1 | PALB2 overlap | | | CNV | copy number change | Detected | Detected |
| 0 | RB1 | Т | TTCAAAA | Insertion | In frame insertion | 54.1% | 53.6% |
| 2 | TSC2 | С | Т | SNV | Stop gained | 49.8% | 47.5% |
| 3 | POLE | С | Т | SNV | Missense variant | 44.1% | 47.0% |
| 4 | CHEK2 | А | G | SNV | Missense variant | 40.8% | 44.9% |
| 5 | MSH6 | GA | G | Deletion | Frameshift variant | 50.9% | 45.0% |
| 6 | BRCA2 | CG | С | Deletion | Frameshift variant | 29.9% | 36.3% |
| 7 | MLH1 | С | Т | SNV | Stop gained | 31.0% | 31.9% |
| 8 | BRCA1 | Т | С | SNV | Missense variant | 39.6% | 35.1% |

a. Sequencing was performed on the MiSeq System.

b. Alignment and variant calling were performed with the DRAGEN Enrichment App. Observed minor allele.

c. frequency (MAF) values are mean values from four technical replicates.

b. Alignment and variant calling with the DRAGEN Enrichment App.

c. Observed variant calls correlate with genotypes previously reported by our collaborator (data not shown).

| Cancer type | Recommended genes for screening |
|-------------------|--|
| Breast | ATM, BARD1, BRCA1, BRCA2, CDH1, CHEK2, NBN, NF1, PALB2, PTEN, STK11, TP53 |
| Colon | APC, AXIN2, BMPR1A, CHEK2 EPCAM, GREM1, MLH1 MSH2, MSH6, PMS2, MSH3, MUTYH, NTLH1, POLD1, POLE, PTEN, SMAD4, STK11, TP53 |
| Ovarian | ATM, BARD1, BRCA1, BRCA2, CDH1, CHEK2, NBN, NF1, PALB2, PTEN, STK11, TP53 |
| Gastric | CDH1 |
| Chir Other | MEN1, NF2, RB1, RET, SDHAF2, SDHB, SDHC, SDHD, TSC1/2, VHL, TP53, WT1 |

Figure 2: Genes included that have known associations with genetic predisposition to specific types of cancers.

Fast library preparation and enrichment workflow

The TruSight Hereditary Cancer Panel uses Illumina DNA Prep with Enrichment, enabling fast library preparation, with sequencing-ready libraries ready in 6.5 hours, including only 2 hours hands-on time. A key component of the Illumina DNA Prep with Enrichment solution is on-bead tagmentation, which uses bead-bound transposomes to mediate a uniform tagmentation reaction (Figure 3). This strategy eliminates the need for separate DNA fragmentation steps. For gDNA inputs between 10-50 ng, saturation-based DNA normalization also eliminates the need for individual library quantification and normalization steps before enrichment. Target enrichment occurs through proven hybrid-capture chemistry, enabling reliable detection of relevant variants for SNVs, indels, and CNVs. Libraries are hybridized to biotin-labeled probes specific for targeted DNA regions. Targets are captured by streptavidin magnetic beads that bind to the biotinylated probes, pulling the bound fragments from solution. After captured fragments are eluted from the beads, the targeted library is ready for sequencing.




Figure 3: Illumina DNA Prep with Enrichment workflow—A uniform tagmentation reaction mediated by eBLTs followed by a single hybridization reaction enables a fast and flexible workflow.

Accurate data

With the ability to assess 113 genes per sample, the TruSight Hereditary Cancer Panel provides a high level of sample throughput while maintaining excellent specificity and uniformity. To demonstrate assay performance, sequencing metrics from two sequencing systems were analyzed using research collaborator samples. Eight samples (in duplicate) with 50 ng DNA input were prepared using Illumina DNA Prep with Enrichment with eight-plex enrichments and sequenced on the MiSeg System and the NextSeq System, and data were evaluated using the BaseSpace Enrichment App v3.1.0. Results showed a high percentage of coverage uniformity (Figure 4).




Figure 4: On-target alignment and coverage uniformity—DNA extracted from collaborator samples were prepared using the TruSight Hereditary Cancer Panel and sequenced on the (left) MiSeq System and (right) NextSeq 550 System. Mean values from two technical replicates are shown for each sample.

Variant calling

To demonstrate variant calling performance at different input levels, sets of 16 samples were prepared with 10 ng, 25 ng, and 50 ng DNA inputs. Sample sets were comprised of four replicates each of Horizon Discovery (HD) samples BRCA Germline I Reference Standard gDNA HD793 and BRCA Germline II Reference Standard gDNA HD794. Each input level was sequenced in 16-plex after preparing with Illumina DNA Prep with Enrichment with 8-plex enrichments. Sequencing was performed on the MiSeq System and data was evaluated using the DRAGEN Enrichment App. Results were concordant to the published list for Horizon Discovery for samples HD793 and HD794, demonstrating reproducible results across all input levels tested. Additional analysis was performed on samples containing unknown variants from research collaborators (Table 4). 50 ng DNA input of eight samples in duplicate were prepared using Illumina DNA Prep with Enrichment with eight-plex enrichments and sequenced on the MiSeg System. Using the DRAGEN Enrichment App for data analysis, variants from different classes (SNV, CNV, and indel) were detected (Table 5), which correlated with genotypes previously reported by our collaborator.

DRAGEN Enrichment App or the BaseSpace Enrichment App can be used for variant calling to provide results in VCF format. Customers can select any third-party tertiary analysis platform to annotate and interpret variants.

For more information, including adjustable assay parameters for the Illumina DNA Prep with Enrichment workflow and the impact on Variant Calling, read the User-definable parameters in the Illumina DNA Prep with Enrichment workflow tech note and the Analyze germline CNVs with TruSight Hereditary Cancer Panel tech note.

Summary

The TruSight Hereditary Cancer Panel enables researchers to access an expert-defined content set for analyzing variation within genes previously linked with a predisposition towards cancer. The optimized probe set provides comprehensive coverage of the targeted regions with high coverage uniformity for identifying many variants. Combining this content with the Illumina DNA Prep with Enrichment method enables a fast, easy workflow with a low sample input requirement, and the flexibility of using any Illumina benchtop sequencing system. The TruSight Hereditary Cancer Panel is a highly efficient targeted sequencing solution to accelerate detection of variants associated with cancer predisposition.

Learn more

TruSight Hereditary Cancer Panel

References

- 1. Illumina. Illumina DNA Prep with Enrichment data sheet. Published 2020. Accessed September 5, 2023.
- 2. Mavaddat N, Pharoah PD, Michailidou K, et al. Prediction of breast cancer risk based on profiling with common genetic variants. J Natl Cancer Inst. 2015;107(5):djv036. Published 2015 Apr 8. doi:10.1093/jnci/djv036
- 3. University of Cambridge, Centre for Cancer Genetic Epidemiology. Breast and Ovarian Analysis of Disease Incidence and Carrier Estimation Algorithm (BOADICEA) ccge.medschl. cam.ac.uk/boadicea/. Accessed September 5, 2023.

Ordering information

| Product | Catalog no. |
|---|-------------|
| Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 Indexes, 96 Samples) | 20091654 |
| Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 Indexes, 96 Samples) | 20091656 |
| Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 Indexes, 96 Samples) | 20091658 |
| Illumina DNA/RNA UD Indexes Set D, Tagmentation (96 Indexes, 96 Samples) | 20091660 |
| Illumina DNA Prep with Enrichment, (S) Tagmentation (96 Samples) | 20025524 |
| Illumina DNA Prep with Enrichment, (S) Tagmentation (16 Samples) | 20025523 |
| Illumina DNA Prep, (S) Tagmentation (96 Samples) | 20025520 |
| Illumina DNA Prep, (S) Tagmentation (16 Samples) | 20025519 |
| iSeq 100 i1 Reagent | 20021533 |
| iSeq 100 i1 Reagent | 20021534 |
| MiSeq Reagent Micro Kit v2 | MS-103-1002 |
| MiSeq Reagent Kit v2 | MS-102-2002 |
| MiSeq Reagent Kit v3 | MS-102-3003 |
| MiniSeq Mid Output Kit | FC-420-1004 |
| MiniSeq High Output Kit | FC-420-1003 |
| NextSeq 500/550 Mid Output Kit v2.5 | 20024905 |
| NextSeq 500/550 High Output Kit v2.5 | 20024908 |
| Flex Lysis Reagent Kit | 20018706 |



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel techsupport@illumina.com | www.illumina.com

© 2023 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html. M-GL-02150 v1.0