

emphatically said that they preferred the iScan System,” Dr. McElfresh added. “With the Experienced User Cards and standard operating procedures already developed, the Illumina platform was an easy choice to make.”

Proof-of-Principle Studies Demonstrate Robust Assay Performance

Casework Genetics conducted proof-of-principle studies using mock evidence samples to perform concordance, sensitivity, and preliminary mixture studies. The company was so impressed with the data generated that it went on to complete a full developmental validation of the HumanOmni1-Quad assay according to the ISO 17025 Standard for international accreditation, and the FBI Quality Assurance document¹, which provides the guidelines for forensic operations throughout the United States.

“After looking at a number of platforms, we found Illumina’s to be the most automated and organized assay; meaning that it would be the easiest to adapt for robust forensic testing.”

“Using the HumanOmni1-Quad BeadChip, we accurately detected trace amounts of genomic DNA—down to less than 1% of the total mixture DNA concentration—at overall DNA concentrations of approximately 0.06 ng,” Dr. McElfresh said. “While this is in line with STR systems for overall sensitivity, it extends the limit of detectable mixture contributors.” One particular advantage of BeadChip technology is the redundant assay of each locus. Coupled with the simple bi-allelic structure of SNP loci, this provides a level of sensitivity that allows analysis of data that would fall under the stochastic threshold of STR systems.

Casework Genetics has used the technology to successfully unravel mixed forensic DNA samples, reviving items that failed to yield interpretable results using STR analysis. “For critical DNA mixtures, STRs often lead to a dead-end of inconclusive data,” Dr. McElfresh said. “SNPs are the light at the end of the tunnel in these cases and can lead directly to solving a crime, getting a violent offender off the streets or innocent persons out of the system.”

Resolution of mixed samples has been an issue for the forensic community, even since the shift from restriction fragment length polymorphism (RFLP) to STR analysis, and the situation has only worsened as sensitivity has increased, making analysis of touch samples routine. According to forensic DNA expert Martin Tracey, Ph.D., Professor of Biological Sciences at Florida International University, “there have been no advances in STR methodologies or interpretive algorithms offering sufficient improvement to resolve complex mixtures. Casework Genetics’ SNP approach using the Illumina BeadChip appears to be the resolution of the mixture comparison dilemma in forensics.”

The assay has also been found to be especially applicable to severely degraded genetic material, including crime scene samples found in less than ideal conditions and environments that are hostile to DNA preservation. The HumanOmni1-Quad BeadChip is able to redundantly assay 50 bp fragments, resulting in more DNA available for analysis and increasing the statistical power of the results. In contrast, DNA fragments of miniSTRs range between 70 bp and ~300 bp. “Given the results we are achieving, we believe it is possible to obtain high-quality results from DNA on commonly handled objects in criminal investigations,” Dr. McElfresh said.

Proprietary Software Analyzes Raw SNP Data

To keep up with the scale of data that Casework Genetics is generating, Dr. McElfresh and his team developed Forensic SNP Analysis Software, an automated analysis tool to handle the large data set produced by this type of sample processing. It was built upon fundamental advancements made by a team of researchers at the Translational Genomics Research Institute (TGen)—led by David Craig, Ph.D., Associate Director, Neurogenomics Division—which provided a straightforward solution to rapidly and sensitively exploit raw SNP allelic intensities from the iScan System². According to Dr. Craig, “the same technologies that sparked a revolution in our understanding of disease have the potential to fundamentally change forensics.”

Based on the data generated by the HumanOmni1-Quad BeadChip, the software targets and compares thousands of SNPs on each of the 23 autosomes, on the X and Y chromosomes, and in the mitochondrial genome. It does not require *a priori* knowledge of the number of individuals in the mixed sample or assume the relative amounts of DNA present from each contributor.

“Casework Genetics’ SNP approach using the Illumina BeadChip appears to be the resolution of the mixture comparison dilemma in forensics.”

To speed processing time, the software breaks the massive data set into subsets that are organized by chromosome, then analyzes those sets individually. This enables the software to be installed and run on portable notebook computers. The end result is a comparison of all samples and a complete forensic analysis of the data required to interpret the relationship of the samples found in the mixture. “Our software is successful in taking the massive amount of data generated by Illumina GenomeStudio[®] software and getting it into the forensic format that officers of the court can use effectively” Dr. McElfresh said.

