ClearPack Lite™

Basic tools for microarray washing and hybridization.

Description

The ClearPack Lite provides laboratories with basic tools that are required to complete the hybridization and washing stages of the 24sure[™] and CytoChip[™] Focus protocols.

Supply

Each pack contains the items listed in Table 1.

Storage

The contents of the ClearPack Lite are shipped at room temperature and should be stored in a cool, dark, environment upon receipt.

Quality Control

Pack contents are sourced from reputable suppliers using the same order codes and batches that are used in Illumina R&D and production laboratories. Although not subject to specific quality control, the parallel use of most consumables should provide a higher degree of confidence compared to equivalent products from alternate suppliers and sources.

Stock Control, Delivery Dates, and Expiration

Illumina operates a sophisticated stock control procedure to ensure reliable supplies of critical components while maximizing useful life. To achieve this, each microarray model number is classified as a volume, high-volume, or special product. For each product, Illumina commits to a target delivery date and a guaranteed delivery date as summarized in Table 2.

Table 2: Stock Control, Delivery, and Expiration

Parameter	Description
Catalog No.	PR-70-431001-00
Description	ClearPack Lite
Classification	High volume
Target Delivery	7 days
Guaranteed Delivery	30 days

Additional Documents

Additional documents available for download include:

- Laboratory Reference Manuals: 24sure Dual Channel, CytoChip
- Laboratory Summary Protocols: 24sure Dual Channel, CytoChip

Ordering Information

Product	Catalog No.
ClearPack Lite	PR-70-431001-00

Table 1: ClearPack Lite Contents

Catalog No.	Product	Quantity	Description
PR-70-431001-00 ClearPack Lite	ClearPack Lite	2	Slide-staining rack, stainless steel with modified handles
		3	400 ml lidded square glass staining dishes
		1	Coverslips No. 1, 22 × 22 mm (200 units)
		3	Magnetic follower with pivot ring, 25 × 6 mm
		3	Microscope slide box

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2014-2015 Illumina, Inc. All rights reserved. Illumina, 24sure, ClearPack Lite, CytoChip, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub. No. 1570-2014-021 Current as of 23 November 2015

