

Illumina Connected Insights

Built with security and
compliance at the core

- Enterprise-level, secure genomics data platform designed with data privacy provisions like Health Insurance Portability and Accountability Act (HIPAA), General Data Protection Regulation (GDPR), and other key regulations in mind
- Compliant with global and local data privacy and requirements with ISO/IEC 27001:2013 certifications and guaranteed data residency
- Multilayered, security-first infrastructure built with encryption, two-factor authentication, role-based management of PHI/PII data, and more

Introduction

Our customers trust us with sensitive information. Analyzing, handling, and storing this large-scale genomics data for a diverse range of applications requires enterprise-level protection. To keep Illumina platforms, products, and web applications secure for everyone, we partnered with top-tier cloud providers around the globe and built Illumina Connected Insights with security at the core.

Connected Insights is a customizable platform that streamlines user-defined tertiary analysis and report generation to help labs address data interpretation bottlenecks as they bring next-generation sequencing (NGS) assays in house or scale existing workflows.

Connected Insights can integrate directly with [Illumina Connected Analytics](#) and shares the same enterprise-level data privacy and security standards, empowering labs performing deep data science and supporting data sharing on a secure and compliant platform. Connected Analytics is a comprehensive, cloud-based data management and secondary analysis platform that enables researchers to manage and process large volumes of genomic data in a secure, scalable, and flexible environment.

Data within Connected Insights and Connected Analytics are hosted on Amazon Web Services (AWS) and maintains compliance with a wide variety of industry-accepted security standards using AWS Well-Architected best practices.¹ By committing to local and global security policies, Illumina aims to reduce roadblocks encountered by researchers and clinicians to realize the true potential of genomics data and solutions.

Key security and privacy features

Availability

In the context of cloud-computing services, internal and external availability risks exist. To address these concerns, Illumina built a business continuity and disaster recovery plan into its business process. Connected Insights is installed on a high-availability cloud infrastructure in ISO/IEC* 27001:2013 that adhere to Uptime Institute Tier

*ISO, International Organization for Standardization;
IEC, International Electrotechnical Commission; AES,
Advanced Encryption Standard

III design standards to guarantee dedicated network connectivity, redundancy, uninterruptible power supply (UPS), and effective data backup strategies.

Record keeping and audit logs

Connected Insights allows for record keeping and audit logs, ensuring IT accountability within the platform for virtually all objects, actions, and activities, including viewing an object.

API protection

Connected Insights was built with application programming interface (API) protection in mind. All service methods require API key signatures, and service is refused to all others. Requests are monitored for abuse.

Encryption for sensitive data

Connected Insights prioritizes confidentiality of data processing activities in the cloud environment. Data uploaded from sequencing instruments undergo encryption both "in transit" and "at rest" using the AES standard and transfer layer security (TLS 1.2 or newer).

Third-party penetration testing

Third-party penetration tests simulate an attack on a system's deployment and are a good way to test defenses actively. Illumina employs an unbiased third party to conduct penetration tests for Connected Insights cloud instances. After the vendor finishes the test, Illumina receives a comprehensive report, detailing the results. Illumina does not release the results of these penetration tests.

Data isolation

Connected Insights offers the highest degree of data isolation by implementing industry-standard data segregation techniques, including the need-to-know principle, enforced through technical and organizational measures, eg, role-based access governed by fine-grained security controls.

Data management and retention

A fully automated data management platform, Connected Insights stores customer data synchronously across multiple availability zones within a geographic region, performs regular data integrity checks, and self-heals to protect against data loss.

Table 1: Illumina Connected Insights certifications and accreditations

Certification	Description
ISO/IEC 27001:2013	International standard for managing risks to the security of information; certification to ISO 27001 proves information management; the standard adopts a process-based approach for establishing, implementing, operating, monitoring, maintaining, and continually improving an ISMS
AWS standards and accreditations	
Service Organization Controls 1/SSAE 16/ISAE 3402	An audit verifying that AWS controls to protect customer data are properly designed and that the individual controls are operating effectively.
Federal Information Security Management Act (FISMA) Moderate	An accreditation granted by the US Government to strengthen federal information system security; for reference, the NIH data centers are rated FISMA moderate
ISMS, information security management system; SSAE, Statements on Standards for Attestation Engagements; ISAE, International Standards for Attestation Engagements	

Global standards and certifications

Connected Insights is ISO/IEC 27001:2013 certified by an independent auditor for the full scope of its activities, including development, management, and support of a cloud-based analysis platform.

Connected Insights was developed in accordance with the Illumina Software Life Cycle (SLC) process under the Illumina Quality Management System (QMS). Additionally, processes within the Illumina QMS have adopted industry best practices and relevant standards.

Additionally, Connected Insights complies with DCB 0129, UK's National Health Service clinical risk management standard attained by manufacturers of health IT systems to ensure clinical safety of products.²

Integrity

Connected Insights uses public key infrastructure (PKI); hashing techniques ensure data flow integrity and origination across the entire solution. Customer database backup occurs up to 24 times per day to decrease the risk of data loss. In addition, the system contains logging that provides notification when data are altered. If improper alteration is detected, rolling back to a previous backed up version is available.

Login policies

Connected Insights enforces strong password requirements, a renewal period, an inactivity timeout, and the option to implement single sign-on (SSO). Individualized logins are also available, enabling multiple users to utilize the platform per instance.

Portability

A vendor lock-in does not exist, removing legal impediments to export client data by only the appropriately permissioned client.

Two-factor authentication

An authentication service is supported by Security Assertion Markup Language (SAML) 2.0 to manage institutional users and passwords (optional). Step-up authentication also ensures protection of sensitive actions. Two-factor authentication is a configuration option available enabling customers to set up their own federated access management process

Privacy by design

Connected Insights supports customers operating in regulated environments and is in accordance with current data protection laws, including GDPR and HIPAA[†].

Illumina facilities that process protected health information (PHI) or personally identifiable information (PII) are in compliance with HIPAA and employ industry best practices such as:

- Buildings are monitored 24 hours a day and keycard accessed
- Offices have a monitored security system
- Computers used to access or store PHI are password protected and have full-disk encryption turned on
- Any access from outside the office is via a secure Virtual Private Network (VPN)

Role-based management of sensitive data

Connected Insights supports customers operating in regulated environments with stringent compliance requirements. Connected Insights includes fine-grained, configurable access controls that govern individual user access and management of sensitive PHI/PII data within the platform. To prevent error, data loss, or tampering, system access is restricted based on which roles require access and the tasks those roles are required to complete.

Transparency

Connected Insights complies with most data residency and privacy requirements; data center regions and providers are disclosed.

Shared responsibilities

Illumina is responsible for protecting the infrastructure that runs all of the services offered in the AWS Cloud. This infrastructure is composed of the hardware, software, networking, and facilities that run AWS Cloud services. Part of this responsibility requires that Illumina performs recurring security patch updates or other updates to protect the environment from emerging threats and supports iterative improvements. Illumina provides these updates during windows defined in the Connected Insights Terms and Conditions. Customers required to comply with HIPAA are responsible for ensuring that they have a HIPAA compliance program in place and that they use Connected Insights in a manner to ensure their compliance.



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com | www.illumina.com

© 2023 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html. M-GL-01599 v2.0.

Guaranteed data residency

Connected Insights provides data residency to address local regulatory and compliance requirements and support clinical customers operating in a regulated environment. Built with the same exceptional security and privacy features as [Connected Analytics](#), Connected Insights employs a region-specific instance model where omics data files, metadata, and health data are stored in the region the user selects. In globally distributed, high-performance computing centers, the platform regulates access to the data; the actual omics data flow, including data download and data view, occurs between the browser and the regional web server directly. When collaborating with partners in different regions, users can implement cross-regional access, reducing latency while ensuring data residency. Data centers supporting Connected Insights include[†]:

- US East (N. Virginia) us-east-1
- UK (London) eu-west-2
- Germany (Frankfurt) eu-central-1

Learn more

[Illumina Connected Insights](#)

Reference

1. Cloud Security—Amazon Web Services (AWS). Amazon website. aws.amazon.com/security. Accessed April 12, 2023.
2. DCB 0129: Clinical Risk Management: Its application in the manufacture of Health IT systems – NHS digital. NHS UK. digital.nhs.uk/data-and-information/information-standards/information-standards-and-data-collections-including-extractions/publications-and-notifications/standards-and-collections/dcb0129-clinical-risk-management-its-application-in-the-manufacture-of-health-it-systems. Accessed June 6, 2023.

[†] Instances available at launch or per first customer request.